- Please unfold the following pages containing the photos quoted in the reference.
- Veuillez déplier les pages suivantes contenant les photos citées en référence.
- Por favor, despliegue las páginas siguientes que contengan las fotos referenciadas.
- Vogliate aprire le seguenti pagine, contenenti le foto in oggetto.
- Falten Sie die folgenden Seiten bitte auf, die die in Verweis zitierten Fotos enthalten.
- Sla volgende pagina's open voor de foto's waar in de tekst naar verwezen wordt.
- Por favor desdobre as páginas seguintes que contêm as fotografias citadas em referência.
- Följande sidor visar bildreferanser.
- Lütfen atıfta bulunan resimleri görmek için sayfaları çevirin.

# Moria

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













#65039-J-05.2023






# CONTENTS

<b>1</b>	<b>DISCLAIMER</b> .....	<b>2</b>
1.1	MAINTENANCE AND WARRANTY .....	2
1.2	USE OF GENERIC PRODUCTS AND REUSE OF SINGLE-USE CONSUMABLES .....	2
1.3	FLAP THICKNESS .....	2
<b>2</b>	<b>EQUIPMENT LIST</b> .....	<b>3</b>
<b>3</b>	<b>REGULATORY INFORMATION</b> .....	<b>3</b>
<b>4</b>	<b>LABELING INFORMATION</b> .....	<b>4</b>
<b>5</b>	<b>IMPORTANT INFORMATION</b> .....	<b>5</b>
5.1	DESCRIPTION .....	5
5.2	INDICATIONS .....	5
5.3	CONTRAINDICATIONS .....	5
5.4	WARNINGS .....	5
5.5	POTENTIAL ADVERSE EVENTS .....	5
5.6	PRECAUTIONS .....	6
5.7	PRODUCT COMPLAINTS .....	6
5.8	PACKAGING .....	6
<b>6</b>	<b>INSTALLATION AND CONNECTION</b> .....	<b>7</b>
<b>7</b>	<b>OPERATION</b> .....	<b>8</b>
<b>8</b>	<b>TROUBLESHOOTING</b> .....	<b>9</b>
<b>9</b>	<b>CARE &amp; MAINTENANCE</b> .....	<b>9</b>
9.1	INTRODUCTION .....	9
9.2	INITIAL TREATMENT AND STORAGE .....	9
9.3	CLEANING-DISINFECTION-STORAGE OF THE HANDPIECE AND ELECTRICAL CABLE .....	10
<b>10</b>	<b>WARRANTY</b> .....	<b>11</b>
10.1	SCOPE OF WARRANTY .....	11
10.2	NON-APPLICATION OF WARRANTY .....	11
10.3	WARRANTY PERIOD .....	11
10.4	LIABILITY .....	11
<b>11</b>	<b>DRAWINGS</b> .....	<b>12</b>

# 1 DISCLAIMER

## 1.1 MAINTENANCE AND WARRANTY

The One Use-*Plus* system has been designed for optimal operation, provided that the recommendations listed in this user manual are followed carefully. If, for any reasons, the system does not perform properly, have it checked immediately by MORIA. In order to maintain the original performance of your microkeratome, MORIA strongly recommends annual maintenance and servicing of all its reusable products.

As only MORIA and its agents are fully expert in MORIA products, servicing and maintenance must be carried out by MORIA or its approved agents.

MORIA shall not be held liable for any malfunction or damage to the apparatus, poor results, or surgical complications due to maintenance being having been carried out by an unqualified operator or third party.

Any such unauthorised intervention shall render the guarantee and any maintenance contract null and void.

## 1.2 USE OF GENERIC PRODUCTS AND REUSE OF SINGLE-USE CONSUMABLES

The materials used in the microkeratome blade, blade holder and head have been selected for their sliding characteristics. The dimensions and tolerances of the blade have been determined in consideration of the dimensions and tolerances of the head of the keratome. MORIA's manufacturing and inspection procedures guarantee there is no conflict in dimensions between head and blade, and that the blade will slide smoothly in the head.

**The reuse of single-use products, or the use of consumables other than those supplied by MORIA, may entail serious surgical consequences for the patient and damage the microkeratome.**

MORIA shall not be held liable in the event of a malfunction or damage to the microkeratome, poor results or surgical complications due to the reuse of a single-use product or the use of consumables other than those supplied by MORIA.

All warranties become null and void if the microkeratome degrades or malfunctions due to such practices.

## 1.3 FLAP THICKNESS

Flap thickness is a key factor for both LASIK and lamellar keratoplasty. Multiple parameters affect flap thickness and standard deviation. Numerous scientific studies have demonstrated that flap thickness is influenced by several patient-related factors such as keratometry (K) readings, corneal anatomy, preoperative pachymetry and refraction, intraocular pressure (IOP), and also by surgery-related factors, such as corneal hydration and the speed of the pass (when using manual microkeratomes).

Ultrasonic pachymetry measurements are not always accurate and reproducible, and results may also vary with surgeons' techniques and calibration of the device. Thus, a MORIA One Use-*Plus* head labeled:

- "90" (SBK head) cuts, on average, a 100-micron flap with variations around this average.
- "130" or "130L" cuts, on average, a 130-micron flap with variations around this average.

MORIA can only guarantee the dimensional characteristics of the head, not the surgical result.

For the EVOLUTION 3 and 3E consoles, please refer to the related instruction manuals (#65038, 65051, 65060).

The most recent version of this user guide and additional information on your keratome are available on the MORIA website: <http://www.moria-surgical.com>.


## 2 EQUIPMENT LIST

Designation	MORIA reference
One Use-Plus: ring -1 and Head SBK (90) or 130	19336/90 19336/130
One Use-Plus: ring 0 and Head SBK (90) or 130	19337/90 19337/130
One Use-Plus: ring +1 and Head SBK (90)	19338/90
One Use-Plus: ring -1 Large Cut and Head Large Cut 110L	19354/110
One Use-Plus: ring -1 Large Cut and Head Large Cut 130L	19354/130
One Use-Plus handpiece (blue color)	19345
Wrench	19345C
Storage box for the One Use-Plus handpiece	22519514
Tonometer	19042
EVOLUTION 3 console	19360
EVOLUTION 3E console	19380
Footswitch	19361
Footswitch Epi-K™	19381
Aspiration tubing	19138
Disposable brush	19149
One Use-Plus single-use plastic ring user manual	65039
EVOLUTION 3 console user manual	65038
EVOLUTION 3E console user manual (serial number below 5000)	65051
EVOLUTION 3E console user manual (serial number 5000 and above)	65060
Annexe "Guidance and manufacturer's declaration: electromagnetic emissions and immunity"	65073
One Use-Plus single-use plastic ring nomogram	65041

## 3 REGULATORY INFORMATION

	MORIA S.A. 27, rue du Pied de Fourche - 03160 Bourbon l'Archambault - France Phone: +33 (0) 4 70 67 09 04 Fax: +33 (0) 4 70 67 06 61 moriam@moriam-int.com <a href="http://www.moria-surgical.com">http://www.moria-surgical.com</a>
Customer Service Information	Contact your local dealer or MORIA
EUROPE	 0459 CLASS IIa according to MDD 93/42/CEE
USA	Product registered at the Food and Drug Administration (FDA): 510(k) K040297 Caution USA only: Federal law restricts the use of this device to physicians or licensed practitioners.
Electrical Safety Standard	IEC 60601
	For EU customers only: this symbol indicates that within the European Union, the product must be discarded in a separate collection bin at the end of its useful life. This applies not only to this device, but also all accessories, including the footpedal and electrical motors, regardless of whether those accessories are marked with the symbol. Do not dispose as unsorted municipal waste. For users out of European Union: please refer to local environmental regulatory regarding waste of electrical and electronic equipment.

## 4 LABELING INFORMATION

 REF XXXXXX	CATALOGUE REFERENCE
	QUANTITY
 YYYY/MM/DD	USE BY
 LOT XXXXXX	BATCH CODE
	STERILIZED USING ETHYLENE OXIDE
	DO NOT REUSE
 XXXXXXX	MANUFACTURER
 YYYY/MM/DD	DATE OF MANUFACTURE
 XXXXXX	CAUTION : CONSULT ACCOMPANYING DOCUMENT(S)
 XXXXXX	CONSULT OPERATING INSTRUCTIONS
	DO NOT USE IF PACKAGE IS DAMAGED
 °C °F %	STORAGE CONDITIONS: <ul style="list-style-type: none"> <li>• TEMPERATURE: XX°C – YY°C / XX°F – YY°F</li> <li>• HUMIDITY RATE: XX% – YY%</li> </ul>
	NON-STERILE
	KEEP DRY
	DISCARD IN A SEPARATE COLLECTION BIN
	CAUTION for USA only: US federal laws restrict this device to sale by, or on the order of, a physician.

## 5 IMPORTANT INFORMATION

### 5.1 DESCRIPTION

The One Use-*Plus* microkeratome is used to create corneal, lamellar, hinged flaps (keratectomy). The One Use-*Plus* microkeratome is an automated, mechanical, linear microkeratome designed to create nasal-hinged flaps. It has three components:

- a handpiece (#19345), containing 2 independent motors: one for advancement, and one for oscillation.
- a single-use suction ring. A dedicated nomogram has been developed for One Use-*Plus* with plastic rings (#65041).
- a single-use plastic head encasing a pre-inserted blade.

The One Use-*Plus* microkeratome operates with the EVOLUTION 3 and 3E control units (#19360, 19380). Please refer to user guide:

- #65038 (EVOLUTION 3)
- #65051 (EVOLUTION 3E serial numbers below 5000)
- #65060 (EVOLUTION 3E serial numbers 5000 and above).

### 5.2 INDICATIONS

The One Use-*Plus* microkeratome is indicated for the creation of lamellar nasal-hinged flaps in corneas with preoperative pachymetry of 500 microns or greater, and keratometry between 39 D and 48 D.

During laser in-situ keratomileusis (LASIK), the flap is then lifted to enable photoablation of the stroma with an excimer laser.

### 5.3 CONTRAINDICATIONS

Patients who are not candidates for LASIK.

Note: Careful preoperative evaluation and sound clinical judgment should be used by the surgeon to determine the risk/benefit ratio. Particular attention must be given before performing a keratectomy on a patient with any of the following conditions:

- preoperative pachymetry under 500 microns
- keratometry under 39 D
- keratometry over 48 D
- inability to withstand a transient rise in intraocular pressure.

### 5.4 WARNINGS

- **Design of this medical device does not allow its reuse. Reusing this medical device affects its clinical performances and exposes the patient to adverse events.**
- Do not mix One Use-*Plus* heads with Epi-K™ components of a brand other than MORIA with the One Use-*Plus* microkeratome.
- Do not mix heads and rings from separate sterile packs. Each disposable head must be used with the disposable ring from the same pack.
- The heads must only be screwed on by hand. Dismantling should be done only with the provided wrench (#19345C) or by hand. Never use any tools or other spanners. Incorrect assembly may cause incomplete or uneven cuts due to lack of blade oscillation.
- Never pull on the cable connected to the handpiece, and never hold the motor by the cable.

### 5.5 POTENTIAL ADVERSE EVENTS

As with any surgical procedure, there is risk involved. LASIK surgery requiring the use of a microkeratome that cuts a corneal flap, potential side effects of laser refractive surgery may include but are not limited to: visual anomalies, dry eye and flap related complications (free cap, incomplete flap, buttonhole, epithelial defect, flap dislocation, flap striae, wrinkles, etc.).

Inappropriate use, deterioration of microkeratome and/or non-respect to contraindications (§ 5.3) and warnings (§ 5.4) expose patient at higher risk to adverse events.

Note: Additional surgery may be necessary to correct some of these potential adverse events.

Document Aletiq : Rév. J-05.2023 (Validée)

#65039-J-05.2023

## **5.6 PRECAUTIONS**

- The keratectomy must be performed only by experienced refractive surgeons with specific training in the use of the One Use-*Plus* microkeratome.
- Preoperative and operative procedures, including knowledge of surgical techniques, proper head and ring selection, and assembly and placement of the microkeratome are important considerations in the successful use of the system by the surgeon. Furthermore, the proper patient selection and compliance will greatly affect the results.

### **For USA only**

**CAUTION: FEDERAL LAW (USA) RESTRICTS THESE DEVICES TO SALE BY OR ON THE ORDER OF A PHYSICIAN.**

Other preoperative, intraoperative and postoperative warnings and precautions are as follows:

#### **• PREOPERATIVE:**

- o Only patients who meet the criteria described in the indications should be selected.
- o Care should be used in the handling and storage of the microkeratome components. They should not be scratched or otherwise damaged. Handpieces should be protected during storage, especially from corrosive environments.
- o Check the label and expiration date on the unopened package of heads and rings.
- o After opening the package, verify that head and ring information is consistent with information on outer package labeling.
- o The head, suction ring, handpiece and control unit should be fully inspected prior to use. Both surfaces of the pre-inserted blade should be inspected under a microscope.
- o The surgeon should be familiar with the various components before using the microkeratome and verify that all parts and necessary instruments are present and properly assembled before the surgery begins. Additional sterile components should be available in case of an unexpected need.
- o The selection of the proper head, ring and stop position for each eye is crucial to the success of the procedure: see MORIA nomogram (#65041).

#### **• INTRAOPERATIVE:**

- o Breakage, slippage, or misuse of microkeratome and its components may cause injury to the patient or operating personnel.
- o Before any usage, lubricate the eye, the suction ring, the head and the blade with a physiological saline solution or another appropriate ophthalmic solution that is compatible with the One Use-*Plus* components.

#### **• POSTOPERATIVE:**

- o The surgeon's postoperative instructions to the patient and the corresponding patient compliance are extremely important.

## **5.7 PRODUCT COMPLAINTS**

Health care professionals (e.g., customers or users of this device) having any complaints about or dissatisfaction with the product quality, identity, durability, reliability, safety, effectiveness and/or performance, should notify MORIA or its distributor by telephone, fax or written correspondence, and have the products checked by MORIA.

When filing a complaint, please provide the component(s) name(s), reference(s), batch number(s), as well as your name and address, the nature of the complaint, and the relevant patient data. Also disinfect and return the components.

## **5.8 PACKAGING**

Packages of each of the components should be intact upon receipt. If a loaner or consignment system is used, all sets should be carefully checked for completeness, and all components should be carefully checked to ensure that there is no damage prior to use. Damaged packages or products should not be used and should be returned to MORIA.

## 6 INSTALLATION AND CONNECTION

Steps	What to do	 <b>IMPORTANT WARNINGS</b>	Related picture
1	Select the suction ring and the head according to the One Use-Plus nomogram	<ul style="list-style-type: none"> <li>• The disposable head and suction ring are supplied sterile and for single use only.</li> <li>• Check that package is undamaged, unopened and within the expiration date.</li> <li>• <b>WARNING: DO NOT MIX HEADS &amp; RINGS OF DIFFERENT PACKS.</b></li> <li>• <b>EACH HEAD MUST BE USED WITH THE RING FROM THE SAME PACK.</b></li> </ul>	A, B
2	Assemble the head onto the handpiece	<ul style="list-style-type: none"> <li>• Open the package containing the disposable head and suction ring.</li> <li>• Screw the head clockwise onto the threaded motor shaft.</li> <li>• Check that the head is screwed on tightly to the motor. There should be no rotation of the head.</li> <li>• If the head is not screwed completely onto the motor, the shaft will not drive the blade holder, and the blade will not move. Poor or irregular resections can result from incorrect assembly.</li> <li>• Do not use the wrench (#19345C) or other tools to assemble it. The wrench should be used for disassembly only.</li> </ul>	E, F
3	Carefully inspect the keratome head after assembly	<ul style="list-style-type: none"> <li>• The disposable head must be perfectly clean inside and outside and free from debris, particles, oxidation and deposits. There should be no scratches on the plates. If any, replace the entire pack (ring and head).</li> <li>• After assembly and using a microscope, carefully inspect the blade under high magnification to ensure that it is not damaged.</li> </ul>	A
4	Check for blade oscillation	<ul style="list-style-type: none"> <li>• Start the motor to check proper blade oscillation</li> <li>• Do not use the One Use-Plus system if the oscillation is not smooth, regular and uninterrupted.</li> </ul>	
5	Inspect the suction ring	<ul style="list-style-type: none"> <li>• The suction ring must be carefully inspected prior to use; it must be perfectly clean, free from debris, particles and deposits.</li> <li>• Carefully check that the suction ports are not obstructed.</li> <li>• Inspect the tubing of the ring, and in case of kinks or obstructions, replace the entire pack (ring and head).</li> </ul>	B
6	Set the stopper	<ul style="list-style-type: none"> <li>• Refer to the indicative nomogram (#65041) for the selection of the ring and stop.</li> <li>• The adjustable stop pin is pre-mounted on the suction ring and used to determine the size of the hinge. The hinge will be placed in a nasal position.</li> <li>• Use the plastic spanner included in the One Use-Plus pack (ring and head) to select the appropriate stop value.</li> <li>• The value of the stop selected must line up with the mark on the suction ring.</li> </ul>	C
7	Connect the suction ring to the aspiration tubing	<ul style="list-style-type: none"> <li>• The aspiration tubing (#19138) is supplied sterile and for single use only.</li> <li>• Check that package is undamaged, unopened and within the expiration date.</li> <li>• Inspect the aspiration tubing and in case of kinks or obstructions, replace it.</li> </ul>	D

## 7 OPERATION

Steps	What to do	 <b>IMPORTANT WARNINGS</b>	Related picture
1	Preliminary check before any usage	<ul style="list-style-type: none"> <li>It is imperative to verify that the devices are perfectly clean and free of any organic or other residues.</li> <li>Also, please refer to EVOLUTION 3 or 3E user manual (#65038, 65051, 65060) for operations prior use. The operator should be familiar with the functioning of the EVOLUTION 3 or 3E control unit (#19360, 19380).</li> <li>If you have additional questions, consult the MORIA website or your MORIA distributor.</li> </ul>	
2	Connect the motor to the control unit, select One Use-Plus and the speed, then check vacuum prior use	<ul style="list-style-type: none"> <li>See EVOLUTION 3 or 3E console user manual (#65038, 65051, 65060).</li> <li>Refer to the indicative nomogram for the selection of speed.</li> </ul>	
3	Insert the head in the rails of the suction ring	<ul style="list-style-type: none"> <li>When inserting the head in the rails of the suction ring, the handpiece must be at a 45° angle from the threaded shaft of the ring, to avoid any blade damage</li> <li>If the bottom of the head (the area where the blade edge protrudes) is in contact with the threaded shaft of the suction ring, inadvertent damage to the blade may occur. Release vacuum and change to a completely new set to avoid the possibility of irregular cut.</li> <li>Once the head properly inserted, slightly reduce the angle from 45° to 0°</li> <li>Carefully engage the head into the ring</li> <li>The shaft of the ring must be aligned with the motor.</li> </ul>	G, H
4	Check for keratome translation prior to use (forward and backward pass)	<ul style="list-style-type: none"> <li>Press vacuum ON, make forward then backward movements using footpedals, then press vacuum OFF. Please refer to EVOLUTION 3 or 3E user manual (#65038, 65051, 65060)</li> <li>Do not use the One Use-Plus system if the forward pass is not smooth, regular and uninterrupted.</li> </ul>	
5	Position the suction ring on the eye. Activate vacuum by pressing the "Vacuum" footswitch once.	<ul style="list-style-type: none"> <li>Verify that the device will not interfere with any other medical devices used during the procedure.</li> </ul>	
6	Check the intraocular pressure with the tonometer	<ul style="list-style-type: none"> <li>The tonometer (#19042) must be perfectly dry and used only on dry eyes. If the pressure is below 65 mm Hg, do not proceed with the surgery.</li> </ul>	
7	<b>Lubricate the suction ring, the head and the blade</b>	<ul style="list-style-type: none"> <li><b>Prior to use, lubricate the rails of suction ring, head, and blade with balanced salt solution or another appropriate solution.</b></li> <li><b>The device should not be placed in contact with solutions for which the potential interactions are unknown.</b></li> </ul>	
8	Once the One Use-Plus system is in position, activate the motor by pressing the "Forward" footswitch. When the microkeratome head touches the stop, immediately release the "Forward" footswitch. Reverse by pressing the "Backward" footswitch.	<ul style="list-style-type: none"> <li>See EVOLUTION 3 or 3E console user manual (#65038, 65051, 65060).</li> <li>Hold the suction ring by its handle; check that your fingers do not hinder the forward movement of the device.</li> <li>Check nothing prevents or modifies head movement in the rails of the suction ring</li> <li>Ensure that there are no obstructions (speculum, eyelashes, eyelids, etc.) in its path.</li> </ul>	
9	Release the "Vacuum" footswitch by pressing once. Remove the suction ring.	<ul style="list-style-type: none"> <li>Discard the disposable head, ring and aspiration tubing after each procedure in an appropriate container.</li> </ul>	

## 8 TROUBLESHOOTING

For any further information and description, please refer to EVOLUTION 3 or 3E console user manual (#65038, 65051, 65060).

In case of defective handpiece, contact your MORIA representative or distributor for repairing.

## 9 CARE & MAINTENANCE

### 9.1 INTRODUCTION

In case of unusual vibrations or noises, do not use the unit and contact your distributor.

Please contact MORIA for any other servicing.

The following are only recommendations. They must be adapted according to the laws in the country of utilization. All the products and solutions for cleaning and disinfection must be used according to the manufacturer's instructions.

The recommendations for cleaning, disinfection and sterilization are updated and available on MORIA's website: <http://www.moria-surgical.com>.

MORIA recommends that you:

- avoid touching potentially contaminated areas,
- use gloves when performing cleaning and decontamination tasks.

### 9.2 INITIAL TREATMENT AND STORAGE

To avoid risks of condensation inside the packaging, unpack and store the devices in a clean, dry environment.

Do not store the instruments in an environment or next to other products that may possibly have a corrosive or magnetic effect. Avoid any contact between devices, especially those comprised of different materials.

The device should be damage free and have no scratches or other surface defects.

Fragile devices or those requiring particular manipulation must be handled separately, taking special care to protect delicate parts.

### 9.3 CLEANING-DISINFECTION-STORAGE OF THE HANDPIECE AND ELECTRICAL CABLE

Steps	What to do	 <b>IMPORTANT WARNINGS</b>	Related picture
1	<p><b>Cleaning</b></p> <ul style="list-style-type: none"> <li>To clean the screw core of the motor and metallic parts, use the disposable brush (#19149) moistened with cleaning solutions.</li> <li>Dip the blade oscillation shaft of the motor in a sterile distilled water bath and run the motor back and forth 10 times.</li> <li>Remove the motor from the cup, unplug it from the console, and wipe the oscillation shaft of the motor with a lint-free cloth.</li> <li>Cable can be clean with a lint-free cloth moistened with cleaning solutions.</li> <li>Air-dry the shaft of the motor with clean, filtered compressed air (medical compressed air).</li> </ul>	<ul style="list-style-type: none"> <li><b>The handpiece must be carefully cleaned after each surgical procedure.</b></li> <li>Do not use abrasives or scrapers to clean any elements of the system. Doing so may result in reduced precision, burrs and/or irregular tissue resection.</li> <li>Do not immerse the motor in the sterile distilled water bath.</li> <li>To prevent damaging the connectors and to ensure proper functioning of the motor, never pull on the cables and never hold the motor by its cable.</li> </ul>	I
2	<p><b>Disinfecting &amp; Drying</b></p> <ul style="list-style-type: none"> <li>Wipe each part of the motor with a lint-free cloth moistened with disinfecting solutions.</li> <li>Use appropriate disinfecting solutions (spray or wipe clean disinfectant) according to the manufacturers' instructions.</li> <li>Dry carefully with disposable instrument cleaner (lint-free) then with clean, filtered, compressed air.</li> </ul>	<ul style="list-style-type: none"> <li>It is imperative to verify that the devices are perfectly clean and free of any organic or other residues.</li> <li>Do not gas sterilize (ETO) the motor.</li> <li>Do not autoclave the motor.</li> <li>No current sterilisation process is compatible with MORIA motors.</li> </ul>	
3	<p><b>Storage</b></p>	<ul style="list-style-type: none"> <li>Failure to completely dry the inside of the motor can result in oxidation.</li> <li>Do not store products that are not completely dry; doing so may create rust and result in irregular corneal resection.</li> <li>When not in use, the One Use-Plus system should be kept in its storage box #22519514 in a dry atmosphere.</li> </ul>	

# 10 WARRANTY

## 10.1 SCOPE OF WARRANTY

Designation	MORIA references
One Use-Plus handpiece	19345
Wrench	19345C
Storage box for the One Use-Plus handpiece	22519514
EVOLUTION 3 console (except the battery)	19360
EVOLUTION 3E console (except the battery)	19380
Footswitch	19361
Footswitch Epi-K™	19381

- The above items as well as spare parts and labor necessary for their repair are covered by warranty. Any items returned must be sent in their original packaging, after having previously been disinfected.
- The maintenance operations and the replacement of spare parts will be exclusively carried out by technicians authorized by MORIA.
- Any defective part exchanged during the period of warranty becomes property of MORIA.
- Instruments and accessories that cannot be reused are excluded from this warranty.
- Future upgrades and/or improvements on the keratome are not implied by this warranty.

## 10.2 NON-APPLICATION OF WARRANTY

The warranty will not be applicable under any of the following conditions:

- Defects or malfunction that occur out of the warranty period (10.3).
- Normal wear and tear.
- Negligence or usage that does not comply with the specifications in the user manual.
- The use of supplies, spare parts, or accessories other than those supplied by MORIA.
- Any disassembly, modification or intervention carried out on the devices by a person not authorized by MORIA.

## 10.3 WARRANTY PERIOD

- The warranty takes effect on the date the material is dispatched.
- The duration of the warranty is 12 months from the date of effect.

## 10.4 LIABILITY

- The liability of MORIA is limited to the supply of the services mentioned in paragraph 10.1. MORIA will not be held responsible for any direct or indirect damage suffered by the client owing to the interventions within the scope of this warranty.
- For any dispute concerning the interpretation or the execution of the present contract or the present general terms and conditions, the Commercial Court of Nanterre (France) will have sole jurisdiction.

# 11 DRAWINGS

## A. DISPOSABLE HEAD

- 1 Threaded shaft
- 2 Head calibration value
- 3 Blade housing
- 4 Slide
- 5 Blade
- 6 Aplanating plate

## B. SUCTION RING

- 1 Vacuum port
- 2 Threaded shaft
- 3 Guiding
- 4 Stopper with engraved stop value
- 5 Handle

## C. SETTING THE STOPPER

- 1 Spanner
- 2 Stopper

## D. CONNECTING THE SUCTION RING TO THE ASPIRATION TUBING

- 1 Aspiration tubing end
- 2 Tubing of the suction ring

## E. MOTOR

- 1 Motor housing
- 2 Motor cord
- 3 Threaded ring
- 4 Oscillation motor shaft
- 5 Head alignment guide
- 6 Advance drive

## F. MOUNTING THE HEAD ON THE MOTOR

- 1 Open the pack containing the disposable head and suction ring
- 2 Screw the head clockwise onto the threaded motor shaft.
- 3 Check that the head is screwed on tightly to the motor. There should be no rotation of the head.

## G. MOUNTING THE HEAD ON THE SUCTION RING

- 1 When inserting the head in the rails of the suction ring, the handpiece must be at a 45° angle from the threaded shaft of the ring
- 2 Once the head properly inserted, slightly reduce the angle from 45° to 0°
- 3 Carefully engage the head into the ring

## H. ALIGNMENT OF THE HEAD, MOTOR AND SUCTION RING

- 1 Proper assembly
- 2 Incorrect alignment
- 3 Alignment of the threaded shaft with the motor

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