EN Instruction for use of the OU LC ring (19186) with the Single-Use artificial chamber (19182) and the infusion set (19181)

IMPORTANT: Check the integrity of the pouches, as well as the expiration date on the labels. The medical device is sterile and for single use only. THE DESIGN OF THESE MEDICAL DEVICES DO NOT ALLOW THEIR REUSE. They should not be resterilized or reused.

To select the appropriate combination of disposable OU LC heads, please refer to the document entitled #65090*.


This document is available on electronic versions on Moria website: www.moria-surgical.com (restricted area - «Library» heading). An individual version is also available upon request (Moria's addresses indicated on back page).

Contents

- 19182** Artificial chamber base with 2 ports (one of which with a one-way stopcock), chamber cover, and white screw
- 19186** OU LC guide ring and blue screw
- 19181** Infusion tubing set consisting of one IV line with drip chamber, extension tubing, and one short line with one-way stopcock and syringe

Short line

IV line

Remove all the sterile parts from the pouches.
Place the base of the chamber on the table and make sure it is stable.

Preparing the lines to obtain 90 mmHg pressure in the artificial chamber:

Attach the short line to the chamber (to the port without stopcock).
Check that the 2 stopcocks are open.


Fill the syringe with the corneal storage medium or physiologic solution (BSS or similar).

Attach the syringe to the short line. Place a finger on the top of the piston and fill the line completely (from the short line to the end of the chamber tubing). Close the chamber stopcock.

* Instructions for use of the One Use Large-Cut (OU LC) heads for artificial chamber.

Height between the liquid in the drip chamber and the top of the artificial chamber: around 1.20 m

Raise the IV bottle to a height of 1.20 m between the drip chamber and the top of the artificial chamber (cornea) to ensure obtaining pressure of 90mmHg.


Attach the IV bottle to the IV line. Open the air valve. Purge the line, and use the roller to close the system.

Check that there is no air bubble in the IV line up to the chamber. Attach the IV line to the chamber stopcock and release the roller of the IV line while keeping the chamber stopcock closed.

Positioning the graft on the artificial chamber:

RECOMMENDATION: Use grafts with a rim having a minimum diameter of 16mm and a well-centered cornea.

If desired, apply cohesive viscoelastic to the endothelium.

Use the syringe to gently infuse the system with corneal storage medium or BSS. Position the donor graft on the top of the piston, endothelium downwards. Slightly lift the graft with forceps to remove the possible air bubbles.

Precisely center the cornea on the stage and place the cover carefully in position (inserting the pin into the slot). Check that the graft is well centered. Seal the chamber by placing the white locking ring over the cover.

Open the chamber stopcock at least three seconds. Then, close the chamber stopcock to maintain pressure of 90 mmHg. **Check the 2 stopcocks are closed.**

If desired, check the pressure with a tonometer.

Donor graft cut:

Main steps:

- remove the epithelium and dry the corneal surface
- measure corneal thickness
- if desired, check the pressure with a tonometer
- place the guide ring for OU LC head on the cover and screw it on with the blue screw
- make the cut(s) with the ONE turbine and the appropriate OU LC heads (refer to document #65090).

Removal of the graft and disposal of the devices:

Once the cuts are completed, unscrew the guide ring.

To prevent a possible collapse of the graft, open the chamber stopcock, unscrew the white screw, and gently lift the cover to grasp the graft.

The graft is now ready to be placed on a Moria punch to obtain the desired diameter.

At the end of the surgical procedure, discard the single-use artificial chamber and its accessories, the OU LC head, and the single-use tubing into an appropriate container for disposal (according to the regulations in force).

CAUTION for USA only: US federal laws restrict this device to sale by, or on the order of, a physician.

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MORIA Inc