EN USE, INSTRUMENT CARE, CLEANING AND STERILISATION OF REUSABLE AND AUTOCLAVABLE INSTRUMENTS
CE 2017

CAUTION

⚠ Moria cannot be held responsible for direct or indirect damage suffered by the customer in the event of inappropriate use, maintenance, cleaning or sterilisation of the medical device.

⚠ Moria disclaims all liability in the event of incidents linked to intervention on the instrument by a person or entity not authorised by Moria.

The Moria guarantee does not apply in the event of dismantling, modification or intervention performed on the instrument (e.g. marking) by a person not authorised by Moria (if the instrument needs repair, it must be decontaminated before shipping, otherwise Moria reserves the right to return the instrument in the condition in which it was received).

Products concerned:

Reusable Moria medical devices carrying the references 20105 to 20113 are compatible with the autoclave (steam steriliser 134°C).

Important:

Devices with the references 20106, 20107, 20109 and 20111 must be used in combination with silicone tubes (Ref. 20120).

Incompatibilities:

- Do not use bleach.
- Do not use soda for instruments with the references 20105, 20109 and 20110.

I. STORAGE CONDITIONS

- Store the instruments in a clean, dry area (temperature between 15 and 35°C and relative humidity between 30 and 75%), to prevent the risk of condensation in the packaging.
- Keep the protective devices on the instruments during storage.
- Do not store the instruments close to products or in an environment that could have a corrosive or magnetic action.
- Prevent any contact between instruments particularly if they are of different materials.

II. INSPECTION BEFORE ANY USE

- As the instruments are not delivered sterile, do not use them without cleaning and sterilising beforehand (see § IV. Cleaning).
- Ensure the selected instrument is appropriate and functional and that it will not interfere with other medical devices during the medical procedure.
- Ensure that the instrument has not suffered impacts and that it is free of marks or signs of wear, particularly for the parts in contact with the patient or medical personnel.
- If in doubt, consult the Moria website or your Moria agent to obtain more information on the indication and restrictions for use of the medical device.

III. PRIOR INSPECTION BEFORE SURGERY

- For distally controlled instruments with references 20105, 20109 and 20110, gently remove the protection keeping the forceps closed.
- Before using the following instruments:
 - IRIS® retinal implant forceps, angled (Ref. 20106)
 - IRIS® retinal implant forceps, straight (Ref. 20107)
 - IRIS® retinal implant vertical forceps (Ref. 20109)
 - IRIS® retinal implant holder (Ref. 20111)
- **Put the silicone tubes (Ref. 20120) on the active parts of the instruments as indicated in the instructions for use #65117.**

IV. CLEANING

IV.1 Preliminary precautions for all instruments

- After use of the following instruments:
 - IRIS® retinal implant forceps, angled (Ref. 20106)
 - IRIS® retinal implant forceps, straight (Ref. 20107)
 - IRIS® retinal implant vertical forceps (Ref. 20109)
 - IRIS® retinal implant holder (Ref. 20111)
- **Remove the silicone tubes (Ref. 20120) from the active parts of the instruments as indicated in the instructions for use #65117.**
- Process articulated instruments in the open position to promote even processing.
- Ensure that the instruments are correctly positioned without overload and that the active parts do not bang together during the process.

IV.2 Protocol for cleaning reusable instruments

Refer to the cleaning protocol IIP001 available on the Moria website (www.moria-surgical.com) – restricted access, 'Library/Documentation' heading) which recommends double cleaning.

Decontamination:

- Check that the silicone tubes on the active parts of the instruments have been removed.
- Immerse immediately after leaving the sterile field in an aldehyde-free detergent-disinfectant (e.g. alkazyme 0.5%) for a minimum of 15 minutes at 20°C in accordance with the manufacturer's recommendations.

1st manual cleaning:

- Mechanical cleaning for a minimum of 2 minutes using soft brushes (metal brushes are prohibited) or using an ultrasound machine (provided that the machine does not recycle the cleaning products).
- Rinse with water.
- Immerse in an aldehyde-free detergent-disinfectant (**different bath from the previous one**).

2nd manual cleaning:

- Proceed as for the 1st manual cleaning.

Final rinse:

- Rinse with distilled, filtered or purified water.

Final drying:

- Wipe the device with a lint-free, non-woven, single use soft cloth.

IV.3 Protocol for cleaning distally controlled forceps
(Refs. 20105, 20109, 20110)

Decontamination:

- Check that the silicone tubes on the active parts of the instruments have been removed.
- Immerse immediately after leaving the sterile field in an aldehyde-free detergent-disinfectant (e.g. alkazyme 0.5%) for a minimum of 15 minutes at 20°C in accordance with the manufacturer's recommendations.

1st manual cleaning:

- Perform mechanical cleaning for a minimum of 2 minutes using soft brushes (metal brushes are prohibited), making sure that you do not damage the active parts.
- Rinse with water.
- Unscrew the active part of the instrument body.

- Screw on the cleaning connector (supplied with the instrument) as indicated in the photos below.

- Fill the syringe with alcohol.
- Connect it to the Luer cone.
- Check that the connection is secure and that impermeability is provided.
- Inject the alcohol.

2nd manual cleaning:

- Proceed as for the 1st manual cleaning without putting the instrument back together.

Final rinse:

- Rinse the components with distilled, filtered or purified water.

Final drying:

- Internal: Dry the inside of the tube with filtered compressed air.
- External: Dry with filtered compressed air and wipe the device with a lint-free, non-woven single use soft cloth.

After drying:

- Rescrew the active part onto the handle of the corresponding instrument. The last 3 figures of the reference engraved on the handle and those on the cone of the active part must correspond. To guarantee good functionality it is prohibited to interchange instrument handles.
- Immediately put back the Teflon® and silicone protection whilst holding the forceps closed.

Important: for effective sterilisation, check that the instruments are perfectly clean and free from all residues.

V. STERILISATION OF THE MEDICAL DEVICE

- Sterilise with an autoclave (steam steriliser) at a minimum of 134°C for 18 minutes. Other cycles and temperatures may be used (to be approved by the user). However, autoclave sterilisation temperatures should not exceed 140°C.
- Adhere to the autoclave drying cycle to prevent oxidation.

Inactivation of UTAs (only for patients at risk of Creutzfeld-Jakob's disease – CJD questionnaire/see French Ministry guide DGS/R13/2011/449 – December 2011).

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