



Care and Maintenance for Reusable Microkeratome Accessories

In accordance with ISO 17664-1 and ISO 17664-2

CLEANING, DISINFECTION, STERILIZATION AND MAINTENANCE

This general reprocessing guide applies to MORIA Reusable Microkeratome Accessories used in ophthalmic surgeries for refractive or keratoplasty procedures. Specific, validated reprocessing instructions are enclosed; they include detailed warnings, applicable accessories, and any special limitations or additional instructions.

Deviations from this reprocessing guide are the sole responsibility of the user/facility. This Guide is written in accordance with ISO 17664-1 and ISO 17664-2.

The Reusable Accessories of the MORIA Microkeratome system is intended for performing lamellar resections on the patient's eye or on donor corneal tissue, and consists of:

- Reusable Suction Ring
- Reusable Keratome Metallic Head
- Handpiece (Electrical or Pneumatic), Hose for Pneumatic Handpiece
- Artificial Chamber (AC)
- Guide Ring for Artificial Chamber (AC)
- ALTK Applanating Lenses

MORIA MICROKERATOME	Suction Ring	Keratome Head	Electrical Handpiece	Pneumatic Handpiece & Flexible Hose	Artificial Chamber (AC)	Guide Ring for AC	ALTK Applanating Lenses
OUP SYSTEM	19391/xx ¹ 19519/-1 19391/xx/OV ²	-	19345	-	-	-	-
M2/M2SU SYSTEM	19325/xx ³ 19379/-1	19327/110 19327/130	19326	-	-	-	-
ONE USE ALTK SYSTEM	-	-	-	19155 and 19353	19161 and 19162	19168	-
OUP ALTK SYSTEM	-	-	19175	-		19173	-
CBm ALTK SYTEM	19309/xx ⁴	19170/xx ⁵	-	19303 and 19353		19171 or 19172	19310, 19311 and 19165, 19166

¹ xx = -1, 0, 1, 2, 3

² xx = 1, 2, 3

³ xx = -1, 0, 1, 2, 3

⁴ xx = -1, 0, 1, 2

⁵ xx = 110, 130, 200, 250, 300, 350, 400

Table of contents

1. GENERAL INFORMATION	3
2. CLEANING, DISINFECTION, STERILIZATION AND MAINTENANCE	3
2.1 Important Information/Recommendations for Use.....	3
2.2 Initial Treatment at the Point of Use	4
2.3 Containment and Transportation	4
2.4 Cleaning	4
2.5 Inspection and Maintenance	4
2.6 Packaging and Storage Between Cleaning and Disinfection or Sterilization	4
2.7 Disinfection, Sterilization	4
2.8 Storage After Disinfection, Sterilization	5
2.9 Limitations Reprocessing.....	5
3. WARRANTY	5
3.1 Exclusions.....	5
3.2 Warranty Period.....	5
4. SAFETY	6
4.1 Symbols used in this Guide.....	6
4.2 Symbols Used on the Device	6
5. CONTACT	7
6. INCIDENTS / ACCIDENTS.....	7
7. APPENDICES	7
7.1 Appendix 1 - Compatibility of Reusable Microkeratome Accessories with Reprocessing Steps	8
7.2 Appendix 2 - Manual Cleaning Instructions for Immersible Devices	9
7.3 Appendix 3 - Wiping Cleaning Instructions for Non-Immersible Devices	10
7.4 Appendix 4 - Lubrication Instructions prior to Sterilization	11
7.5 Appendix 5 - Packaging Instructions for Sterilization and Storage Box	12
7.6 Appendix 6 - Steam Sterilization Instructions for Autoclavable Devices	13
7.7 Appendix 7 - Wiping Disinfection Instructions for Non-Autoclavable Devices	14
7.8 Appendix 8 - VH ₂ O ₂ Sterilization Instructions for Non-Autoclavable Devices	15

1. GENERAL INFORMATION

The Reusable Microkeratome Accessories have been specifically designed to ensure optimal performance and safety, provided that the cleaning, disinfection, sterilization, and maintenance instructions are strictly followed in accordance with ISO 17664-1 and ISO 17664-2.

In the event of a malfunction, abnormal wear, or any doubt regarding the proper condition of the device, it is imperative to discontinue its use and have it inspected without delay by MORIA or an authorized representative.

MORIA strongly recommends performing annual maintenance and inspection of all reusable products to maintain their original performance and ensure continued compliance with the microkeratome's initial specifications.

2. CLEANING, DISINFECTION, STERILIZATION and MAINTENANCE

The following sections outline the reprocessing procedures applicable to MORIA Reusable Microkeratome Accessories.

As not all accessories are compatible with conventional reprocessing methods (e.g., immersion cleaning or autoclave steam sterilization), MORIA has developed and validated dedicated reprocessing protocols for each accessory category, in compliance with the Spaulding classification (non-critical, semi-critical, and critical devices).

For details, refer to **Appendix 1 - Compatibility of Microkeratomes Accessories with Reprocessing Steps**

2.1 Important Information/Recommendations for Use

Reusable devices are not delivered sterile and must be cleaned and sterilized or disinfected before use.



Recommendations for Use:

- **Detergent:**

It is recommended to use cleaning solutions containing low-foaming detergent agents (enzymatic, non-enzymatic, or alkaline) with a pH range between 6.0 and 8.0.
Detergents should be used at the concentration and temperature specified by the detergent manufacturer to ensure optimal cleaning performance.
Rinsing procedures must be strictly followed to ensure the complete removal of all detergent residues. Inadequate rinsing may result in chemical residues that could adversely affect subsequent processing steps or compromise device integrity.
Certain alkaline detergents are specifically formulated to be compatible with medical devices. However, users must consult the detergent manufacturer's instructions, limitations, and warnings to verify material compatibility and to identify any components that may be adversely affected by the cleaning agents.
- **Water:**

The use of hard water should be avoided. Softened tap water may be used for most rinsing steps; however, purified water should be used for the final rinse to prevent mineral deposits.
- **Manual Cleaning Tools:**

Tools required for manual cleaning include surgical scrub brushes, soft low-lint cloths, cotton-tipped applicators, and brushes of various sizes and lengths.
Do not use abrasive cleaning tools (e.g., scouring pads or metal brushes), as these may reduce precision, cause scratches, or damage metal surfaces.
Cleaning tools must be cleaned and inspected after each use. The cloth should be clean, lint-free, and replaced frequently. Brushes must be kept clean; worn brushes and disposable cleaning tools should be discarded.
- **Disassembling before treatment/cleaning**

All components of the device must be disassembled prior to cleaning. The suction ring, metallic head, artificial chamber, guide ring for the AC, handpiece, flexible hose, and ALTK applanation lenses must each be cleaned separately.
- **Warnings:**
 - Do not immerse handpieces (electrical or pneumatic) or the flexible hose in water.
 - Do not steam sterilize electrical handpieces or ALTK applanation lenses in an autoclave.
 - To prevent damage to the connectors and ensure proper electrical handpiece function, never pull on the cables or hold the motor by its cable.
 - Do not use abrasive materials or scrapers to clean any device components; use soft brushes instead.

2.2 Initial Treatment at the Point of Use

Cleaning begins at the point of use to prevent soil and contaminants from drying on the devices.

2.2.1 Initial Treatment for Immersible Devices

This procedure applies only to the suction ring, keratome head, artificial chamber, guide ring for the AC, and ALTK applanation lenses.

For detailed instructions, refer to **Appendix 2 - Manual Cleaning Instructions for Immersible Devices**

2.2.2 Initial Treatment for Non-Immersible Devices

This procedure applies only to the handpiece (electrical or pneumatic) and the flexible hose.

For detailed instructions, refer to **Appendix 3 - Wiping Cleaning Instructions for Non-Immersible Devices**.

2.3 Containment and Transportation

Used devices must be transported to the reprocessing location (Center Sterile Reprocessing) in closed or covered containers to prevent contamination risks.

2.4 Cleaning

2.4.1 Manual Cleaning Instruction for Immersible Devices

This procedure applies only to the suction ring, keratome head, artificial chamber, guide ring for the AC, and ALTK applanation lenses.

For detailed instructions, refer to **Appendix 2 - Manual Cleaning Instructions for Immersible Devices**

2.4.2 Wiping Cleaning for Non-Immersible Devices

This procedure applies only to the handpiece (electrical or pneumatic) and the flexible hose.

For detailed instructions, refer to **Appendix 3 - Wiping Cleaning Instructions for Non-Immersible Devices**.

2.5 Inspection and Maintenance

After cleaning and drying:

- Check the general condition of the devices (e.g., no broken parts, free and functional movement).
- Ensure that there is no oxidation or residue.
- Verify that the device remains identifiable (e.g., markings are visible).
- Make sure the device has not been subjected to impact and is free from marks or signs of wear, particularly on parts that come into contact with the patient or medical staff.

Prior to each steam sterilization cycle, it is recommended to lubricate the pneumatic handpiece. Proper lubrication ensures optimal performance and extends the service life of the device.

For detailed instructions, refer to **Appendix 4 - Lubrication Instructions prior to Sterilization**

2.6 Packaging and Storage Between Cleaning and Disinfection or Sterilization

Reusable devices that are stored between cleaning and sterilization or disinfection must be thoroughly dried to prevent microbial contamination that may result from moisture.

Containment devices may be stacked for storage. Store in a clean, dry, and dust-free environment.

For details on packaging for terminally sterilized medical devices and storage boxes, refer to **Appendix 5 - Packaging Instructions for Sterilization and Storage**.

2.7 Disinfection, Sterilization

2.7.1 Steam Sterilization for Autoclavable Devices

This procedure applies only to the suction rings, keratome heads, artificial chambers, guide rings for the AC, pneumatic handpieces, and flexible hoses.

For detailed instructions, refer to **Appendix 6 - Steam Sterilization Instructions for Autoclavable Devices**.

2.7.2 Disinfection for Non-Autoclavable Devices

This procedure applies only to electrical handpieces and ALTK applanating lenses.

For detailed instructions, refer to **Appendix 7 - Disinfection Instructions for Non-Autoclavable Devices**

2.7.3 VH₂HO₂ Sterilization for Non-Autoclavable Devices

This procedure applies only to electrical handpieces (except for 19326) and ALTK applanating lenses.

For detailed instructions, refer to **Appendix 8 - VH₂O₂ Sterilization Instructions for Non-Autoclavable Devices**.

2.8 Storage After Disinfection, Sterilization

Sterile-packaged devices should be stored in a designated, limited-access area that is well ventilated and provides protection from dust, moisture, insects, vermin, and extremes of temperature and humidity.

The shelf life depends on the type of sterile barrier employed, the storage method, and the environmental and handling conditions. Each health care facility should define a maximum shelf life for sterilized medical devices prior to use.

2.9 Limitations Reprocessing

MORIA does not specify a maximum number of uses for its reusable devices.

The service life of each device depends on several factors, including the method and duration of use, the procedures performed between uses, and the reprocessing treatments applied. The lifespan is also influenced by wear and potential damage during use, the type of sterile barrier employed, storage conditions, environmental factors, handling, and reprocessing parameters.

When MORIA's recommended reprocessing procedures are properly followed, repeated processing has minimal impact on the integrity and performance of the instruments.

A thorough inspection, followed by a functional test prior to use, is the most reliable method to determine whether a device should be withdrawn from service or sent for repair.

Improper, ineffective, or insufficient maintenance may shorten the device's service life and will void the instrument warranty.

3. WARRANTY

The warranty covers any malfunction, material damage, or manufacturing defect in products used in accordance with the instructions provided in this leaflet.

This warranty does not affect the application of statutory guarantees under national legislation governing the sale of consumer goods.

Any equipment to be returned must first be cleaned and disinfected or sterilized, then shipped in its original packaging.

3.1 Exclusions


- Any damage resulting from failure to follow the instructions described in this guide.
- Any damage resulting from transport, improper use or negligence, incorrect handling, system modification, inadequate maintenance, or the use of parts or accessories not supplied or recommended by MORIA in this guide.
- Any defect or malfunction of the system occurring outside the warranty period (as defined in the "Warranty Period" section).
- Normal wear and tear.

3.2 Warranty Period

- The warranty takes effect on the date the equipment is shipped.
- The warranty is valid for a period of 12 months from the effective date.
- In addition to the standard one-year warranty, any repairs performed are guaranteed for three months from the date of invoicing.

4. Safety


















4.1 Symbols used in this Guide

Symbol	Meaning	Note
	Warning	Improper operation may result in serious injury ¹ or death to the user, patient
	Caution	Improper operation may result in bodily injury ² or property damage

¹Serious injury means an injury that results in vision loss, causes subsequent complications, or requires hospitalization or long-term outpatient treatment.

²Bodily injury means an injury, burn, or similar condition that does not require hospitalization or long-term outpatient treatment.

4.2 Symbols Used on the Device

SYMBOL	MEANING
	Catalogue number
	Batch code
	Serial number
	Quantity
	Manufacturer
	Date of manufacture
	Do not use if package is damaged and consult instructions for use
	Keep dry
	Keep away from sunlight
	Consult instruction for use or consult electronic instructions for use
	Caution
	Discard in a separate collection bin
	Patient information web site
	Medical device
	CAUTION for USA only: US federal laws restrict this device to sale by, or on the order of, a physician.
	Unique device identifier
	Product compliant with MDD 93/42/CEE or MDR 2017/745

5. Contact

If you encounter any technical problems with our products, please contact MORIA headquarters or your local distributor.

The latest versions of the user manuals and additional information about our systems are available on the MORIA website:

<http://www.moria-surgical.com>

6. INCIDENTS / ACCIDENTS



IMPORTANT WARNINGS

If any serious incident occurs in relation to the device or its accessories, please contact:

MORIA: correspondant-materiovigilance@moria-int.com

or your Local Distributor

or the Competent Authority of the Member State in which you are located.

7. APPENDICES

COMPATIBILITY OF MORIA REUSABLE MICROKERATOME ACCESSORIES WITH REPROCESSING STEPS



Certain MORIA Reusable Microkeratome Accessories may be incompatible with conventional reprocessing methods, such as immersion cleaning or steam sterilization in an autoclave.

MORIA has developed and validated appropriate reprocessing methods for each accessory category, in accordance with the Spaulding classification (Non-Critical, Semi-Critical, and Critical Devices).

7.1.1 General Compatibility

- The Suction Ring, Keratome Head, Artificial Chamber, and Guide Ring are suitable for immersion cleaning.
- The Suction Ring, Keratome Head, Artificial Chamber, Guide Ring, and Pneumatic Handpieces are suitable for steam sterilization in an autoclave.
- Electrical and Pneumatic Handpieces must not be immersed.
- Electrical Handpieces and ALTK Applanating Lenses must not be autoclaved.

7.1.2 Validated Reprocessing Steps by Device Category

Device Type	Validated Reprocessing Steps
Immersible Devices	Manual Cleaning in an Ultrasonic Bath
Non-Immersible Devices	Wiping Cleaning
Autoclavable Devices	Steam Sterilization
Non-Autoclavable Devices	Intermediate-Level Wiping Disinfection
Non-Autoclavable Devices	VH ₂ O ₂ (Vaporized Hydrogen Peroxide) Sterilization

7.1.3 Table of Recommended Reprocessing Steps

MORIA Accessories	Manual Cleaning	Wiping Cleaning	Wiping Disinfection	Steam Sterilization	VH202 Sterilization
Suction Ring	YES	-	-	YES	-
Keratome Head	YES	-	-	YES	-
Electrical Handpiece	NO	YES	YES	NO	YES ⁶
Pneumatic Turbine and Flexible Hose	NO	YES	-	YES	-
Artificial Chamber (Ac)	YES	-	-	YES	-
Guide Ring for AC	YES	-	-	YES	-
ALTK Applanating Lenses	YES	-	YES	NO	YES

⁶ Only for the V-PRO system from STERIS and for handpieces 19345 and 19175

Legend:

YES = Recommended and validated

NO = Not compatible

- = Not applicable

Notes

- All reprocessing procedures must follow validated parameters defined by MORIA in accordance with ISO 17664-1 and ISO 17664-2.
- Only use cleaning and disinfecting agents recommended for ophthalmic surgical instruments.
- Ensure complete drying before assembly or sterilization.
- Regular inspection is required to verify absence of wear, corrosion, or damage.

MANUAL CLEANING FOR IMMERSIBLE DEVICES



APPLICABLE ONLY TO SUCTION RINGS, KERATOME HEADS, ARTIFICIAL CHAMBERS, GUIDE RINGS FOR AC, AND ALTK APPLANATING LENSES.

Accessories	MORIA References
Suction Ring	19391/xx, 19519/-1, 19391/xx/OV, 19325/xx, 19379/-1, 19309/xx
Keratome Head	19327/xx, 19170/xx
Artificial Chamber	19161, 19162
Guide Ring for AC	19168, 19173, 19171, 19172
ALTK Applanating Lenses	19310, 19311, 19165, 19166

7.2.1 Initial Treatment for Immersible Devices

Keep the devices moistened after use and before cleaning.

1. Disassemble the devices
2. Immerse the devices in the cleaning solution¹ for at least 20 min within 30 min after surgery.

7.2.2 Cleaning for Immersible Devices

Pre-cleaning

1. Rinse with cold tap water (<30°C) for at least 30 seconds².
2. Immerse the device in the cleaning solution¹ (30-40°C) for at least 15 min to moisten and loosen soil.
3. Scrub all device surfaces with the cleaning solution using a soft cleaning brush, taking care not to damage the device, until all visible soil has been removed.
4. For devices with lumen: flush the cleaning solution through the lumen³ five times, then flush with cold tap water (<30°C)³ five additional times.
5. Rinse with cold tap water² (<30°C) for at least 30 seconds.

Manual Cleaning

1. Immerse the device in the cleaning solution¹ (30-40°C) and sonicate for at least 10 min in an ultrasonic bath.
2. Rinse with cold tap water² (<30°C) for at least 30 seconds.
3. For devices with lumen: flush demineralized water (30-40°C) through the lumen³ five times.
4. Immerse in demineralized water (30-40°C) for at least 5min.
5. For devices with lumen: flush 70° alcohol through the lumen³ five times.
6. Wipe and dry with filtered air.



¹Refer to manufacturer's detergent recommendation of use regarding concentration, T°.

²Expose all surfaces and cavities to flowing water.

³Use a 5mL syringe to flush through the lumen.

Manual Cleaning Validation Study Information:

- Detergent Neodisher Mediclean Forte at 0.5% solutions (Dr Weigert, Hamburg, Germany) were used during the cleaning process validation.
- Trial Report reference: CLE_002

WIPING CLEANING FOR NON-IMMERSIBLE DEVICES



APPLICABLE ONLY TO ELECTRICAL HANDPIECES, PNEUMATIC HANDPIECES, AND FLEXIBLE HOSES.

Accessories	MORIA References
Electrical handpiece	19175, 19345, 19326
Pneumatic handpiece	19155, 19303
Flexible hose	19353

7.3.1 Initial Treatment for Not Immersible Devices

1. Wipe the entire device⁴ with a lint-free surgical wipe moistened with the cleaning solution^{1,2}.
2. Wipe the entire device⁴ with a lint-free surgical wipe moistened with distilled water².

7.3.2 Cleaning for Not Immersible Devices

1. Spray the cleaning solution¹ on the entire device⁴ until all parts are thoroughly wetted on moving parts.
2. Wipe the entire device⁴ with a lint-free surgical wipe moistened with the cleaning solution^{1,2} until soil and organic matter have been visibly removed.
3. For electrical handpieces: rinse the screw core of the motor with a brush⁵ moistened with distilled water (see Figure 1).
4. Wipe the entire device⁴ with a lint-free surgical wipe moistened with distilled water²
5. Dry carefully with a lint-free surgical wipe³ then with clean, filtered, compressed air.




Figure 1: Rinsing of the screw core with a brush⁵



¹Refer to manufacturer's detergent recommendation of use.

²Replace the wipe when it becomes dry.

³Replace the wipe when it becomes too wet.

⁴For handpieces: including cable and oscillation shaft. The metallic parts must be strongly wiped.

⁵Use MORIA disposable brush 19149X25

Wiping Cleaning Validation Study Information:

- CaviCide spray and CaviWipes wipes (cleaning and disinfection products from Metrex) were used during the wiping cleaning validation.
- Trial Report reference: CLE_003

LUBRICATION PRIOR TO STERILIZATION



APPLICABLE ONLY TO PNEUMATIC HANDPIECES.

Accessories	MORIA References
Pneumatic Handpiece	19155, 19303

Pneumatic handpieces must be lubricated before each steam sterilization cycle using an appropriate lubrication spray and nozzle^{1,2}.

1. Insert the tip nozzle into the smaller of the two large tubes (see Figure 2), which is the air intake hole, and always hold the spray can in an upright position.
2. Cover the larger of the two air tubes (the air exhaust hole) and the turbine head with a clean tissue (see Figure 3), then spray for 1 second³.
3. Keep the pneumatic handpiece in an upside-down position for 2 min (see Figure 4) prior to steam sterilization to allow excess lubricant to drain onto a lint-free cloth.




Figure 2



Figure 3



Figure 4



¹The use of KaVo lubrication spray is recommended.

²Refer to the manufacturer's instructions for proper spray application..

³A single burst of KaVo spray is sufficient to lubricate the turbine while preventing excess lubricant inside the mechanism.

PACKAGING FOR STERILIZATION AND STORAGE

Accessories	Sterilization Tray MORIA References	Storage Box MORIA References
Suction Ring	22519513 ¹ , 255193301, 2519343 ¹	-
Keratome Head	255193301, 2519343 ¹	-
Artificial Chamber	22519167 ¹	-
Guide Ring for AC	22519167 ¹	-
Pneumatic Handpiece and Hose	22519139 ¹	-
Electrical Handpiece	STERIS tray ³ only for 19345, 19175	225195142, 225191742, 225191742
ALTK Applanating Lenses	STERIS tray ³	225191692

7.5.1 General Information



The packaging for terminally sterilized medical devices should meet the following requirements:

- be suitable for steam sterilization and comply with EN ISO 11607-1
- provide adequate protection for both the devices and the sterilization packaging against mechanical damage.

Packaging of devices in tray:

- The set of devices shall be placed in dedicated MORIA sterilization trays or cases to ensure adequate protection of the instruments. Trays or cases with lids shall be wrapped in a standard medical-grade steam sterilization wrap, using the AAMI double-wrap method or an equivalent validated process.
- A rigid container system may also be used for sterilization, provided it is validated for the intended sterilization method.

Packaging of devices individually:

- Medical-grade steam sterilization pouches of appropriate sizes shall be used to double-package individual devices. The inner pouch shall be sufficiently large to contain the instrument without stressing the seals or tearing the material, yet small enough to fit within a secondary pouch without compromising the integrity of the overall package.

7.5.2 List of Dedicated Packaging per System

Accessories	OUP system	M2/M2SU system	ONE USE ALTK SYSTEM	OUP ALTK SYSTEM	CBm ALTK SYTEM
Suction Ring	22519513 ¹	255193301	-	-	22519343 ¹
Keratome Head	-		-	-	
Electrical Handpiece	22519514 ² STERIS tray ³	22519931 ²	-	22519174 ² STERIS tray ³	-
Pneumatic Handpiece and Flexible Hose	-	-	22519139 ¹	-	22519139 ¹
Artificial Chamber (AC)	-	-	22519167 ¹		
Guide Ring for AC	-	-			
ALTK Applanating Lenses	-	-	-	-	22519169 ² STERIS tray ³

¹Steam sterilization tray

²Storage box for non-autoclavable devices - do not use for sterilization.

³VH₂O₂ sterilization tray: VP004 or P0033 from STERIS - do not use an aluminum tray.

STEAM STERILIZATION FOR AUTOCLAVABLE DEVICES

⚠ APPLICABLE ONLY TO SUCTION RINGS, KERATOME HEADS, PNEUMATIC HANDPIECES, FLEXIBLE HOSES, ARTIFICIAL CHAMBERS, AND GUIDE RINGS FOR AC.

Accessories	MORIA References
Suction Ring	19391/xx, 19519/-1, 19391/xx/OV, 19325/xx, 19379/-1, 19309/xx
Keratome Head	19327/xx, 19170/xx
Pneumatic Handpiece	19155, 19303
Flexible Hose	19353
Artificial Chamber	19161, 19162
Guide Ring for AC	19168, 19173, 19171, 19172

- Steam autoclave sterilization using a pre-vacuum (forced air removal) cycle is recommended by MORIA for reusable devices. The minimum sterilization parameters are provided in the table below.
- The sterilizer manufacturer's recommendations shall always be followed. When sterilizing multiple devices set in a single cycle, ensure that the maximum load specified by the manufacturer is not exceeded.
- Device sets shall be properly prepared and packaged in trays and/or cases that allow adequate steam penetration and direct contact with all device surfaces.
- The devices and/or device trays shall undergo a complete sterilization drying cycle, as residual moisture from autoclaves may lead to staining, discoloration, or corrosion.

Steam Sterilization Cycle	Preparation	Exposure Time	Temperature ⁵	Minimum Drying Time, in Room Chamber ⁶
Pre-vacuum ^{1 2}	Double-wrapped	18 min	134°C / 273°F	30 min
Pre-vacuum (UK) ³	Double-wrapped	3 min	134°C / 273°F	30 min
Pre-vacuum (USA) ⁴	Double-wrapped	4 min	132°C / 270°F	30 min



- Gravity displacement sterilization cycles are not recommended, as the required cycle times are excessively long and not suitable for devices with lumens.
- Flash (immediate-use) steam sterilization in unwrapped packaging and without the recommended drying time, with exposure at 132-134 °C (270-273 °F) for the times listed in the table, should be used only in emergency situations. Devices must be thoroughly cleaned and disassembled prior to sterilization.

¹Recommended by the World Health Organization (WHO) and the French Ministry of Health (Guide DGS/RI3/2011/449 - December 2011) regarding the risk of TSE/CJD contamination.

²The Robert Koch Institute (RKI) recommends 134 °C (273 °F) for 5 minutes, extended to 18 minutes if an alkaline detergent was not used during cleaning, due to TSE/CJD contamination risks.

³Local or national specifications shall be followed where steam sterilization requirements are stricter or more conservative than those listed in the table.

⁴United States only: per ANSI/AAMI ST79.

⁵The sterilization temperature shall not exceed 140 °C (284 °F).

⁶Use the autoclave drying cycle to prevent oxidation or corrosion of devices.

Sterilization Validation Study Information:

- Devices packaged in dedicated MORIA trays, and wrapped using the double-wrapping method according to AAMI standards
- Trial Report reference: STE_002

WIPING DISINFECTION FOR NON-AUTOCLAVABLE DEVICES



APPLICABLE ONLY TO ELECTRICAL HANDPIECES AND ALTK APPLANATOR LENSES.

Accessories	MORIA References
Electrical Handpiece	19175, 19345, 19326
ALTK Applanating Lenses	19310, 19311, 19165, 19166

1. Spray disinfection solution¹ on the entire device⁴ until all parts are thoroughly wetted. Allow the solution to act for 3 min (see Figure 5 for handpieces).
2. Wipe the entire device⁴ with a lint-free surgical wipe moistened with the disinfection solution^{1,2} for at least 1 min. Repeat the process using a new moistened wipe (see Figure 6 for handpieces).
3. Wrap the device with a new wipe moistened with disinfection solution¹ and allow it to act for at least 10 min.
4. Wipe the entire device⁴ with a lint-free surgical wipe moistened with distilled water². Repeat the process using a new moistened wipe.
5. Dry carefully with a lint-free surgical wipe³, then with clean, filtered, compressed air.



Figure 5: Application of cleaning and disinfection spray



Figure 6: Wiping with the cleaning and disinfection product, with emphasis on the metallic parts



¹Refer to manufacturer's disinfectant recommendation of use.

²Replace the wipe when it becomes dry.

³Replace the wipe when it becomes wet.

⁴For handpieces: including cable and oscillation shaft. The metallic parts must be strongly wiped

Wiping Disinfection Validation Study Information:

- CaviCide spray and CaviWipes wipes (cleaning and disinfection products from Metrex) were used during the wiping disinfection validation.
- Trial Report reference: DIS_001

VH₂O₂ STERILIZATION FOR NON-AUTOCLAVABLE DEVICES



APPLICABLE ONLY TO ELECTRICAL HANDPIECES AND ALTK APPLANATOR LENSES.

Accessories	MORIA References
Electrical Handpiece ¹	19175, 19345
ALTK Applanating Lenses	19310, 19311, 19165, 19166

The listed Electrical Handpieces¹ and ALTK Applanator Lenses are compatible with the low-temperature hydrogen peroxide (VH₂O₂) sterilization² using the V-PRO System from STERIS³



¹Excluding 19326

²Do not use trays made of aluminum or the dedicated handpiece storage box. Aluminum is not compatible with VH₂O₂ sterilization. The dedicated storage boxes are not perforated.

³Sterilizers V-PRO max, V-PRO max2, V-PRO 60, V-PRO s2 - compatible with the Lumen, Non-Lumen and Flexible cycles

Hydrogen Peroxide Sterilization Validation Study Information:

- Devices packaged in STERIS trays VP004 and VP0033, and wrapped using the double-wrapping method according to AAMI standards
- Trial Report reference: STE_003



MORIA SA
27, rue du Pied de Fourche
03160 Bourbon L'Archambault
FRANCE
Phone: +33 (0)1 46 74 46 74
Fax: +33 (0)1 46 74 46 00
moria@moria-int.com
www.moria-surgical.com

MORIA Inc
1050 Cross Keys Drive
Daylestown, PA 18902
USA
Phone: (800) 441 1314
Fax: +1 (215) 230 7670
orders@moriausa.com
www.moria-surgical.com

Moria Surgical Spain
C. de Pablo Iglesias, 58,
08908 L'Hospitalet de Llobregat
Barcelona - SPAIN
Phone: +34932896270
info@moria-spain.com
www.offtaltech.com

Moria Japan K.K.
Arcadia Building 6F
1-12-3 Kanda
Sudacho Chiyoda-Ku
Tokyo 101-0041
JAPAN
Phone: 81-3-6260-8309
Fax: 81-3-6260-8310
moria@moriajapan.com
www.moriajapan.com

ALCHIMIA S.R.L.
(a Moria company)
Viale Austria, 14
35020 - Ponte San Nicolò (PD)
ITALY
Phone: +39 049 89.62.074
Fax: +39 049 89.62.071
info@alchimiasrl.com
www.alchimiasrl.com

**Moria COMMERCIAL
(CHINA) CO., LTD.**
上海利亚贸易有限公司
Room H 6F Kai Building NO.432
West Huai Hai Road, Changning
district, Shanghai, 200052, P.R.C
中国上海市长宁区淮海西路432号
凯利大厦6楼08室
Phone/Fax: +86 021 52586095
moria@moria-int.com

Moria Limited
Suite 1904, 19th Floor,
ST. George's Building,
No. 2 Ice House Street,
Central - HONG KONG
moria@moria-int.com
www.moria-surgical.com