## EN

#65034 - INSTRUCTION FOR USE  
INSTRUMENT CARE, CLEANING AND STERILISATION OF REUSABLE AND AUTOCLAVABLE INSTRUMENTS according ISO 17664-1

These instructions are recommended for the care, cleaning, maintenance and sterilization of reusable MORIA surgical instruments. This document is intended to assist health care personnel in safe handling practices, effective reprocessing, inspection and maintenance of MORIA reusable devices.

## 1. PRE-REQUISITES BEFORE USE AND DIRECTIONS FOR USE

**Warning:** Reusable Instruments are not delivered sterile and must be cleaned and sterilized before first use.

From a functional perspective, a visual inspection must be carried out prior to use to check for any burrs or debris that could damage tissue or personal protective equipment. Furthermore, the integrity of cutting tools must be verified. Any device deemed to be dull or non-functional in any manner should be returned to MORIA for maintenance or exchange.

- Preliminary precautions:
  - Process the instruments in open articulated position, without undressing at the rear, to promote through even processing.
  - Make sure that the instruments are positioned correctly without overloading and that the active parts do not touch each other during the process.
  - Regularly check for the absence of organic or metallic debris/deposits under a microscope.
  - If necessary, ensure that the MORIA protective tip is properly positioned and it accompanies the instrument during all stages of cleaning / disinfection / sterilization.
  - The use of metallic brushes, scrub pads and other articles likely to damage the instruments must be avoided.

## 2. POINT OF USE

Point of use processing aims to make subsequent cleaning easier. All reusable devices must undergo immediate point-of-use processing.

The use of soft-bristled brushes is authorized to clean the parts from all soils that can potentially alter the action of detergents and decontaminants.

## 3. CLEANING AND DRYING

See table 'Cleaning and drying'.

## 4. INSPECTION AND MAINTENANCE

After cleaning, and drying, please follow the instructions hereunder:

- Check the general aspect of the active parts (e.g no twisted or broken parts, free movement) and check that there is no oxidation of residues.
- Make sure that the instrument has not been subjected to impact and that it is free of marks or signs of wear, primarily the parts in contact with the patient or with the medical staff.
- Check the readability of instrument identification.
- Refer to 65034-Appendix for complete inspection instructions of instrument (available on MORIA Web Site: <http://www.moria-surgical.com>)

## 5. PACKAGING

The packaging for terminally sterilized medical devices should fulfill the following requirements:

- Suitable for steam sterilization and EN ISO 11607-1 and 11607-2
- Sufficient protection of the instruments as well as of the sterilization packaging to mechanical damage

For the packaging instruments in :

- Set of instruments should be placed in generic or dedicated sterilization tray/cases to protect the instruments. Trays and cases with lids should be wrapped in standard medical grade, steam sterilization wrap using the AAMI ST79 (2017) double wrap or equivalent method.

Container cases could also be used for sterilization.

For the packaging of instruments individually:

- Medical grade steam sterilization pouches of the appropriate sizes should be used to double package instrument. Ensure that the inner pouch is large enough to contain the instrument without stressing the seals or tearing the packaging but small enough to be placed in a secondary pouch without compromising the integrity of the total package.

The use of protection pouches made of soft silicone tubing of an appropriate size and thickness is recommended in order to protect the active parts.

## 6. STERILIZATION

Steam autoclave (moist heat) sterilization using a pre-vacuum (forced air removal) cycle is recommended per MORIA for reusable devices, according to minimum parameters sterilization in table below.

Sterilizer manufacturer recommendations should always be followed. When sterilizing multiple instrument sets in one sterilization cycle, ensure that the manufacturer's maximum load is not exceeded.

Instrument sets should be properly prepared and packaged in trays and/or cases that will allow steam to penetrate and make direct contact with all surfaces.

The instrument and/or instrument tray should be processed through a complete sterilization drying cycle as residual moisture from autoclaves can promote staining, discoloration, and rust.

## 12 (EN) SYMBOLS (FR) SYMBOLES (ES) SÍMBOLOS (IT) SIMBOLI (PT) SÍMBOLOS (SV) SYMBOLER (SK) SYMBOLY (SL) SIMBOLI (HU) SZIMBÓLUMOK (BG) СИМВОЛИ (CZ) SYMBOLY (LT) SIMBOLIA

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## 3. CLEANING and DRYING

	Automated Cleaning Process	Manual Cleaning process
<b>Preparation at the Point Use</b>	Immerse in enzymatic/cleaning solution <sup>1</sup>	20 min within 30 min after surgery
<b>Transportation</b>	Safe storage and transportation in a closed container to the reprocessing area to avoid any damage and contamination to the environment.	
<b>Pre-Cleaning</b>	1. Rinse with <b>tap cold water</b> 2. Immerse in <b>enzymatic/cleaning solution</b> <sup>1</sup> to moisten and loosen soil.	<30°C for 30 s
	3. Scrub all surfaces of instrument with the <b>cleaning solution</b> using a soft cleaning brush, without damaging the tip, until all visible soil has been removed.	30 to 40°C for 15 min
	4. For instrument with lumen: flush the <b>cleaning solution</b> through lumen, then with <b>cold tap water</b> <sup>1</sup>	5 times for cleaning solution 5 times for water, <30°C
	5. Rinse with <b>cold tap water</b> <sup>1</sup>	<30°C for 30 s
<b>Cleaning Type</b>	<b>Automated Cleaning (Washer-disinfector)*</b>	<b>Manual cleaning</b>
<b>Cleaning</b>	1. Pre-cleaning with cold water 2. Cleaning with enzymatic/cleaning solution <sup>1</sup> and sonicate in ultrasonic bath	<30°C for 2 min 30-40°C for 10 min
<b>Rinsing</b>	3. First rinse with hot water for neutralization <sup>4</sup> 4. Second rinse with cold water	40-60°C for 2 min <30°C for 1 min
	5. Final rinse with demineralized water <sup>5</sup>	90°C for 3 min
<b>Heating</b>	6. Thermal Heating Disinfection Cycle	90°C for 5 min
<b>Drying</b>	7. Dry	110-120°C for 20 min

<sup>1</sup> Recommended by the World Health Organization (WHO), and French Ministry guide DGS/R3/2011/449 – December 2011 in regards of TSE/CJD risks contaminations.

<sup>2</sup> Robert Koch Institute (RKI) recommends 134°C / 273°F 5 min and to extend to 18 min in case of alkaline detergent not used in cleaning steps and in regards of TSE/CJD risks contaminations.

<sup>3</sup> Local or national specifications must be followed where steam sterilization requirements are stricter or more conservative than those listed in the table

<sup>4</sup> US only : ANSI/AAMI ST79

<sup>5</sup> Sterilization temperature must not exceed 140°C/284°F

<sup>6</sup> Use the autoclave drying cycle to avoid oxidation.

<sup>7</sup> Gravity displacement sterilization cycles are not recommended because cycle times are too long to be practical and not adapted for lumen instruments.

<sup>8</sup> Flash (immediate-use) steam sterilization in unwrapped packaging and without the recommended dry time, by exposure at 132-134°C / 270-273°F for the exposure times listed in table should only be used as an emergency procedure. Instruments must be cleaned and disassembled.

<sup>9</sup> Refer to manufacturer's enzymatic/cleaning solution recommendation of use regarding concentration, T°.

<sup>10</sup> Expose all surfaces and cavities to flowing water.

<sup>11</sup> Use a 5 mL syringe to flush through lumen.

<sup>12</sup> If an alkaline detergent is used, adequate neutralization in accordance with the manufacturer's enzymatic/cleaning solution instructions is to be performed.

## Cleaning Validation study information:

See Trial Reports references: STE\_001 according to NF EN ISO 17665-1: 2006, AAMI TIR 12 (2020), ANSI/AAMI ST 79 (2017)

7. STORAGE

Store instruments in clean dry conditions in order to prevent risk of condensation inside the packaging.

Keep the protective cover or the MORIA protective tip on the active part of the instrument during storage.

Do not store instruments near products or in an environment that may have a corrosive or magnetic effect.

Avoid contact between instruments especially if they are made of different materials.

## 8. TRANSPORT

Transportation inside healthcare center :

Used devices must be transported to the processing location (i.e. sterile Processing) in closed or covered containers to prevent contamination risks.

Transportation outside the healthcare center:

Used devices must be transported to the original manufacturer (i.e. Repair processing) in closed or covered containers to prevent contamination risks.

Instrument sets should be transported to the original manufacturer (i.e. Repair processing) in closed or covered containers to prevent contamination risks.

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