## EN

65125 - INSTRUCTION FOR USE  
EN PIM UU

## 1. DESCRIPTION / FEATURES

PIM UU instruments are EO (Ethylene Oxide) sterilized and ready for use.

## 2. INTENDED AND INDICATION FOR USE

PIM UU instruments are intended as a single-use instrument used in ophthalmic surgery for anterior segment procedures.

## 3. INTENDED USER PROFILE

PIM UU instruments must be used only by experienced ophthalmic surgeon/qualified physician

## 4. TARGET PATIENT POPULATION

PIM UU instruments are indicated for use without restriction to the target population (age, sex and health status).

## 5. CONTRAINDICATION

No known contraindication

## 6. WARNING / CAUTION / PRECAUTION / ADVERSE EVENTS

a) General  
• Federal (USA) law restricts this device sale by, or on order of, a physician.

• Do not use if expiration date shown on the label has passed.

• Do not use if storage conditions have failed (humidity traces). This device should be stored in a clean, dry location at room temperature and away from direct sunlight.

• Do not use if the packaging or device is open, damaged, or wet (check the integrity of pouch).

• Eliminate contact with instrument tip while not in surgical use as damage to the tip could occur.

• Follow medical operating procedures during surgery.

• Stop the use of the device in case of any malfunction during the procedure, including potential complication described in part 5.

• This product is intended for single-use, one patient only. Reuse or re-sterilization may cause diminished product performance, contamination or infection.

## b) Specific to the medical device

• No specific warning to the medical device.

## 7. PRECAUTION

The design of these medical devices are « single use only ». Reuse of this device may affect its mechanical and biological features, and its clinical performances; and may cause device failures or may expose the patient to adverse events such as contamination, bacterial infections, or allergic reactions.

This instruction for use refers to the use of medical device but does not describe surgical steps to perform ophthalmic procedures

## 8. INTENDED PERFORMANCE AND CLAIMS

• Intended performances of:

- Capsulorhexis forceps: grasping tissue to realize capsulorhexis.

- Eyelid speculums: maintain eyelids open.

- Tissue forceps: easily manipulating tissues, adapted manipulation of corneal tissues,

- Sutures forceps: nylon filament holding, good lens manipulation,

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- Schaar: snijden van weefsels,

- Ringen: eye holding,

- Needle holders: good needle holding,

- Markers: to mark the chosen angle axis, to mark cornea or sclera,

- Hooks for lens: manipulate the crystalline and IOL lens, divide the nucleus.

## 9. PRE-REQUISITES BEFORE USE AND DIRECTIONS FOR USE

a) Checking the medical device  
• Determine package integrity. If package has been compromised, breached wet and/or damaged in any way, do not use this device.

b) Preparing the medical device

• Affect sterile transfer of the product to the sterile field.

• Remove product from protective packaging.

• Verify tip actuation by squeezing the handle and determine if tip is damaged. If tip is damaged, do not use.

• The instrument is now prepared for use.

c) Using the medical device

• Follow established surgical procedures.

d) End of procedure

• This single use surgical instrument should be considered clinical waste and should be disposed in accordance with clinical waste laws applicable in your country.

## 10. RESIDUAL RISK

• To date, the risks identified by the manufacturer associated to the medical device

• There are no adverse reactions associated with the use of this product other than those related to the specific procedure being performed (to be assessed by the physician).

## 11. POSSIBLE ADVERSE EVENT

• Any serious incident occurring in conjunction with this device must be reported to the manufacturer and the competent authority of the user's and/or patient's Member State.

## 12. DISPOSAL

• Discard into an appropriate container for destruction (according to local regulations).

## 13. LIST OF REFERENCES

• All single use instruments. Refer to MORIA instrumentation catalog.

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