- Please unfold the following pages containing the photos quoted in the reference.
- Veuillez d  plier les pages suivantes contenant les photos cit  es en r  f  rence.
- Por favor, despliegue las p  ginas siguientes que contengan las fotos referenciadas.
- Vogliate aprire le seguenti pagine, contenenti le foto in oggetto.
- Falten Sie die folgenden Seiten bitte auf, die die in Verweis zitierten Fotos enthalten.
- Sla volgende pagina's open voor de foto's waar in de tekst naar verwezen wordt.
- Por favor desdobre as p  ginas seguintes que cont  m as fotografias citadas em refer  ncia.
- F  ljande sidor visar bildreferanser.
- L  tten atıfta bulunan resimleri g  rmek i  in sayfaları   evirin.

Moria

MORIA S.A. – 27, rue du Pied de Fourche - 03160 BOURBON L'ARCHAMBAULT – France
#65078-C-05.2023

EN - CBm-ALTK CBSU for artificial chamber	
USER GUIDE	1-16
FR - CBm-ALTK CBSU pour chambre artificielle	
MANUEL UTILISATEUR.....	17-34
ES - CBm-ALTK CBSU para cámara artificial	
MANUAL DE USUARIO	35-50
IT - CBm-ALTK CBSU per camera artificiale	
MANUALE DI ISTRUZIONI	51-66
DE - CBm-ALTK CBSU für künstliche Kammer	
BENUTZERHANDBUCH.....	67-82

NL - CBm-ALTK CBSU voor kunstmatige kamer	
HANDLEIDING.....	83-98
PT - CBm-ALTK CBSU para câmara artificial	
MANUAL DO UTILIZADOR.....	99-116
SV - CBm-ALTK CBSU för artificiell kammare	
ANVÄNDARHANDBOKL.....	117-132
TR - CBm-ALTK Vakum kamarası için CBSU	
KULLANIM KILAVUZU	133-149

Ref.: 65078

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















Document Aletiq : Rév. C-05.2023 (Validée)

CONTENTS

1	WARNINGS	2
1.1	MAINTENANCE.....	2
1.2	USE OF GENERIC PRODUCTS OR REUSE OF SINGLE-USE CONSUMABLES	2
2	LIST OF EQUIPMENT AND ASSOCIATED ACCESSORIES	3
3	REGULATORY INFORMATION	4
4	MANUFACTURER INFORMATION ON LABELS.....	5
5	IMPORTANT INFORMATION	5
5.1	DESCRIPTION	5
5.2	INDICATIONS.....	5
5.3	PRODUCT PERFORMANCE – THICKNESS OF CORNEAL FLAPS	6
5.4	CONTRAINDICATIONS	6
5.5	WARNINGS	6
5.6	POTENTIAL ADVERSE EFFECTS	7
5.7	PRECAUTIONS	7
5.8	PRODUCT COMPLAINT	8
5.9	PACKAGING.....	8
6	INSTALLATION AND CONNECTION	8
7	OPERATION.....	10
7.1	LAMELLAR RESECTION ON A DONOR EYE (GRAFT).....	10
7.2	RESECTION ON THE PATIENT	11
8	IN THE EVENT OF A PROBLEM	11
9	MAINTENANCE	11
9.1	INTRODUCTION	11
9.2	INITIAL TREATMENT AND STORAGE	12
9.3	PROTOCOL FOR CLEANING / DISINFECTING / STERILIZING / STORING THE TURBINE	12
9.4	PROTOCOL FOR CLEANING / DISINFECTING / STERILIZING / STORING THE TURBINE HOSE.....	13
9.5	PROTOCOL FOR CLEANING / DISINFECTING / STERILIZING / STORING CBm-ALTK HEADS AND SUCTION RINGS	13
10	WARRANTY.....	14
10.1	SCOPE OF WARRANTY	14
10.2	WARRANTY EXCLUSIONS.....	14
10.3	WARRANTY PERIOD	14
10.4	LIABILITY.....	14
11	PHOTOS	15

This user guide presents all the information necessary for using and maintaining the CB microkeratome (associated with the ALTK artificial chamber and the single-use artificial chamber), its heads (reusable or single-use) and its accessories (see list below).

The latest version of these user guides and additional information about your microkeratome are available on Moria's restricted-access website: <http://www.moria-surgical.com>. You must register to access it.

1 WARNINGS

1.1 MAINTENANCE

The CBm-ALTK and CBSU system for artificial chambers was designed to ensure optimal functionality, provided that the maintenance recommendations are scrupulously followed. If for any reason the system does not function correctly, please have it checked immediately by Moria. Moria strongly recommends annual maintenance and servicing of all its reusable products in order to maintain the original performance of your microkeratome.

As only Moria and its agents are fully expert in Moria products, servicing and maintenance must be carried out by Moria or its approved agents.

Moria is not liable for malfunction or damage to apparatus, poor results, or surgical complications because maintenance was carried out by an unqualified operator or third party.

Any such unauthorized intervention renders the warranty and any maintenance contract null and void.

1.2 USE OF GENERIC PRODUCTS OR REUSE OF SINGLE-USE CONSUMABLES

The materials for the blade, blade holder and microkeratome head have been selected for their biocompatibility and drag-free characteristics. Blade dimensions and tolerances have been designed to match the dimensions and tolerances of the keratome head. Moria's manufacturing and control procedures guarantee that there is no dimensional incompatibility between the head and the blade, and ensures that the blade slides freely in the head.

Products designed for single-use must not be reused under any circumstances. The reuse of single-use products, or the use of consumables other than those supplied by Moria, may entail serious surgical consequences and damage the microkeratome.

Moria is not liable for malfunction or damage to the microkeratome, poor results or surgical complications caused by the reuse of a single-use product or products not supplied by Moria. The elements of this system should be connected only to Moria equipment (console, heads, rings, etc.). All warranties become null and void if the microkeratome degrades or malfunctions due to such practices.

2 LIST OF EQUIPMENT AND ASSOCIATED ACCESSORIES

Equipment described in this user guide	MORIA reference
CBm-ALTK head	19170/130 19170/150 19170/200 19170/250 19170/300 19170/350 19170/400
Box of 10 sterile CBm blades	19333
Box of 5 sterile CBSU heads for artificial chamber	19178/50 19178/90 19178/110 19178/130 19178/200 19178/250 19178/300 19178/350

Devices associated with the equipment described in this user guide	MORIA reference
Base of the ALTK artificial chamber	19161
Universal cover for ALTK artificial chamber	19162
Guide ring for CBm microkeratome	19171
Cover for ALTK + "CBm-ALTK" artificial chamber	19172
Single-use artificial chamber with guide ring for CBSU head	19179
CBSU guide ring for disposable artificial chamber	19180
Infusion tubing kit for artificial chamber	19181
Single-use artificial chamber	19182
Infusion set including syringe (for U.S. only)	19183
Pneumatic turbine for CB microkeratome	19303
CB suction ring (size +2, +1, 0, -1)	19309/2 19309/1 19309/0 19309/-1
Turbine Hose	19353
EVOLUTION 3E console (SN >5000)	19380
EVOLUTION 3E pedal	19381
EVOLUTION 3E pedal	19381J
EVOLUTION 3E pedal (for China)	19381C
EVOLUTION 3 pedal (also for EVOLUTION 3E SN>7000)	19361


Accessories	MORIA reference
Sterilization box for the turbine and turbine hose	22519139
ALTK sterilization box	22519167
ALTK flattening lens (7, 9, 11 mm)	19165
ALTK flattening lens (6, 8, 10 mm)	19166
ALTK flattening lens storage box	22519169
Tonometer	19042
Tonometer storage box	22519047
CB 7-9 flattening lens for recipient	19310
CB 8-10 flattening lens for recipient	19311

User information	MORIA reference
User guide for ALTK for artificial chambers for microkeratomes ONE & CBm-ALTK	65009
Selecting CBSU heads to prepare ultra-thin lamellar grafts	65082
User guide for EVOLUTION 3E (versions 5000 and later)	65060/INTL
User guide for EVOLUTION 3 (French/English version)	65037/65038
User guide for EVOLUTION 2 (French/English version)	65015/65016

3 REGULATORY INFORMATION

	MORIA S.A. 27, rue du Pied de Fourche – 03160 Bourbon l'Archambault – France Telephone: +33 (0) 4 70 67 09 04 Fax: +33 (0) 4 70 67 06 61 moria@moria-int.com http://www.moria-surgical.com
Customer Service Information	Contact your local distributor or Moria representative.
EUROPE	 0459 in compliance with European Directive MDD 93/42/EEC
USA	Registered with the Food and Drug Administration
Electrical Safety Standard	IEC 60601 Class II BF
	Exclusively for users in the Europe Union: this symbol indicates that, in the EU, the product at the end of its life must be discarded in an appropriate container. This applies not only to the product but also to all accessories, including electric motors and pedal driven machines, regardless of whether this symbol is displayed on all components or not. Do not discard into a public waste bin that does not indicate that it is part of a sorting program. Exclusively for users in the Europe Union: see the environmental standards for electrical and electronic waste for the countries concerned.

4 MANUFACTURER INFORMATION ON LABELS

 XXXXXX	CATALOGUE REFERENCE
	QUANTITY
 YYYY/MM/DD	USE BY
 XXXXXX	BATCH CODE
	STERILIZED USING ETHYLENE OXIDE
	DO NOT REUSE
 XXXXXXX	MANUFACTURER
 YYYY/MM/DD	DATE OF MANUFACTURE
 XXXXXX	CAUTION : CONSULT ACCOMPANYING DOCUMENT(S)
 XXXXXX	CONSULT OPERATING INSTRUCTIONS
	DO NOT USE IF PACKAGE IS DAMAGED
	STORAGE CONDITIONS: • TEMPERATURE: XX°C – YY°C / XX°F – YY°F • HUMIDITY RATE: XX% – YY%
	NON-STERILE
	KEEP DRY
	DISCARD IN A SEPARATE COLLECTION BIN
	CAUTION for USA only: US federal laws restrict this device to sale by, or on the order of, a physician.

5 IMPORTANT INFORMATION

5.1 DESCRIPTION

The CB is a mechanical microkeratome with a single-use or metallic (reusable) head.
It consists of:

- a CB pneumatic turbine (#19303)
- a single-use CBSU for artificial chamber including a pre-assembled blade, or reusable CBm-ALTK metallic head (to which a sterile blade needs to be attached).

The CB microkeratome is used with the EVOLUTION 2, 3 et 3E console. Please refer to the corresponding user guide for more information.

5.2 INDICATIONS

- Resection of donor cornea (graft mounted on the artificial chamber):
CBm-ALTK or CBSU heads for artificial chambers are indicated for a complete corneal lamellar (hingeless) resection on the graft.
- Corneal resection of the recipient's eye:
The CBm-ALTK head is indicated for performing a complete (hingeless) lamellar resection on the recipient's eye.

Important: CBm (#19332/xxx), CBSU (#19312/xxx) and CBSU heads for artificial chambers (#19178/xxx) have a stop that prevents a complete (hingeless) resection CB (#19309/xx) on the patient when used with CB rings (#19309/xx). The use of these devices entails creating a hinge on the recipient's eye. Only CBm-ALTK heads have no stops and allow a complete resection on the patient.

5.3 PRODUCT PERFORMANCE – THICKNESS OF CORNEAL FLAPS

The thickness of the corneal flap is a key factor in lamellar keratoplasty. Many parameters can affect the thickness of the flap and its predictability (standard deviation). Numerous scientific studies have shown the effect of many patient-related factors on the thickness of the corneal flap such as keratometry (K), the anatomy of the cornea, the preoperative pachymetry, intraocular pressure (and the pressure in the artificial chamber for preparing the donor graft), as well as factors related to the surgery itself such as corneal hydration and the travel speed of the mechanical microkeratome.

Ultrasound pachymetric measurements are not always accurate or reproducible, so results can vary depending on the technique and calibration of the equipment.

For CBm-ALTK heads, for example, a MORIA head calibrated as “130” cuts, on average, a corneal flap of 160 microns:

- “130” corresponds to a measurable dimension in the head of the microkeratome (distance between the blade and plateau).
- 160 microns is the average thickness that this head cuts with variations around this average.

Head (nominal)	Increment	Average depth of cut
130	+30	160
200	+50	250
250	+60	310
300	+70	370
350	+85	435
400	+100	500

For CBSU heads for artificial chambers, see the document «Selecting CBSU heads to prepare ultrathin lamellar grafts» (#65082).

Moria can only guarantee the dimensional characteristics of the head, not the surgical result. The surgeon is strongly advised to carry out tests beforehand to determine the thickness and diameters of the expected cuts based on his/her own technique.

5.4 CONTRAINDICATIONS

CBm-ALTK ((#19170/xxx) heads and CBSU heads for artificial chambers (#19178/xxx) must not be used in any way for a LASIK procedure due to the thickness of the resections obtained.

Important: CBm (#19332/xxx), CBSU (#19312/xxx) and CBSU heads for artificial chambers (#19178/xxx) have a stop that prevents a complete (hingeless) resection CB (#19309/xx) when used with CB rings (#19309/xx). The use of these devices entails creating a hinge on the recipient's eye.

5.5 WARNINGS

- **The design of CBm blades and CBSU heads for artificial chambers does not allow them to be used under any circumstances. Reusing these products degrades their clinical performance and exposes the patient to adverse reactions.**
- Do not use materials and/or components with the CB microkeratome that are not made by Moria.
- CBm blades must be inserted by hand. Do not use tools or pliers. Incorrect assembly can produce an incomplete or defective cut due to improper oscillation of the blade.

- The heads must be screwed by hand. Do not use tools or pliers. Incorrect assembly can produce an incomplete or defective cut due to improper oscillation of the blade.
- CB blades (#19333) are compatible with CBm-ALTK heads for artificial chambers.

Head	Reference	CBm blade (#19333)
CBm-ALTK	19170/xxx	X
CBSU for artificial chamber	19178/xxx	N/A

5.6 POTENTIAL ADVERSE EFFECTS

As with any surgical procedure, there is a risk. As lamellar keratectomy requires the use of a microkeratome, potential complications include and are not limited to those connected with corneal flaps (incomplete flap, hole in the flap, epithelial defect, micro-striations on the flap, folds in the flap, etc.).

Improper use and/or deterioration of the product and/or non-observance of contraindications (§ 5.4) or of warnings (§ 5.5) expose the patient to further adverse effects.

5.7 PRECAUTIONS

- A keratectomy should only be performed by an experienced ophthalmic surgeon who has received training in the CB microkeratome.
- Pre- and peroperative requirements, including knowledge of surgical techniques, of the head used, of the choice of ring, assembly and positioning of the microkeratome are important factors in the optimal use of the system by the surgeon.

Other pre- and peroperative precautions must also be taken:

• PREOPERATIVE:

- o Special care must be taken when handling, storing and maintaining the components of the microkeratome. They must not be knocked or damaged. The CB turbine must be protected during storage, especially in a corrosive environment.
- o Check the label and expiry date on the unopened peelable pack of CBSU heads for artificial chamber and CBm blades.
- o After opening a sterile pack of CBSU heads for artificial chambers, check that the markings on the heads match the information on the pack label.
- o The head, ring, CB turbine and console must be carefully inspected before use. The CBm blades (after insertion into the head #19170/xxx) and the CBSU blades must always be inspected under a microscope before use.
- o The surgeon must be familiar with the components before using the microkeratome, and must check that all the necessary parts and instruments are present and correctly assembled before beginning the surgical procedure. Spare sterile components should be available in case they are needed.
- o The selection of the head (and potentially the ring) is an important factor in the success of the surgery.

• PEROPERATIVE:

- o Breakage or misuse of the microkeratome and its components can injure the patient, operating room staff or damage the donor graft.
- o Before using, lubricate the suction ring, the rotating cover, the head and the blade with a physiological saline solution or any other appropriate ophthalmic solution that is compatible with the components of the CB.

• POSTOPERATIVE:

- o The surgeon's recommendations to the patient should be followed.

5.8 PRODUCT COMPLAINT

Any health professional (specifically, customer or user of the system) who has a complaint or has experienced dissatisfaction with the quality of the product, its description, durability, reliability, safety, effectiveness and/or performance must communicate it to Moria or to its distributor by telephone and confirm it in writing (indicating the product name, reference, serial or batch number, name and address, nature of the complaint, and patient characteristics). The faulty product(s), once decontaminated, must also be checked by Moria.

5.9 PACKAGING

The packaging of every component must be intact when received. All elements must be carefully checked to ensure that there is no damage before use. Products that are damaged or in damaged packaging must not be used but must be returned to Moria.

Do not use a single-use product if the packet is open or torn.

6 INSTALLATION AND CONNECTION

ACTIONS	 WARNINGS / PRECAUTIONS	PHOTOS
1. Select the head <u>1.1 CBm-ALTK:</u> select the desired pre-calibrated head. The calibration value is shown on the front of the head and determines the depth of the cut. <u>1.2 CBSU for artificial chamber:</u> see the selecting guide for CBSU heads for artificial chambers (#65082) for the selection of the CBSU head. The calibration value is shown on the front of the head and determines the depth of the cut.	Disposable CBSU heads for artificial chambers are supplied sterile and are single-use only. Check that the packaging is not damaged or open and that the use-by date has not expired.	A (CBm-ALTK) B (CBSU)
2. Insert the blade into the head <u>2.1 CBm-ALTK:</u> Check the general condition of the blade, particularly its cutting edge before and after inserting the blade into the head. The shape of the blade only allows it to be inserted in one way thus eliminating any risk of error. - The blade can only be inserted when the face is in line with the head. - Insert the blade gently into the slot provided for the purpose in the head. - Guide the blade so that the cutting edge does not come in contact with the metal of the head. - Push the blade in gently and completely. Light pressure should allow the blade to slip in easily. If that is not the case, check that there are no debris or particles in the slot. If there is resistance, use another blade. - Check that the blade holder is correctly centered in the head (the throat of the holder must be visible when viewed from top). - Moisten the CBm-ALTK thoroughly with the physiological serum or any other appropriate ophthalmic solution. <u>2.2 CBSU for artificial chamber:</u> - The CBSU heads for artificial chambers are supplied sterile with their blades pre-inserted to avoid any handling of the blade. - Open the peelable pack and place the blister on a clean surface without breaking the asepsis chain. Place the blister on a flat surface so that the cylindrical base is visible, then remove the top of the blister carefully.	- The head must be perfectly clean and must present no debris, particle, oxidation, deposit or scuffmark. - CBm blades are supplied sterile and are single-use only. Check that the packaging is not damaged or open and that the use-by date has not expired. - Check the general condition of the blade. Blades with damaged cutting edges must not be used.	C, D

3. Assemble the head of the CB turbine 3.1 CBm-ALTK: the CB turbine screws onto the head. Hold the head firmly by its two lateral faces and screw the CB turbine clockwise using finger strength only. 3.2 CBSU for artificial chamber: Holding the head in the blister, screw it onto the turbine. After assembly, remove the head from the blister and moisten it with physiological serum or any other appropriate ophthalmic solution.	• Before any use, check that the CB turbine is clean and in working order. • Check that the head is solidly screwed into the CB turbine. The head must not be able to move once mounted. • If the head is not correctly screwed to the CB turbine, the blade will not oscillate. A faulty cut may result from incorrect assembly.	E, F
4. Carefully inspect the head after assembling • Carefully inspect the blade using a microscope under high magnification to ensure that it is not damaged. If it is, discard it into an appropriate container and replace it.	• The head must be perfectly clean and present no debris, particle, oxidation, deposit or scuffmark.	
5. Connect the supply hose to the CB turbine and console • Connect the supply hose to the rear of the CB pneumatic turbine. The shape of the connector only allows it to be connected in one way (the two orifices are different dimensions). Screw it in using finger strength only. • Connect the turbine supply hose connector to the front of the console. The shape of the connector only allows it to be connected in one way (the two orifices are different dimensions). Screw it in using finger strength only.	• Check that the gaskets are present, as well as their appearance, on the connector on the CB pneumatic turbine and on the connector on the console. • Check the integrity of the hose (not cracked, pinched, etc.). Replace it if necessary. • Make sure that the tubing does not develop a kink.	G
6. Check the oscillation of the blade • Activate the vacuum then the turbine to check that the blade is properly engaged and oscillates correctly.	• See the user guide for the console used before using it and the user guide for artificial chamber (#65009). • Do not use the turbine if the blade does not oscillate smoothly and continuously.	
For anterior lamellar grafts (preparation of the recipient eye)		
7. Select the suction ring then inspect it • The rings need to have been sterilized in accordance with the recommendations in this user guide. • Select the appropriate suction ring (4 sizes available). The size of the ring determines the size of the cut.	• The suction rings must be carefully checked before every use. They must be perfectly clean, and present no debris, particle, oxidation, deposit or scuffmark. • Carefully check the vacuum connection. If obstructed, use another ring. • Do not use CBSU heads on the recipient eye (because they do not permit a complete resection).	H
8. Connect the suction ring to the vacuum tube • Connect the suction ring to the vacuum tube • Position the ring on the eye • Activate the vacuum and check that the ring is properly centered and that the eye is held by the suction.	• The single-use vacuum tube (#19138) is supplied sterile. • Check that the packaging is not damaged or open and that the use-by date has not expired. • Inspect the vacuum tube. If it is damaged or obstructed, use another vacuum tube.	I, J
9. Select the flattening lens then check the size of the resection • The flattening lenses (transparent and mounted on a sleeve) contain two reticles indicating two different diameters. They are used with the suction rings to check the size of the resection obtained with the CBm-ALTK head. • Select the theoretically appropriate lens for the ring used. • Change the size of ring, if necessary, to suit the size of the resection. Check again using the lenses.	• The flattening lenses must be carefully checked before every use. They must be perfectly clean, and must present no debris, particle, oxidation, deposit or scuffmark. • They need to have been decontaminated and disinfected in accordance with the recommendations in this user guide.	K

7 OPERATION

7.1 LAMELLAR RESECTION ON A DONOR EYE (GRAFT)

ACTIONS	 WARNINGS / PRECAUTIONS
1. Check the equipment before every use It is essential to make sure that the equipment is perfectly clean and that it presents no debris, particle, oxidation or deposit.	
2. Connect the CB turbine to the console Check the pressure delivered by the console at the turbine (operating range between 2.5 and 3.5 bars).	See the user guide for the console employed before use.
4. To prepare the donor eye: Mount the CB turbine on the ALTK artificial chamber or on the single-use artificial chamber.	- See the ALTK artificial chamber User Guide for microkeratomes ONE & CBm-ALTK (#65009) before use. - See the selecting guide for CBSU heads for artificial chambers (#65082).
5. Activate the vacuum by pressing once on the "Vide / Vacuum" pedal.	Check that the equipment does not interfere with other medical equipment used during the surgery.
6. Check the intraocular pressure with a tonometer Moria recommends a pressure of 90mm Hg.	For CBSU heads for artificial chambers: see the selecting guide for CBSU heads for artificial chambers (#65082).
7. To prepare the donor eye: Lubricate the sliding guides for the ring for the artificial chamber, the head and the blade.	- Before use, lubricate the various elements with physiological serum or any other appropriate ophthalmic solution. - The equipment must not be in contact with a solution with unknown interactions.
8. Activate the oscillation of the blade and perform the resection - Once the microkeratome is in position, activate the vacuum and the oscillation of the blade by pressing and holding down the "Avance / Forward" pedal. - Rotate the microkeratome gently and smoothly by hand until a free corneal flap (or "free cap") is created. - Depending on the thickness of the desired graft, a second cut may be made. Note that on the ALTK chamber (#19161) it is possible to adjust the diameter by using lenses #19165 & #19166. The adjustment ring on chamber #19162 allows the diameter of the cut to be adjusted.	- See the user guide for the console used. - See the ALTK artificial chamber user guide for microkeratomes ONE & CBm-ALTK (#65009) before use. - See the selecting guide for CBSU heads for artificial chambers (#65082). - Check that nothing interferes with the movement of the head in the ring sliding guides of the artificial chamber.
8. Release the vacuum by pressing once on the "Vide / Vacuum" pedal.	At the end of the procedure, discard the CBm blade or single-use head (CBSU) as well as the vacuum tube, into an appropriate container.

7.2 RESECTION ON THE PATIENT

ACTIONS	 WARNINGS / PRECAUTIONS
1. Check the equipment before every use It is essential to make sure that the equipment is perfectly clean, and that it presents no debris, particle, oxidation or deposit.	
2. Connect the CB turbine to the console Check the pressure delivered by the console at the turbine (operating range between 2.5 and 3.5 bars).	See the user guide for the console employed before use.
3. Check the oscillation of the blade - Activate the turbine to check that the blade is properly engaged and oscillates correctly. - Check that the CBm-ALTK or CBSU head slides freely in the ring (do one test pass with the turbine in the ring of the chamber before placing the graft, or on the suction ring of the recipient eye before placing it on the eye).	Do not use the turbine if the blade does not oscillate smoothly and continuously.
4. Position the suction ring on the patient - Once the steps have been carried out to install and connect the elements, position the suction ring on the patient's eye. - Activate the vacuum and check that the ring is properly centered and held on the eye. - Check the intraocular pressure with a tonometer. Moria recommends a pressure of 60mmHg. - Lubricate the suction ring, the head and the blade.	- Before use, lubricate the various elements with physiological serum or any other appropriate ophthalmic solution. - The equipment must not be in contact with a solution with unknown interactions. - Check that the equipment does not interfere with other medical equipment used during the surgery. - The tonometer (#19042) must be perfectly dry and must be used on a dry surface.
5. Activate the oscillation of the blade and perform the resection - Once the microkeratome is in position, activate the vacuum and the oscillation of the blade by pressing and holding down the "Avance / Forward" pedal. - Rotate the microkeratome gently and smoothly by hand until a free corneal flap (or "free cap") is created.	- See the user guide for the console used. - See the ALTK artificial chamber user guide for microkeratomes ONE & CBm-ALTK (#65009) before use. - See the selecting guide for CBSU heads for artificial chambers (#65082). - Check that nothing interferes with the movement of the head in the ring sliding guides of the artificial chamber. - For the recipient, make sure that there are no obstructions to the movement of the head (speculum, eyelash, eyelid, etc.).
6. Release the vacuum by pressing once on the "Vide / Vacuum" pedal.	At the end of the procedure, discard the CBm blade or single-use head (CBSU) as well as the vacuum tube, into an appropriate container.

8 IN THE EVENT OF A PROBLEM

For further information and addition descriptions, see the user guide for EVOLUTION 2, 3 or 3E and/or the user guide for the ALTK artificial chamber for microkeratomes ONE & CBm-ALTK and single-use artificial chamber.

9 MAINTENANCE

The recommendations for cleaning, disinfection and sterilization are updated and available on Moria's restricted-access website: <http://www.moria-surgical.com>. You must register to access it.

9.1 INTRODUCTION

If there is unusual vibration or abnormal noise, do not use the equipment and contact the Moria representative. The following are only recommendations. They must be adapted to the regulations in the country in which

the equipment is used. All cleaning and disinfection products and solutions must be used in accordance with the manufacturer's instructions.

Moria recommends:

- Do not touch potentially contaminated zones.
- Use gloves during cleaning and decontamination operations.

9.2 INITIAL TREATMENT AND STORAGE

To avoid the risk of condensation inside packaging, close it and store the equipment in a clean dry environment. Do not store instruments in an environment, or near other products, that could cause magnetic effects or create corrosion. Avoid all contact between items of equipment.

The equipment must be damage and defect free.

Equipment that is fragile or requires special handling must be handled separately taking care not to damage delicate parts.

For instructions on disinfecting, cleaning and sterilizing the ALTK artificial chamber, tonometer and flattening lenses for the donor, see the user guide for the "ALTK Artificial Chamber for microkeratomes ONE & CBm-ALTK".

9.3 PROTOCOL FOR CLEANING / DISINFECTING / STERILIZING / STORING THE TURBINE

The pneumatic turbine must be disconnected from the turbine supply hose and from the head.
It must be carefully cleaned, disinfected, dried and sterilized after every surgical procedure.

ACTIONS	 WARNINGS / PRECAUTIONS
1. Manual external cleaning • After use, clean the outside of the turbine with a lint-free cloth dampened with a cleaning solution (eg: Alkazyme® / Laboratoires Alkapharm) to remove any residue.	• Do not use abrasives or scrapers to clean items. It could damage their usability. • Do not immerse the turbine in any solution whatsoever, including water. • Cleaning is not disinfecting! Disinfection is a separate step described in point 2.
2. Manual external disinfection • After cleaning, dry every part of the turbine with a lint-free cloth dampened with a disinfecting solution.	• Use appropriate disinfection solutions in accordance with the manufacturer's instructions. • Do not immerse the turbine in any solution whatsoever, including water.
3. Lubrication then cleaning before sterilization • Lubricate the handpiece using a KaVo spray only and the appropriate nozzle before every autoclave sterilization cycle. A single jet of KaVo spray is sufficient to lubricate the turbine while avoiding excess lubricant inside the turbine. • Clean the turbine with a lint-free cloth then keep it upside down for 20 minutes before the autoclave sterilization to expel excess lubricant.	• Use the recommended spray. Do not use other lubricant sprays, vaseline or grease. • Shake the spray before using. • Place the spray nozzle on the smallest hole of the turbine and press once only. • KaVo spray (ref 411-9630) is not sold by Moria. It can be purchased from Kavo website (www.kavo.com).
4. Sterilization • Place the turbine and its hose in the dedicated sterilization box. • Sterilize the turbine in the autoclave (134°C – 273°F) for 18 minutes. Sterilization in USA only: sterilize with an autoclave (steam sterilizer with a vacuum phase) at a minimum of 132°C for at least 4 minutes.	• Remove the turbine from the sterilizer immediately after the sterilization procedure. • Never use the turbine immediately after sterilization. Let it cool down to room temperature first.
5. Storage • Once this sterilization cycle is completed, remove the handpiece as promptly as possible and store it upright until it cools down and is dry. • Once cooled down and dry, the handpiece can be stored in an enclosed dry place.	• To avoid oxidation, make sure that the turbine is perfectly dry before putting it in the storage box. • While it is not being used, the CB turbine must be kept in its storage box in a dry atmosphere (Photo L). • Before use, check that the handpiece has no lubricant on it. If it has, do not use it. Do not lubricate before use.

9.4 PROTOCOL FOR CLEANING / DISINFECTING / STERILIZING / STORING THE TURBINE HOSE

The turbine supply hose must be disconnected from the pneumatic turbine and console, cleaned, disinfected, dried and sterilized after every use.

ACTIONS	 WARNINGS / PRECAUTIONS
1. Cleaning Clean the outside of the turbine supply hose with a lint-free cloth dampened with a disinfecting solution.	Do not immerse.
2. Disinfection Disinfect the outside of the turbine supply hose with a disinfectant (disinfectant spray or pad).	Do not immerse.
3. Rinsing and drying - Rinse the outside of the turbine supply hose with a disposable lint-free cloth soaked in sterile water. - Dry the turbine supply hose with a disposable lint-free cloth.	
4. Sterilization - Sterilize the turbine supply hose in the autoclave (134°C – 273°F) for 18 minutes. Sterilization in USA only: sterilize with an autoclave (steam sterilizer with a vacuum phase) at a minimum of 132°C for at least 4 minutes. - Remove it from the sterilizer immediately after the sterilization procedure. - Dry the turbine supply hose meticulously after sterilization with a disposable lint-free cloth.	Never use the turbine supply hose immediately after the sterilization is completed. Let it cool down to room temperature first.
5. Storage Store the turbine supply hose in a dry, dust-free environment.	To prevent the turbine supply hose from deforming, it is important to coil it carefully in its box (Photo L).

9.5 PROTOCOL FOR CLEANING / DISINFECTING / STERILIZING / STORING CBm-ALTK HEADS AND SUCTION RINGS

ACTIONS	 WARNINGS / PRECAUTIONS
1. Double cleaning - After taking it out of the block, immerse the device immediately in an aldehyde-free detergent-disinfectant solution (eg: Alkazyne® / Laboratoires Alkapharm) for at least 15 minutes in accordance with the manufacturer's instructions. - Clean by hand (or by machine provided that the machine does not recycle the cleaning products) after disassembling if necessary. - Rinse with water. - Immerse in an aldehyde-free disinfectant-detergent solution (different bath from the previous). - Clean by hand (or by machine provided that the machine does not recycle the cleaning products). - Rinse with filtered or osmotic water. - Dry the device with a single-use, lint-free non-woven soft cloth.	- Manual cleaning: clean the medical device mechanically using soft brushes (do not use metal brushes or scouring pads). - IMPORTANT: Change the bath water after every use. - Incompatibility: - Do not use bleach. - Do not use this protocol for products made of aluminum alloy. - Do not use this protocol on heat-sensitive microkeratomes (motor and turbine). - Inactivation of Non-Conventional Transmissible Agents (NCTA) only for patients at risk for Creutzfeld-Jacob disease. For France, see instruction no. DGS/ RI3/2011/449 of 1st December 2011.
2. Sterilization - Sterilize with an autoclave (steam sterilizer) at a minimum of 134°C – 273°F for 18 minutes sterilization plateau. - Sterilization in USA only: sterilize with an autoclave (steam sterilizer with a vacuum phase) at a minimum of 132°C for at least 4 minutes.	

French bibliography:

- December 2011 – Instruction no. DGS/RI3/2011/449 of 1st December 2011 updating recommendations aimed at reducing the risk of transmitting non-conventional transmissible agents during invasive procedures.
- November 2005 – Ministerial Guide to the treatment of medical devices in ophthalmology and contactology - <http://www.sante.gouv.fr>.

10 WARRANTY

10.1 SCOPE OF WARRANTY

Designation	MORIA references
CB suction ring (size -1, 0, +1, +2)	19309/-1 19309/0 19309/+1 19309/+2
CBm-ALTK head	19170/130 19170/150 19170/200 19170/250 19170/300 19170/350 19170/400
Turbine hose	19353
Pneumatic turbine for CB microkeratome	19303
Sterilization box for the turbine and turbine hose	22519139
EVOLUTION 3E console (excluding battery)	19380
EVOLUTION 3E pedal	19381
EVOLUTION 3E pedal (for China)	19381C
EVOLUTION 3E pedal	19381J
EVOLUTION 3 pedal (also for EVOLUTION 3E SN>7000)	19361

- Warranty covers parts and labor necessary for the proper operation of the equipment cited above. All returned equipment must be sent in its original packaging, after first having been decontaminated.
- The warranty does not cover single-use instruments or single-use accessories.
- Parts must be serviced and replaced exclusively by personnel approved by Moria.
- All defective parts exchanged during the warranty period remain the exclusive property of Moria.

10.2 WARRANTY EXCLUSIONS

- Fault or malfunction of the system that appears outside the warranty period (as defined in section 10.3).
- Normal equipment wear and tear.
- Negligence, misuse or any use that does not comply with the specifications in the user guide.
- The use of consumables, parts or accessories not supplied by Moria.
- Dismantling, modification or intervention on the equipment by a person not authorized by Moria.

10.3 WARRANTY PERIOD

- The warranty takes effect on the day that the equipment is shipped.
- The warranty is for 12 months from the effective date.

10.4 LIABILITY

- Moria's liability is limited to the supply of the services specified in section 10.1. Moria is not liable for direct or consequential damage, including financial, suffered by the customer due to interventions under this warranty.
- Any dispute over the interpretation or performance of this contract or over the general terms and conditions shall be within the competence of the Nanterre Commercial Court (France).

11 PHOTOS

A. CBm-ALTK HEAD

- 1 Head calibration value
- 2 Female collar (base)
- 3 Blade slot
- 4 Sliding guide

B. CBSU HEAD FOR ARTIFICIAL CHAMBER

- 1 Head calibration value
- 2 Female collar (base)

C. CBm BLADE

- 1 Cutting edge
- 2 Blade holder

D. CBm-ALTK HEAD WITH CBm BLADE

- 1 CBm-ALTK head
- 2 CBm blade

E. ASSEMBLY: CBm-ALTK HEAD ON THE CB TURBINE

- 1 CB turbine
- 2 CBm-ALTK head

F. ASSEMBLY: CBSU HEAD FOR ARTIFICIAL CHAMBER ON THE CB TURBINE

- 1 CB turbine
- 2 CBSU head for artificial chamber
- 3 CBSU head blister for artificial chamber

G. ASSEMBLY: TURBINE AIR SUPPLY HOSE

- 1 Turbine air supply hose – turbine side
- 2 CB turbine
- 3 CBm-ALTK head
- 4 Front face of the console
- 5 Turbine air supply hose – console side

H. CBSU MANUAL SUCTION RING

- 1 Handle
- 2 Stop pin
- 3 Smooth pivot (manual ring)
- 4 Sliding guide
- 5 Aperture

I. CONNECTING THE SUCTION RING TO THE ASPIRATION TUBING

- 1 Aspiration tubing end
- 2 Suction ring

J. MOUNTING CB MICROKERATOME ON THE SUCTION RING

- 1 CB turbine with CBm-ALTK head
- 2 Suction ring

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