
- Please unfold the following pages containing the photos quoted in the reference.
- Veuillez déplier les pages suivantes contenant les photos citées en référence.
- Por favor, despliegue las páginas siguientes que contengan las fotos referenciadas.
- Vogliate aprire le seguenti pagine, contenenti le foto in oggetto.
- Falten Sie die folgenden Seiten bitte auf, die die in Verweis zitierten Fotos enthalten.
- Sla volgende pagina's open voor de foto's waar in de tekst naar verwezen wordt.
- Por favor desdobre as páginas seguintes que contêm as fotografias citadas em referência.
- Följande sidor visar bildreferanser.
- Lütfen atıfta bulunan resimleri görmek için sayfaları çevirin.

# Moria

**MORIA S.A. – 27, rue du Pied de Fourche - 03160 BOURBON L'ARCHAMBAULT – France**  
#65044-C-05.2023

Document Aletiq : Rév. C-05.2023 (Validée)

|  |       |   |       |
|--|-------|---|-------|
| <b>EN - Epi-K™ Microkeratome</b><br><b>With Single-use plastic ring</b><br>INSTRUCTION MANUAL.....   | 1-12  | <b>NL - Epi-K™</b><br><b>Plastic ring voor eenmalig gebruik</b><br>HANDLEIDING.....                   | 61-72 |
| <b>FR - Microkératome Epi-K™</b><br><b>avec anneau plastique</b><br>MANUEL UTILISATEUR.....  | 13-24 | <b>PT - Epi-K™</b><br><b>Anel plástico descartável</b><br>MANUAL DE INSTRUÇÕES.....                   | 73-84 |
| <b>ES - Epi-K™</b><br><b>Anillo de plástico de uso único</b><br>MANUAL DE INSTRUCCIONES.....   | 25-36 | <b>SV - Epi-K™</b><br><b>Plastring för engångsanvändning</b><br>INSTRUKTIONSMANUAL.....               | 85-96 |
| <b>IT - Epi-K™</b><br><b>Testina monouso in plastica con lama</b><br><b>premontata e anello di suzione in plastica</b><br><b>monouso</b><br>MANUALE DELL'UTENTE..... | 37-48 | <b>TR - Epi-K™</b><br><b>Tek kullanımlık plastik halka</b><br>KULLANIM TALİMATLARI EL KİTABI...97-108 |       |
| <b>DE - Epi-K™</b><br><b>Einweg-Kunststoffring</b><br>BENUTZERHANDBUCH.....  | 49-60 |   |       |

**Ref.: 65044**

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













#65044-C-05.2023






# CONTENTS

|           |  |           |
|-----------|--|-----------|
| <b>1</b>  | <b>DISCLAIMER</b> .....  | <b>2</b>  |
| 1.1       | CARE AND MAINTENANCE .....   | 2         |
| 1.2       | USE OF GENERIC PRODUCTS AND REUSE OF SINGLE-USE CONSUMABLES .....          | 2         |
| 1.3       | CREATION OF AN EPITHELIAL FLAP .....                                       | 2         |
| <b>2</b>  | <b>EQUIPMENT LIST</b> .....  | <b>3</b>  |
| <b>3</b>  | <b>REGULATORY INFORMATION</b> .....  | <b>3</b>  |
| <b>4</b>  | <b>LABELING INFORMATION</b> .....  | <b>4</b>  |
| <b>5</b>  | <b>IMPORTANT INFORMATION</b> .....   | <b>5</b>  |
| 5.1       | DESCRIPTION .....  | 5         |
| 5.2       | INDICATIONS .....  | 5         |
| 5.3       | CONTRAINDICATIONS .....  | 5         |
| 5.4       | WARNINGS .....   | 5         |
| 5.5       | POTENTIAL ADVERSE EVENTS .....   | 5         |
| 5.6       | PRECAUTIONS .....  | 6         |
| 5.7       | PRODUCTS COMPLAINTS .....  | 6         |
| 5.8       | PACKAGING .....  | 6         |
| <b>6</b>  | <b>INSTALLATION AND CONNECTION</b> .....                                   | <b>7</b>  |
| <b>7</b>  | <b>OPERATION</b> .....   | <b>8</b>  |
| <b>8</b>  | <b>TROUBLESHOOTING</b> .....   | <b>9</b>  |
| <b>9</b>  | <b>CARE &amp; MAINTENANCE</b> .....  | <b>9</b>  |
| 9.1       | INTRODUCTION .....   | 9         |
| 9.2       | INITIAL TREATMENT AND STORAGE .....  | 9         |
| 9.3       | CLEANING- DISINFECTION-STORAGE OF THE HANDPIECE AND ELECTRICAL CABLE ..... | 10        |
| <b>10</b> | <b>WARRANTY</b> .....  | <b>11</b> |
| 10.1      | SCOPE OF WARRANTY .....  | 11        |
| 10.2      | NON-APPLICATION OF WARRANTY .....  | 11        |
| 10.3      | WARRANTY PERIOD .....  | 11        |
| 10.4      | LIABILITY .....  | 11        |
| <b>11</b> | <b>DRAWINGS</b> .....  | <b>12</b> |

# 1 DISCLAIMER

## 1.1 CARE AND MAINTENANCE

The Epi-K™ system has been designed for optimal operation, provided that the recommendations listed in this user manual are followed carefully. If, for any reasons, the system does not perform properly, have it checked immediately by MORIA. In order to maintain the original performance of your microkeratome, MORIA strongly recommends annual maintenance and servicing of all its reusable products.

As only MORIA and its agents are fully expert in MORIA products, servicing and maintenance must be carried out by MORIA or its approved agents.

MORIA shall not be held liable for any malfunction or damage to the apparatus, poor results, or surgical complications due to maintenance being having been carried out by an unqualified operator or third party.

Any such unauthorised intervention shall render the guarantee and any maintenance contract null and void.

## 1.2 USE OF GENERIC PRODUCTS AND REUSE OF SINGLE-USE CONSUMABLES

The materials used in the microkeratome blade, blade holder, and head have been selected for their biocompatibility and sliding characteristics. The dimensions and tolerances of the separator have been determined in consideration of the dimensions and tolerances of the head of the epikeratome. MORIA's manufacturing and inspection procedures guarantee there is no conflict in dimensions between head and separator, and that the separator will slide smoothly in the head.

**Single-use devices should not be re-used. Doing so will negatively affect their clinical performance and increase the potential for adverse events.**

**The reuse of single-use products, or the use of consumables other than those supplied by MORIA, may entail serious surgical consequences for the patient and damage the microkeratome.**

MORIA shall not be held liable in the event of a malfunction or damage to the microkeratome, poor results or surgical complications due to the reuse of a single-use product or the use of consumables other than those supplied by MORIA.

**MORIA handpieces must only be connected to MORIA devices (console unit, heads, suction rings, etc.).**

All warranties become null and void if the microkeratome degrades or malfunctions due to such practices.

## 1.3 CREATION OF AN EPITHELIAL FLAP

Creation of an epithelial flap depends on a number of factors relating to the product (separator integrity, and advance rate) as well as patient characteristics and surgical technique. MORIA strongly recommend that surgeons attend an Epi-K™ training session to gain familiarity with the device and the technique. MORIA can only guarantee the dimensional characteristics of the system, not the surgical result.

For the EVOLUTION 3E console with serial number below 5000, please refer to instruction manual #65051. For the EVOLUTION 3E console with serial number 5000 and above, please refer to instruction manual #65060/INTL.

The most recent version of this user guide and additional information on your keratome are available on the MORIA website: <http://www.moria-surgical.com>.


## 2 EQUIPMENT LIST

| Designation  | MORIA reference |
|--|-----------------|
| Box of 10 Epi-K™ ring -1 and 10 Epi-K™ heads   | 19384           |
| Box of 10 Epi-K™ ring 0 and 10 Epi-K™ heads  | 19385           |
| Epi-K™ handpiece (green color)   | 19342           |
| Wrench   | 19345C          |
| Storage box for the Epi-K™ handpiece   | 22519512        |
| Tonometer  | 19042           |
| EVOLUTION 3E console / footswitch Epi-K™   | 19380/19381     |
| Aspiration tubing  | 19138           |
| Disposable brush   | 19149           |
| Epi-K™ single use plastic ring user manual   | 65044           |
| EVOLUTION 3E console user manual (serial number below 5000)                              | 65051           |
| EVOLUTION 3E console user manual (serial number 5000 and above)                          | 65060/INTL      |
| Annexe "Guidance and manufacturer's declaration: electromagnetic emissions and immunity" | 65073           |
| Guidelines (A4 format)   | 65068           |
| Guidelines (pocket format)   | 65069           |
| Epi-K™ protocol, drug regimen and contact lens with or without the flap                  | 65071           |

## 3 REGULATORY INFORMATION

|  |  |
|--|--|
|    | MORIA S.A.<br>27, rue du Pied de Fourche - 03160 Bourbon l'Archambault - France<br>Téléphone : +33 (0) 4 70 67 09 04<br>Fax : +33 (0) 4 70 67 06 61<br>moria@moria-int.com<br><a href="http://www.moria-surgical.com">http://www.moria-surgical.com</a>  |
| Customer Service Information   | Contact your local dealer or MORIA   |
| EUROPE   |  0459 CLASS IIA according to MDD 93/42/CEE  |
| USA  | Product registered at the Food and Drug Administration (FDA): 510(k) K043183.<br>Caution USA only: Federal law restricts the use of this device to physicians or licensed practitioners.   |
| Electrical Safety Standard   | IEC 60601  |
|  | For EU customers only: this symbol indicates that within the European Union, the product must be discarded in a separate collection bin at the end of its useful life. This applies not only to this device, but also all accessories, including the footpedal and electrical motors, regardless of whether those accessories are marked with the symbol. Do not dispose as unsorted municipal waste.<br>For users out of European Union: please refer to local environmental regulatory regarding waste of electrical and electronic equipment. |

## 4 LABELING INFORMATION

|   |  |
|---|--|
|  REF XXXXXX | CATALOGUE REFERENCE  |
|             | QUANTITY   |
|  YYYY/MM/DD | USE BY   |
|  LOT XXXXXX | BATCH CODE   |
|             | STERILIZED USING ETHYLENE OXIDE  |
|             | DO NOT REUSE   |
|  XXXXXXX    | MANUFACTURER   |
|  YYYY/MM/DD | DATE OF MANUFACTURE  |
|  XXXXXX     | CAUTION : CONSULT ACCOMPANYING DOCUMENT(S)   |
|  XXXXXX     | CONSULT OPERATING INSTRUCTIONS   |
|             | DO NOT USE IF PACKAGE IS DAMAGED   |
|  °C °F %    | STORAGE CONDITIONS: <ul style="list-style-type: none"> <li>• TEMPERATURE: XX°C – YY°C / XX°F – YY°F</li> <li>• HUMIDITY RATE: XX% – YY%</li> </ul> |
|             | NON-STERILE  |
|            | KEEP DRY   |
|           | DISCARD IN A SEPARATE COLLECTION BIN   |
|           | CAUTION for USA only: US federal laws restrict this device to sale by, or on the order of, a physician.  |

## 5 IMPORTANT INFORMATION

### 5.1 DESCRIPTION

The Epi-K™ epikeratome is used to create a corneal epithelial sheet during an Epi-LASIK procedure, prior to photoablation of the stroma with an excimer laser.

The Epi-K™ epikeratome is an automated, mechanical, linear epikeratome designed to create epithelial flaps. It has three components:

- a handpiece (#19342), containing two independent motors: one for advancement, and one for oscillation.
- a single-use suction ring.
- a single use plastic head encasing a pre-inserted separator.

The Epi-K™ epikeratome operates with the EVOLUTION 3E control unit (#19380).

Please refer to user guide:

- #65051 (EVOLUTION 3E serial number below 5000)
- #65060/INTL (EVOLUTION 3E serial number 5000 and above).

### 5.2 INDICATIONS

The Epi-K™ epikeratome is indicated for the cleavage of the epithelium from the Bowman's membrane in patients suitable for surface ablation with intact Bowman's membrane, and keratomeries between 39 D and 48 D.

### 5.3 CONTRAINDICATIONS

- Patients who have had any previous surgery that could have damaged Bowman's membrane, including but not limited to RK, PRK, LASEK, PTK, or Epi-LASIK.
- Patients with any pathologies such as pterygium, keratoconus, corneal ulcers, corneal scars etc. that could have damaged Bowman's membrane.
- Patients who have had a foreign body in cornea.

### 5.4 WARNINGS

- **Design of this medical device does not allow its reuse. Reusing this medical device affects its clinical performance and exposes the patient to adverse events.**
- Do not mix Epi-K™ heads with One Use-Plus heads.
- Do not use disposable materials and/or components or a brand other than MORIA with the Epi-K™ epikeratome.
- Do not mix heads and rings from separate sterile packs. Each disposable head must be used with the disposable ring from the same pack.
- The heads must only be screwed on by hand. Dismantling should be done only with the provided wrench (#19345C) or by hand. Never use any tools or other spanners. Incorrect assembly may cause incomplete or uneven cuts due to lack of separator oscillation.
- Never pull on the cable connected to the handpiece, and never hold the motor by the cable.

### 5.5 POTENTIAL ADVERSE EVENTS

As with any surgical procedure, there is risk involved. Epi-LASIK surgery requiring the use of an epikeratome that cuts an epithelial flap, potential side effects of laser refractive surgery may include but are not limited to: visual anomalies, dry eye and flap related complications (incomplete flap, buttonhole, epithelial defect, flap dislocation, flap striae, wrinkles, etc.).

Inappropriate use, deterioration of epikeratome and/or non-respect to contraindications (§ 5.3) and warnings (§ 5.4) expose patient at higher risk to adverse events.

Note: Additional surgery may be necessary to correct some of these potential adverse events.

## **5.6 PRECAUTIONS**

- Epi-LASIK must be performed only by experienced refractive surgeons with specific training of the Epi-K™ epikeratome.
- Preoperative and operating procedures, including knowledge of surgical techniques, suction ring and stop selection, assembly and placement of the epikeratome are important considerations in the successful use of the system by the surgeon. Furthermore, proper patient selection and his compliance will greatly affect the results.

### **For USA only**

**CAUTION: FEDERAL LAW (USA) RESTRICTS THESE DEVICES TO SALE BY OR ON THE ORDER OF A PHYSICIAN.**

Other preoperative, intraoperative and postoperative warnings and precautions are as follows:

#### **• PREOPERATIVE**

- o Only patients who meet the criteria described in the indications should be selected.
- o Care should be used in the handling and storage of the epikeratome components. They should not be scratched or otherwise damaged. Handpieces should be protected during storage, especially from corrosive environments.
- o Check the label and expiration date on the unopened package of heads and rings.
- o After opening the package, verify that head and ring information is consistent with information on outer package labelling.
- o The head, suction ring, handpiece and control unit should be fully inspected prior to use. Both surfaces of the pre-inserted separator should be inspected under a microscope.
- o The surgeon should be familiar with the various components before using the epikeratome and verify that all parts and necessary instruments are present and properly assembled before the surgery begins. Additional sterile components should be available in case of an unexpected need.
- o The selection of the proper ring and stop position for each eye is crucial to the success of the procedure.

#### **• INTRAOPERATIVE**

- o Breakage, slippage, or misuse of epikeratome and its components may cause injury to the patient or operating personnel.
- o Before any usage and during the whole epithelial flap creation, lubricate the eye, the suction ring, the head and separator with a physiological saline solution or another appropriate ophthalmic solution that is compatible with the Epi-K™ components.

#### **• POSTOPERATIVE**

- o The surgeon's postoperative instructions to the patient and the corresponding patient compliance are extremely important.

## **5.7 PRODUCTS COMPLAINTS**

Health care professionals (e.g., customers or users of this device) having any complaints about or dissatisfaction with the product quality, identity, durability, reliability, safety, effectiveness and/or performance, should notify MORIA or its distributor by telephone, fax or written correspondence, and have the products checked by MORIA.

When filing a complaint, please provide the component(s) name(s), reference(s), batch number(s), as well as your name and address, the nature of the complaint, and the relevant patient data. Also disinfect and return the components.

## **5.8 PACKAGING**

Packages of each of the components should be intact upon receipt. If a loaner or consignment system is used, all sets should be carefully checked for completeness, and all components should be carefully checked to ensure that there is no damage prior to use.

Damaged packages or products should not be used and should be returned to MORIA.

## 6 INSTALLATION AND CONNECTION

| Steps | What to do  |  <b>IMPORTANT WARNINGS</b>  | Related picture |
|-------|---|--|-----------------|
| 1     | Select the suction ring and the head                  | <ul style="list-style-type: none"> <li>• The disposable head and suction ring are supplied sterile and for single use only.</li> <li>• Check that package is undamaged, unopened and within the expiration date.</li> <li>• <b>WARNING: DO NOT MIX HEADS &amp; RINGS OF DIFFERENT PACKS.</b></li> <li>• <b>EACH HEAD MUST BE USED WITH THE RING FROM THE SAME PACK.</b></li> </ul>   | A, B            |
| 2     | Assemble the head onto the handpiece                  | <ul style="list-style-type: none"> <li>• Open the package containing the disposable head and suction ring.</li> <li>• Screw the head clockwise onto the threaded motor shaft.</li> <li>• Check that the head is screwed on tightly to the motor. There should be no rotation of the head.</li> <li>• If the head is not screwed completely onto the motor, the shaft will not drive the blade holder, and the blade will not move. Poor or irregular resections can result from incorrect assembly.</li> <li>• Do not use the wrench (#19345C) or other tools to assemble it. The wrench should be used for disassembly only.</li> </ul> | E, F            |
| 3     | Carefully inspect the epikeratome head after assembly | <ul style="list-style-type: none"> <li>• The disposable head must be perfectly clean inside and outside and free from debris, particles, oxidation and deposits. There should be no scratches on the plates. If any, replace the entire pack (ring and head).</li> <li>• After assembly and using a microscope, carefully inspect the separator under high magnification to ensure that it is not damaged.</li> </ul>  | A               |
| 4     | Check for separator oscillation                       | <ul style="list-style-type: none"> <li>• Start the motor to check proper blade oscillation</li> <li>• Do not use the Epi-K™ system if the oscillation is not smooth, regular and uninterrupted.</li> </ul>   |                 |
| 5     | Inspect the suction ring                              | <ul style="list-style-type: none"> <li>• The suction ring must be carefully inspected prior to use; it must be perfectly clean, free from debris, particles and deposits.</li> <li>• Carefully check that the suction ports are not obstructed.</li> <li>• Inspect the tubing of the ring, and in case of kinks or obstructions, replace the entire pack (ring and head).</li> </ul>   | B               |
| 6     | Set the stopper                                       | <ul style="list-style-type: none"> <li>• Refer to the indicative guidelines for the selection of the ring and stop.</li> <li>• The adjustable stop pin is pre-mounted on the suction ring and used to determine the size of the hinge. The hinge will be placed in a nasal position.</li> <li>• Use the plastic spanner included in the Epi-K™ pack (ring and head) to select the appropriate stop value.</li> <li>• The value of the stop selected must line up with the mark on the suction ring.</li> </ul>   | C               |
| 7     | Connect the suction ring to the aspiration tubing     | <ul style="list-style-type: none"> <li>• The aspiration tubing (#19138) is supplied sterile and for single use only.</li> <li>• Check that package is undamaged, unopened and within the expiration date.</li> <li>• Inspect the aspiration tubing and in case of kinks or obstructions, replace it.</li> </ul>  | D               |

## 7 OPERATION

| Steps | What to do   |  <b>IMPORTANT WARNINGS</b>   | Related picture |
|-------|--|---|-----------------|
| 1     | Preliminary check before any usage   | <ul style="list-style-type: none"> <li>It is imperative to verify that the devices are perfectly clean and free of any organic or other residues.</li> <li>Also, please refer to EVOLUTION 3E user manual (#65051, 65060) for operations prior use. The operator should be familiar with the functioning of the EVOLUTION 3E control unit (#19380).</li> <li>If you have additional questions, consult the MORIA website or your MORIA distributor.</li> </ul>  |                 |
| 2     | Connect the motor to the control unit, select the Epi-K™ and check vacuum prior use  | <ul style="list-style-type: none"> <li>See EVOLUTION 3E console user manual (#65051, 65060).</li> </ul>   |                 |
| 3     | Insert the head in the rails of the suction ring   | <ul style="list-style-type: none"> <li>When inserting the head in the rails of the suction ring, the handpiece must be at a 45° angle from the threaded shaft of the ring, to avoid any blade damage</li> <li>If the bottom of the head (the area where the separator edge protrudes) is in contact with the threaded shaft of the suction ring, inadvertent damage to the separator may occur. Release vacuum and change to a completely new set to avoid the possibility of irregular cut.</li> <li>Once the head properly inserted, slightly reduce the angle from 45° to 0°</li> <li>Carefully engage the head into the ring<br/>The shaft of the ring must be aligned with the motor.</li> </ul> | G, H            |
| 4     | Check for epikeratome translation prior to use (forward and backward pass)   | <ul style="list-style-type: none"> <li>Press vacuum ON, make forward then backward movements using footpedals, then press vacuum OFF. Please refer to EVOLUTION 3E user manual (#65051, 65060)</li> <li>Do not use the Epi-K™ system if the forward pass is not smooth, regular and uninterrupted.</li> </ul>   |                 |
| 5     | Position the suction ring on the eye. Activate vacuum by pressing the "Vacuum" footswitch once.  | <ul style="list-style-type: none"> <li>Verify that the device will not interfere with any other medical devices used during the procedure.</li> </ul>   |                 |
| 6     | Check the intraocular pressure with the tonometer  | <ul style="list-style-type: none"> <li>The tonometer (#19042) must be perfectly dry and used only on dry eyes. If the pressure is below 65 mm Hg, do not proceed with the surgery.</li> </ul>   |                 |
| 7     | <b>Lubricate the suction ring, the head and the separator</b>  | <ul style="list-style-type: none"> <li><b>Prior to use and during procedure, lubricate the rails of suction ring, the, head and separator with balanced salt solution or another appropriate solution.</b></li> <li>The device should not be placed in contact with solutions for which the potential interactions are unknown.</li> </ul>  |                 |
| 8     | Once the Epi-K™ system is in position, activate the motor by pressing the "Forward" footswitch (Speed 1, then Speed 2, then Speed 3). When the microkeratome head touches the stop, immediately release the "Forward" footswitch. Reverse by pressing the "Backward" footswitch. | <ul style="list-style-type: none"> <li>See EVOLUTION 3E console user manual (#65051, 65060).</li> <li>Hold the suction ring by its handle; check that your fingers do not hinder the forward movement of the device.</li> <li>Check nothing prevents or modifies head movement in the rails of the suction ring</li> <li>Ensure that there are no obstructions (speculum, eyelashes, eyelids, etc.) in its path</li> </ul>  |                 |
| 9     | Release the "Vacuum" footswitch by pressing once. Remove the suction ring.   | <ul style="list-style-type: none"> <li>Discard the disposable head, suction ring and aspiration tubing after each procedure in an appropriate container.</li> </ul>   |                 |

## **8 TROUBLESHOOTING**

For any further information and description, please refer to EVOLUTION 3E console user manual (#65051, 65060).

In case of defective handpiece, contact your MORIA representative or distributor for repairing.

## **9 CARE & MAINTENANCE**

### **9.1 INTRODUCTION**

In case of unusual vibrations or noises, do not use the unit and contact your distributor.

Please contact MORIA for any other servicing.

The following are only recommendations. They must be adapted according to the laws in the country of utilization. All the products and solutions for cleaning and disinfection must be used according to the manufacturer's instructions.

The recommendations for cleaning, disinfection and sterilization are updated and available on MORIA's website: <http://www.moria-surgical.com>.

MORIA recommends that you:

- avoid touching potentially contaminated areas,
- use gloves when performing cleaning and decontamination tasks.

### **9.2 INITIAL TREATMENT AND STORAGE**

To avoid risks of condensation inside the packaging, unpack and store the devices in a clean, dry environment. Do not store the instruments in an environment or next to other products that may possibly have a corrosive or magnetic effect. Avoid any contact between devices, especially those comprised of different materials. The device should be damage free and have no scratches or other surface defects.

Fragile devices or those requiring particular manipulation must be handled separately, taking special care to protect delicate parts.

### 9.3 CLEANING- DISINFECTION-STORAGE OF THE HANDPIECE AND ELECTRICAL CABLE

| Steps | What to do   |  <b>IMPORTANT WARNINGS</b>   | Related picture |
|-------|--|---|-----------------|
| 1     | <b>Cleaning</b> <ul style="list-style-type: none"> <li>• To clean the screw core of the motor and metallic parts, use the disposable brush (#19149) moistened with cleaning solutions.</li> <li>• Dip the blade oscillation shaft of the motor in a sterile distilled water bath and run the motor back and forth 10 times.</li> <li>• Remove the motor from the cup, unplug it from the console, and wipe the oscillation shaft of the motor with a lint-free cloth.</li> <li>• Cable can be cleaned with a lint-free cloth moistened with cleaning solutions.</li> <li>• Air-dry the shaft of the motor with clean, filtered compressed air (medical compressed air).</li> </ul> | <ul style="list-style-type: none"> <li>• <b>The handpiece must be carefully cleaned after each surgical procedure.</b></li> <li>• Do not use abrasives or scrapers to clean any elements of the system. Doing so may result in reduced precision, burrs and/or irregular tissue resection.</li> <li>• Do not immerse the motor in the sterile distilled water bath.</li> <li>• To prevent damaging the connectors and to ensure proper functioning of the motor, never pull on the cables and never hold the motor by its cable.</li> </ul> | I               |
| 2     | <b>Disinfecting &amp; Drying</b> <ul style="list-style-type: none"> <li>• Wipe each part of the motor with a lint-free cloth moistened with disinfecting solutions.</li> <li>• Use appropriate disinfecting solutions (spray or wipe clean disinfectant) according to the manufacturers' instructions.</li> <li>• Dry carefully with disposable instrument cleaner (lint-free) then with clean, filtered, compressed air.</li> </ul>   | <ul style="list-style-type: none"> <li>• It is imperative to verify that the devices are perfectly clean and free of any organic or other residues.</li> <li>• Do not gas sterilize (ETO) the motor.</li> <li>• Do not autoclave the motor.</li> <li>• No current sterilisation process is compatible with MORIA motors.</li> </ul>   |                 |
| 3     | <b>Storage</b>   | <ul style="list-style-type: none"> <li>• Failure to completely dry the inside of the motor can result in oxidation.</li> <li>• Do not store products that are not completely dry; doing so may create rust and result in irregular corneal resection.</li> <li>• When not in use, the Epi-K™ system should be kept in its storage box (#22519512) in a dry atmosphere.</li> </ul>   |                 |

## **10 WARRANTY**

### **10.1 SCOPE OF WARRANTY**

| <b>Designation</b>                        | <b>MORIA references</b> |
|---|-------------------------|
| Epi-K™ handpiece (green color)            | 19342                   |
| Wrench                                    | 19345C                  |
| Storage box for the Epi-K™ handpiece      | 22519512                |
| EVOLUTION 3E console (except the battery) | 19380                   |
| EVOLUTION 3E footswitch Epi-K™            | 19381                   |

- The above items as well as spare parts and labor necessary for their repair are covered by warranty. Any items returned must be sent in their original packaging, after having previously been disinfected.
- The maintenance operations and the replacement of spare parts will be exclusively carried out by technicians authorized by MORIA.
- Any defective part exchanged during the period of warranty becomes property of MORIA.
- Instruments and accessories that cannot be reused are excluded from this warranty.
- Future upgrades and/or improvements on the keratome are not implied by this warranty.

### **10.2 NON-APPLICATION OF WARRANTY**

The warranty will not be applicable under any of the following conditions:

- Defects or malfunction that occur out of the warranty period (10.3).
- Normal wear and tear.
- Negligence or usage that does not comply with the specifications in the user manual.
- The use of supplies, spare parts, or accessories other than those supplied by MORIA.
- Any disassembly, modification or intervention carried out on the devices by a person not authorized by MORIA.

### **10.3 WARRANTY PERIOD**

- The warranty takes effect on the date the material is dispatched.
- The duration of the warranty is 12 months from the date of effect.

### **10.4 LIABILITY**

- The liability of MORIA is limited to the supply of the services mentioned in paragraph 10.1. MORIA will not be held responsible for any direct or indirect damage suffered by the client owing to the interventions within the scope of this warranty.
- For any dispute concerning the interpretation or the execution of the present contract or the present general terms and conditions, the Commercial Court of Nanterre (France) will have sole jurisdiction.

## **11 DRAWINGS**

### **A. Disposable head**

- 1 Threaded shaft
- 2 ONE EPI
- 3 Separator housing
- 4 Slide
- 5 Separator
- 6 Aplanating plate

### **B. Suction ring**

- 1 Vacuum port
- 2 Threaded shaft
- 3 Guiding
- 4 Stopper with engraved stop value
- 5 Handle

### **C. Setting the stopper**

- 1 Spanner
- 2 Stopper

### **D. Connecting the suction ring to the aspiration tubing**

- 1 Aspiration tubing end
- 2 Tubing of the suction ring

### **E. Motor**

- 1 Motor housing
- 2 Motor cord
- 3 Threaded ring
- 4 Oscillation motor shaft
- 5 Head alignment guide
- 6 Advance drive

### **F. Mounting the head on the motor**

- 1 Open the pack containing the disposable head and suction ring
- 2 Screw the head clockwise onto the threaded motor shaft.
- 3 Check that the head is screwed on tightly to the motor. There should be no rotation of the head.

### **G. Mounting the head on the suction ring**

- 1 When inserting the head in the rails of the suction ring, the handpiece must be at a 45° angle from the threaded shaft of the ring
- 2 Once the head properly inserted, slightly reduce the angle from 45° to 0°.

### **H. Alignment of the head, motor and suction ring**

- 1 Proper assembly
- 2 Incorrect alignment
- 3 Alignment of the threaded shaft with the motor

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