



## Epi-LASIK retaining the epithelial flap

### Ring Selection

- For myopic treatments on normal or small eyes (corneal diameter up to 11.5 mm), select the **-1, 0 or +1** suction ring. Use of a higher ring value will decrease diameter of the epithelial resection by 0.1 mm.
- For hyperopic treatments, treatments requiring a large ablation zone (wavefront-guided, large pupils), or large eyes (corneal diameter > 11.5 mm), select the **LD** (Large Diameter) suction ring.

### Stop Selection for either ring

- Select the **7.5** stop for average  $K < 43.5$  D
- Select the **8.0** stop for average  $K \geq 43.5$  D

## Epi-LASIK discarding the epithelial flap

### Ring and Stop Selection

- For myopic treatments on normal or small eyes (corneal diameter up to 11.5 mm), select the **-1, 0 or +1** suction ring with the **9.0** stop for all cases decentering nasally in order to get a free cap. To decenter correctly, allow 0.5 mm of sclera to show through the ring. Use of a higher ring value will decrease diameter of the epithelial resection by 0.1 mm.
- For hyperopic treatments, treatments requiring a large ablation zone, or large eyes (corneal diameter > 11.5 mm), select the **LD** (Large Diameter) suction ring with the **9.5** stop.

Performance may deviate from the guidelines depending on several patient-related factors (IOP, degree of myopia, age, gender, ethnicity, keratometry, etc.) and surgery-related factors (corneal hydration, suction time, etc.). Guidelines are a general guide only. It is strongly recommended that every surgeon establish his own nomogram. Moria shall not be responsible for any direct, incidental, consequential or exemplary damage suffered by any party, even if that party has not been advised of the possibility of such damage.

To obtain more information

**MORIA SA**  
27, rue du Pied de Fourche  
03160 Bourbon Archambault  
FRANCE  
Phone: +33 (0) 1 46 74 46 74  
Fax: +33 (0) 1 46 74 46 70  
[www.moria-surgical.com](http://www.moria-surgical.com)

**MORIA JAPAN KK**  
Inami Bldg. 6F  
3-24-2 Hongo Bunkyo-ku  
Tokyo 113-0033  
JAPAN  
Phone: +81 3 5842 6761  
Fax: +81 3 5842 6762  
[www.moriajapan.com](http://www.moriajapan.com)

**Moria Shanghai Office**  
RM4, 10 FL, Kai Li Building  
N°432 Huai Hai West Road  
200052 Shanghai  
CHINA  
Phone: +86 21 5258 5068  
Fax: +86 21 5258 5067  
[www.moria-surgical.com.cn](http://www.moria-surgical.com.cn)

**Moria**



## Epi-LASIK retaining the epithelial flap

For myopic treatments on normal or small eyes (corneal diameter up to 11.5 mm), select the **-1** or **0** suction ring with the following stops:

- **7.5** stop for average K < 43.5 D
- **8.0** stop for average K ≥ 43.5 D

Use of a "0" ring value will decrease diameter of the epithelial resection by 0.1 mm compared to use of a "-1" ring on same eye.

For large eyes (corneal diameter > 11.5 mm), select the **-1** suction ring with the **7.5** stop.

## Epi-LASIK discarding the epithelial flap

Select the **-1** or **0** suction ring with the **9.5** stop for all cases decentering nasally in order to get a free cap.

Use of a "0" ring value will decrease diameter of the epithelial resection by 0.1 mm compared to use of a "-1" ring on same eye.

Performance may deviate from the guidelines depending on several patient-related factors (IOP, degree of myopia, age, gender, ethnicity, keratometry, etc.) and surgery-related factors (corneal hydration, suction time, etc.). Guidelines are a general guide only. It is strongly recommended that every surgeon establish his own nomogram. Moria shall not be responsible for any direct, incidental, consequential or exemplary damage suffered by any party, even if that party has not been advised of the possibility of such damage.

To obtain more information



**MORIA SA**  
27, rue du Pied de Fourche  
03160 Bourbon Archambault  
FRANCE  
Phone: +33 (0) 1 46 74 46 74  
Fax: +33 (0) 1 46 74 46 70  
[www.moria-surgical.com](http://www.moria-surgical.com)

**MORIA JAPAN KK**  
Inami Bldg. 6F  
3-24-2 Hongo Bunkyo-ku  
Tokyo 113-0033  
JAPAN  
Phone: + 81 3 5842 6761  
Fax: + 81 3 5842 6762  
[www.moriajapan.com](http://www.moriajapan.com)

**Moria Shanghai Office**  
RM4, 10 FL, Kai Li Building  
N°432 Huai Hai West Road  
200052 Shanghai  
CHINA  
Phone: +86 21 5258 5068  
Fax: +86 21 5258 5067  
[www.moria-surgical.com.cn](http://www.moria-surgical.com.cn)

# Moria