EN - ANNEXE - USER MANUAL EVOLUTION 3E	3
FR - ANNEXE - MANUEL UTILISATEUR EVOLUTION 3E	7
ES - ANEXO – INSTRUCCIONES DE USO EVOLUTION 3E.....	11
IT - ALLEGATO - MANUALE UTENTE EVOLUTION 3E	15
DE - ANHANG – BENUTZERHANDBUCH EVOLUTION 3 E.....	19
NL - BIJLAGE - HANDLEIDING EVOLUTION 3E	23
PT - ANEXO – MANUAL DO UTILIZADOR EVOLUTION 3E.....	27
SV - ANNEX – BRUKSANVISNING FÖR EVOLUTION 3E.....	31
TR - EK – EVOLUTION 3E KULLANIM KILAVUZU.....	35

Ref.: 65073

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Document Aletiq : Rév. C-05.2023 (Validée)

ANNEXE - USER MANUAL EVOLUTION 3E

Guidance and manufacturer's declaration :

**electromagnetic emissions
electromagnetic immunity
recommended separation distances**

WARNING AGAINST ELECTROMAGNETIC COMPATIBILITY (EMC)

The EVOLUTION 3E needs special precautions regarding EMC.

The EVOLUTION 3E has to be installed and put into service according to the EMC requirements.

The EVOLUTION 3E can be affected by portable and mobile RF communications.

WARNING AGAINST ELECTROSTATIC DISCHARGES (ESD)

Pins of connectors identified with the ESD warning symbol should not be touched and connections of these connectors should not be made unless ESD precautionary procedures are used.

When a cable is connected to a connector identified with the ESD warning symbol, the EVOLUTION 3E shall be switch off. Then, touch the main earth terminal or metal part which is connected the earth in order to unload your body and connected the wished cable. Follow the starting procedure (refer to user manual #65060).

All staff involved must receive an explanation of the ESD warning symbol and a training in ESD precautionary procedure.

The training's content must include :

- the introduction of ESD warning
- a practical example of cable connection when ESD warning is present
- the insurance that the staff has completely understood the procedure and the risk if the staff doesn't apply it.

WARNING AGAINST INTERFERENCE RISK

The use of accessories, transducers and cables other than those specified, with the exception of transducers and cables sold by the manufacturer as replacement parts for internal components, may result in increased emissions or decreased immunity EVOLUTION 3E.

The EVOLUTION 3E should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, the EVOLUTION 3E should be observed to verify normal operation in the configuration in which it will be used.

The use of the accessory, transducer or cable with EVOLUTION 3E other than those specified may result in increased emissions or decreased immunity of the EVOLUTION 3E.

Table 1 - Guidance and manufacturer's declaration – electromagnetic emissions – for all ME equipment and ME systems

Guidance and manufacturer's declaration – electromagnetic emissions		
EVOLUTION 3E is intended for use in the electromagnetic environment specified below. The customer or the user of EVOLUTION 3E product should assure that it is used in such an electromagnetic environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The EVOLUTION 3E uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	The EVOLUTION 3E is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Table 2 - Guidance and manufacturer's declaration – electromagnetic immunity – for all ME equipment and ME systems

Guidance and manufacturer's declaration – electromagnetic immunity			
The EVOLUTION 3E is intended for use in the electromagnetic environment specified below. The customer or the user of the EVOLUTION 3E should assure that it is used in such an electromagnetic environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	±1 kV mode differential ±2 kV mode common	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 s	0 Vac during 10 ms 92 Vac during 100 ms 161 Vac during 1 s 0 Vac during 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the EVOLUTION 3E requires continued operation during power mains interruptions, it is recommended that the EVOLUTION 3E be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m for 50Hz and for 60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE UT is the a.c. mains voltage prior to application of the test level

Table 3 - Guidance and manufacturer's declaration – electromagnetic immunity – for ME equipment and ME systems that are not life-supporting

Guidance and manufacturer's declaration – electromagnetic immunity			
The EVOLUTION 3E is intended for use in the electromagnetic environment specified below. The customer or the user of the EVOLUTION 3E should assure that it is used in such an electromagnetic environment.			
Immunity test	Immunity test	Immunity test	Immunity test
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of the EVOLUTION 3E, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 1.17 \sqrt{P}$
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2,5 GHz	3 V/m	$d = 1.17 \sqrt{P} \text{ 80 MHz to 800 MHz}$ $d = 2.33 \sqrt{P} \text{ 800 MHz to 2,5 GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey (1) should be less than the compliance level in each frequency range. (2)</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
(1) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the EVOLUTION 3E is used exceeds the applicable RF compliance level above, the EVOLUTION 3E should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the EVOLUTION 3E.			
(2) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Table 4 – Recommended separation distances between portable and mobile RF communications equipment and the ME equipment or ME system – for ME equipment and ME systems that are not life-supporting

Recommended separation distances between portable and mobile RF communications equipment and the EVOLUTION 3E			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.17 \sqrt{P}$	80 MHz to 800 MHz $d = 1.17 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2.33 \sqrt{P}$
0,01	0.12	0.12	0.23
0,1	0.37	0.37	0.74
1	1.16	1.16	2.33
10	3.69	3.69	7.37
100	11.6	11.6	23.3

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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MORIA S.A. – 27, rue du Pied de Fourche - 03160 BOURBON L'ARCHAMBAULT – France