# #65085

EN

## Instructions for use for the ONE® Single-Use Vacuum Trephine # 17201Dxxx

### Indication

Moria Single-Use trephines are devices intended to incise a circular disc with partial or full thickness on the patient's cornea or on donor corneal tissue.

### Device components

The device consists of a trephine kit (Fig.1) and a suction unit (Fig.2).

The most recent version of this document is available on Moria website (<http://www.moria-surgical.com>) via the restricted access heading «Documentation/Library». Translations of this document are also available in other languages under that same heading. Registration required.

### Warning

- These instructions explain only how to use this medical device and do not describe the steps necessary to perform a corneal graft. The surgical procedure remains the responsibility of the surgeon experienced in corneal grafts.
- The trephine is sterile and ready for use. Do not use it if the pouch in which it is packaged is damaged.
- Do not dismantle the trephine, as this may damage the device.
- The blade is sharp. Do not touch it.
- The trephine is for single-use only: it is not designed to be resterilized or reused under any circumstances.
- Single Use instruments are indicated for use with restriction (Child, pregnant and/or breast-feeding are not eligible for Ophthalmic surgery without doctor's consent)
- At-risk populations (pediatric population, pregnant and breastfeeding women etc...) eligibility is subject to surgeon's judgment and experience.

### Directions

#### 1/ Select the appropriate diameter and assemble the device

The diameter is indicated on the packaging and on the trephine (Fig.3).

- Connect the trephine tubing (E, F) to the syringe (Fig.4) before use.
- The blade is delivered retracted at 250µm above the cornea (i.e. one full turn).

#### 2/ Check before using

- Test that the blade (A) descends fully by turning the collar clockwise (Fig. 5) until it meets the depth stop, and check that the movement is smooth. Check under the microscope (minimum 20x) that the cutting edge of the blade has no visible defect. Defects can compromise the performance of the cut.
- Retract the blade fully in the trephine before placing it on the patient's eye.

#### 3/ To use the device

**Center the trephine on the cornea:** Mark the center of the cornea, then position the trephine using the centration indicator (C). Check that the device is clear of interference from other devices such as speculums. Check that the cornea is fully flattened under the centration indicator. In case of severe keratoconus with marked protrusion of the cone, the trephination can be improved by flattening the apex of the cone.

**Apply suction:** Depress the syringe piston completely (D), then press the trephine gently and uniformly onto the cornea. Slowly release the syringe piston.

**Confirm suction:** Pull upward on the trephine very gently. It should remain firmly attached to the cornea.

##### Trephining:

- Hold the trephine between thumb and index finger on the grip area (H) (Fig. 6).
- Slowly turn the collar (G) clockwise (Fig. 5) to cause the blade to descend.
- If desired, the reference on the collar can be used to count the number of turns made (Fig. 7). Each complete turn corresponds to a vertical blade descent of 250µm. Pachymetry performed by the surgeon gives an approximate value and should be taken into account when the blade descends (Fig 8).

**Removing the trephine:** Depress the syringe piston to release suction and remove the trephine vertically.

The base of the trephine in contact with the cornea has 16 equidistant marks (B) that can be used for the remainder of the procedure.

At the end of the surgical procedure, retract the blade (A) and place the trephine and its accessories in an appropriate container for destruction (in compliance with applicable regulations).

Explanation of symbols (General list of symbols used by Moria)

|                                                                                              |                                                                                                         |                                                                                                      |
|----------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------|
|  XXXXXX     | CATALOGUE REFERENCE                                                                                     |                                                                                                      |
|             | QUANTITY                                                                                                |                                                                                                      |
|  YYYY/MM/DD | USE BY                                                                                                  |                                                                                                      |
|  XXXXXX     | BATCH CODE                                                                                              |                                                                                                      |
|  EO         | STERILIZED USING ETHYLENE OXIDE                                                                         |                                                                                                      |
|             | DO NOT REUSE                                                                                            |  DO NOT RESTERILIZE |
|  XXXXXXXX   | MANUFACTURER                                                                                            |                                                                                                      |
|  YYYY/MM/DD | DATE OF MANUFACTURE                                                                                     |                                                                                                      |
|  MD         | MEDICAL DEVICE                                                                                          |                                                                                                      |
|             | CONSULT OPERATING INSTRUCTIONS                                                                          |                                                                                                      |
|             | DO NOT USE IF PACKAGE IS DAMAGED                                                                        |                                                                                                      |
|             | KEEP DRY                                                                                                |                                                                                                      |
|  only       | CAUTION for USA only: US federal laws restrict this device to sale by, or on the order of, a physician. |                                                                                                      |

List of references

|                      |                              |
|----------------------|------------------------------|
| Reference: 17201Dxxx |                              |
| 17201Dxxx            | VACUUM TREPHINE Dxxx         |
| Packaging            | 1 per box/Sterile/Single-Use |

Accessories

|             |                                                                                                                                      |
|-------------|--------------------------------------------------------------------------------------------------------------------------------------|
| Reference   | 17200Dxxx                                                                                                                            |
| Description | ONE® – Disposable corneal vacuum punch<br>Diameters: 6.50 to 9.00 mm with 0.25 mm increments. Also 6.00 and 9.50 mm diameters        |
| Packaging   | 1 per box/Sterile/Single-Use                                                                                                         |
| Reference   | 17204                                                                                                                                |
| Description | ONE® – Artificial Chamber for Single-Use Trephine<br>This trephine is compatible only with the single-use artificial chamber #17204. |
| Packaging   | 1 per box/Sterile/Single-Use                                                                                                         |

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**MORIA Inc**