EN - ACP SYSTEM	
INSTRUCTION MANUAL	1-11
FR - SYSTÈME ACP	
MANUEL UTILISATEUR	12-23
ES - SISTEMA DE ACP	
MANUAL DEL USUARIO.....	24-35
IT - SISTEMA ACP	
MANUALE D'USO	36-47
DE - ACP-SYSTEM	
BENUTZERHANDBUCH	48-58
NL - ACP-SYSTEEM	
HANDLEIDING	59-70
PT - SISTEMA ACP	
MANUAL DO UTILIZADOR	71-82
SV - ACP-SYSTEM	
BRUKSANVISNING	83-93
TR - ACP SİSTEMİ	
KULLANICI EL KİTABI.....	94-105

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Moria

Contents

I.	DISCLAIMER	2
II.	GUIDANCE AND MANUFACTURER’S DECLARATION: ELECTROMAGNETIC EMISSIONS AND IMMUNITY	2
III.	EQUIPMENT AND ACCESSORIES LIST	2
A.	EQUIPMENT LIST	2
B.	ACCESSORIES LIST	2
C.	USER INFORMATION	2
IV.	TECHNICAL SPECIFICATIONS.....	3
V.	REGULATORY INFORMATION	3
VI.	INTENDED USE	3
VII.	LABELING INFORMATION	4
VIII.	INSTALLATION AND CONNECTION	5
A.	RECOMMENDATIONS.....	5
B.	ELECTRICAL SETTING.....	5
C.	CHOOSING THE PRESSURE.....	5
D.	CONNECTING ACP TUBING SET.....	5
IX.	OPERATION	6
A.	MOUNT THE CORNEAL TISSUE ON THE ARTIFICIAL CHAMBER.....	6
B.	ACTIVATE PRESSURE.....	6
C.	MODES OF USE	6
	1. “OPEN SYSTEM” = PRESSURE DIRECTLY COMING FROM ACP SYSTEM TO THE ARTIFICIAL CHAMBER	6
	2. “CLOSED SYSTEM” = PRESSURE PREVIOUSLY ESTABLISHED BY ACP SYSTEM	7
D.	REMOVE THE CORNEAL TISSUE AFTER OPERATION.....	7
	1. “OPEN SYSTEM”	7
	2. “CLOSED SYSTEM”	7
E.	DISCONNECT ACP TUBING SET.....	8
F.	END OF PROCEDURE	8
X.	COMPATIBILITY FOOTSWITCH – ACP SYSTEM	8
XI.	TROUBLE SHOOTING.....	8
A.	ACP SYSTEM.....	8
B.	DIGITAL DISPLAY	8
C.	PRESSURE SETTING	8
D.	ABSENCE OF AIR DELIVERY.....	9
E.	LEAKAGE AFTER CORNEAL TISSUE SETTING ON THE ARTIFICIAL CHAMBER.....	9
XII.	CARE & MAINTENANCE.....	9
A.	CARE.....	9
B.	MAINTENANCE.....	9
	1.REPLACING FUSE	9
XIII.	WARRANTY.....	10
A.	SCOPE OF WARRANTY.....	10
B.	NON-APPLICATION OF WARRANTY	10
C.	WARRANTY PERIOD	10
D.	LIABILITY.....	10
XIV.	DRAWINGS	11
A.	ACP SYSTEM (FRONT).....	11
B.	ACP SYSTEM (BACK)	11
C.	FOOTSWITCH	11
D.	POWER SUPPLY	11
E.	CHANGING FUSE	11

The most recent version of this user guide is available on MORIA website: <http://www.moria-surgical.com>.

I. DISCLAIMER

ACP system has been designed for optimal operation, provided that the recommendations listed in this user manual are followed carefully. If, for any reasons, ACP system does not perform properly, have it checked immediately by MORIA. MORIA strongly recommends having the system thoroughly inspected by MORIA on a routine basis every year.

The use of materials and/or components of a brand other than MORIA with ACP system will immediately nullify the MORIA warranty. MORIA may not be held responsible for any damage resulting from the use of materials and/or components of a brand other than MORIA.

II. GUIDANCE AND MANUFACTURER'S DECLARATION: ELECTROMAGNETIC EMISSIONS AND IMMUNITY

Refer to annex document (#65114).

III. EQUIPMENT AND ACCESSORIES LIST

A. EQUIPMENT LIST

Designation	MORIA reference
ACP pack	G0
ACP system	19188
ACP footswitch	19189
ACP power supply	19190
Supply cord EEC	19450
Supply cord USA	19451
Supply cord UK	19452
Supply cord China	19453
Supply cord Brazil	19454
ACP carrying case	19191

B. ACCESSORIES LIST

Designation	MORIA reference
ACP tubing set x10	19192
Artificial chamber	-

C. USER INFORMATION

Designation	MORIA reference
ACP system user manual	65112
ACP system user manual (XX=other)(MORIA website)	65112XX
Annex "Guidance and manufacturer's declaration: electromagnetic emissions and immunity"	65114
Annex "Guidance and manufacturer's declaration: electromagnetic emissions and immunity" (XX=other)(MORIA website)	65114XX
ACP tubing set instruction for use	65113
ACP tubing set instruction for use (XX=other)(MORIA website)	65113XX

IV. TECHNICAL SPECIFICATIONS

Reference	mm	kg	°C	HR	Pressure	Elect. Supply
19188	235 x 165 x 85	1,54	15°C to 35°C 59°F to 95°F	45% to 75%	Pressure range : From 30 ± 10 mmHg to 250 mmHg Overpressure: 250 mmHg Precision of displayed values : ± 5 mmHg	DC 15V 24W
19189	98 x 66 x 28	0,34	15°C to 35°C 59°F to 95°F	45% to 75%		
19190	139,5 x 62,3 x 33,8 (case) 1200 ± 30 for the cord	0,36	15°C to 35°C 59°F to 95°F	45% to 75%		AC input : 100-240V AC 47-63 Hz DC output 15V 3,67A

The operating atmospheric pressure must be in between 70 and 106 kPa (525 and 795 mmHg).

V. REGULATORY INFORMATION

	MORIA S.A. 27, rue du Pied de Fourche - 03160 Bourbon l'Archambault - France Phone +33 (0) 4 70 67 09 04 Fax +33 (0) 4 70 67 06 61 moria@moria-int.com http://www.moria-surgical.com
Electrical Safety Standard, electromagnetic compatibility	IEC 60601-1-1 Class I IEC 60601-1-2
Annex "Guidance and manufacturer's declaration: electromagnetic emissions and immunity"	Refer to annex document (#65114).
EUROPE	CE according to MDD 93/42/CEE
USA	Product registered at the Food and Drug Administration (FDA). Caution USA only: Federal law restricts the use of this device to physicians or licensed practitioners.
Customer Service Information	Contact your local dealer or MORIA.
	For European Union customers only: this symbol indicates that within the European Union, the product must be discarded in a separate collection bin at the end of its useful life. This applies not only to this device, but also all accessories, including the footswitch and power supply, regardless of whether those accessories are marked with the symbol. Do not dispose as unsorted municipal waste. For users outside of the European Union: please refer to local environmental regulatory regarding waste of electrical and electronic equipment.

VI. INTENDED USE

ACP system is intended to pressurize an artificial chamber at a desired value with a tissue previously mounted.

Note: for use with cadaver tissue only.

It is recommended the ACP system be used only by experienced operators with specific knowledge in corneal preparation.

VII. LABELING INFORMATION

 XXXXXXXX	CATALOGUE NUMBER OR SKU
 XXXXXXXX	SERIAL NUMBER
 XXXXXXXX	DATE OF MANUFACTURE
	MANUFACTURER
	CAUTION : CONSULT ACCOMPANYING DOCUMENT(S)
	CONSULT OPERATING INSTRUCTIONS
	DISCARD IN A SEPARATE COLLECTION BIN
	WARNING AGAINST ELECTROSTATIC DISCHARGES
	WARNING AGAINST INTERFERENCE RISK
	REFER TO USER MANUAL
	SIGN OF GENERAL SAFETY
	ALTERNATIVE CURRENT
	DIRECT CURRENT
IP X6	FOOTSWITCH PROTECTED AGAINST POWERFUL WATER JETS
	THIS SIDE UP
	FRAGILE
	KEEP DRY
	STORAGE CONDITIONS: <ul style="list-style-type: none"> • TEMPERATURE: 15°C – 35°C / 59°F – 95°F • HUMIDITY RATE: 45% – 75%

VIII. INSTALLATION AND CONNECTION

A. RECOMMENDATIONS

Following is a list of general recommendations when installing ACP system:

- Make sure ACP system is installed no more than 2 meters (78.7 inches) away from user.
- ACP system must remain in front of user during operation to visualize every indication (value, warning).
- User must be able to hear signals from ACP system properly.
- User must be cautious about ACP system displacement due to its connections.
- ACP system must not be installed above or close to a radiating device.

B. ELECTRICAL SETTING

Steps	What to do	 IMPORTANT WARNINGS	Related pictures
1	Using the screw connector, plug in the ACP footswitch (#19189) to the rear panel of ACP system (#19188).	Finger tight only.	B5-C
2	Connect ACP power supply (#19190) to the supply cord (#19450-19451-19452-19453-19454). Connect ACP power supply to the rear panel of ACP system. Connect ACP supply cord to the power socket.	Visually check integrity of ACP power supply and ACP supply cord prior use. To change fuse, see chapter XII. A. To avoid electrical shock, ACP system must only be connected to a grounded power system. Always allow an easy access to the electric supply and power socket. Position ACP system so that displays, connections and buttons remain visible and accessible.	B3-D
3	Switch on ACP system turning button "O/I" to "I". Warm-up 5 min before use.	Visually check the yellow LED is ON and fixed. The display indicates the pressure in mmHg.	B2-A3-A5

C. CHOOSING THE PRESSURE

Steps	What to do	 IMPORTANT WARNINGS	Related pictures
4	Turn the knob to choose the desired pressure value on the digital display.	Visually check that : <ul style="list-style-type: none"> • the display indicates the pressure in real time. • the value increases or decreases in function of the sense of the turning knob. • the yellow LED is still ON and fixed. MORIA recommends to gradually increasing the pressure.	A3-A4-A5

D. CONNECTING ACP TUBING SET

Steps	What to do	 IMPORTANT WARNINGS	Related pictures
5	Visually check the integrity of ACP tubing set (#19192) before set up.	MORIA strongly advises against the reuse of ACP tubing set as it may damage connector.	
6	Connect ACP tubing set to ACP system and to the artificial chamber.	Refer to instructions for use of ACP tubing set (#65113). Poorly connected tubing, either on ACP system or on the artificial chamber, may also cause insufficient pressure. It is important to check the integrity of all connections. Warning: <ul style="list-style-type: none"> • do not reverse the position of ACP tubing set. • if liquid appears in ACP tubing set (excluding infusion part), replace it immediately. • visually check integrity of ACP tubing set after set up. 	A2

IX. OPERATION

A. MOUNT THE CORNEAL TISSUE ON THE ARTIFICIAL CHAMBER

Steps	What to do	 IMPORTANT WARNINGS	Related pictures
7	Fill the syringe coming from ACP tubing set with saline salt solution or corneal storage medium.	Refer to instructions for use of ACP tubing set (#65113). Check there is no air bubble in the syringe.	
8	Connect the syringe to the stopcock of ACP tubing set.		
9	Make sure stopcock is open to the syringe and artificial chamber, while being closed to ACP system.		
10	Press the syringe to infuse the artificial chamber.	Make sure there is no air bubble in the tubing.	
11	Mount the corneal tissue on the artificial chamber.	Make sure there is no air bubble below the corneal tissue.	

B. ACTIVATE PRESSURE

Steps	What to do	 IMPORTANT WARNINGS	Related pictures
12	Push on the footswitch to open the output valve and to release air in ACP tubing set.	Visually check that: <ul style="list-style-type: none"> • the yellow LED is OFF. • the green LED is ON. 	A5-A6-C

C. MODES OF USE

1. "OPEN SYSTEM" = PRESSURE DIRECTLY COMING FROM ACP SYSTEM TO THE ARTIFICIAL CHAMBER

Steps	What to do	 IMPORTANT WARNINGS	Related pictures
13-1a	Before proceeding to the operation, check the pressure display.	In case pressure display is different from desired value, readjust the pressure value on the digital display by turning the knob. Note: ACP system can deviate during a few minutes of ± 15 mmHg. After warm-up time, ACP system becomes stable.	A3-A4
13-1b	Proceed to the operation: shift stopcock of ACP tubing set to be open to both ACP system and artificial chamber, while being closed to the syringe.		
13-1c	Pressure is established in the artificial chamber and applied directly from ACP system. Corneal tissue is ready for operation.	Visually check that the displayed value is similar to the desired one. If not, readjust the pressure value on the digital display by turning the knob.	A3-A4

Note: In open system, the pressure displayed on ACP system is the same as established in the artificial chamber.

2. “CLOSED SYSTEM” = PRESSURE PREVIOUSLY ESTABLISHED BY ACP SYSTEM

Steps	What to do	 IMPORTANT WARNINGS	Related pictures
13-2a	Before proceeding to the operation, check the pressure display.	In case pressure display is different from desired value, readjust the pressure value on the digital display by turning the knob.	A3-A4
13.2b	Establish pressure in the artificial chamber by shifting the stopcock of ACP tubing set to be open to both ACP system and artificial chamber, while being closed to the syringe.	Visually check that the displayed value is similar to the desired one. If not, readjust the pressure value on the digital display by turning the knob. Note: ACP system can deviate during a few minutes of ± 15 mmHg. After warm-up time, ACP system becomes stable.	A3-A4
13-2c	Close the system by shifting the stopcock of ACP tubing set to be open to both ACP system and syringe, while being closed to the artificial chamber.	Note: once the system is closed, syringe can gently be pushed away.	
13-2d	Corneal tissue is ready for operation.		

D. REMOVE THE CORNEAL TISSUE AFTER OPERATION

1. “OPEN SYSTEM”

Steps	What to do	 IMPORTANT WARNINGS	Related pictures
14-1a	Shift stopcock of ACP tubing set to be open to both artificial chamber and syringe, while being closed to ACP system.		
14-1b	Push on the footswitch to close the output valve.	Visually check that: <ul style="list-style-type: none"> • the yellow LED is ON and fixed. • the green LED is OFF. 	A5-A6-C
14-1c	Remove the corneal tissue.		

2. “CLOSED SYSTEM”

Steps	What to do	 IMPORTANT WARNINGS	Related pictures
14-2a	Proceed to the operation by shifting the stopcock of ACP tubing set to be open to both artificial chamber and syringe, while being closed to ACP system.		
14-2b	Push on the footswitch to close the output valve.	Visually check that: <ul style="list-style-type: none"> • the yellow LED is ON and fixed. • the green LED is OFF. 	A5-A6-C
14-2c	Remove the corneal tissue.		

E. DISCONNECT ACP TUBING SET

Steps	What to do	 IMPORTANT WARNINGS	Related pictures
15	Disconnect ACP tubing set.	Finger tight only. Discard ACP tubing set after each procedure in an appropriate container.	A2

F. END OF PROCEDURE

Steps	What to do	 IMPORTANT WARNINGS	Related pictures
16	Switch off ACP system turning button "O/I" to "O".	Visually check that: <ul style="list-style-type: none">• the yellow LED is OFF.• the green LED is OFF.• the display is OFF.	B2 A3-A5-A6

X. COMPATIBILITY FOOTSWITCH – ACP SYSTEM

ACP footswitch (#19189) must be used only with ACP system (#19188).

XI. TROUBLE SHOOTING

A. ACP SYSTEM

Message	Reason	Intervention
ACP system doesn't switch on	No power. Failure of the power supply, the supply cord connection or the fuse.	Check the connection of the main ACP power supply, the green LED of power supply must be ON. If not, unplug the system and check the integrity of the fuse and replace it if damaged. If unable to solve the problem: do not use ACP system and contact your distributor or MORIA to have it fixed.
ACP system switches off during operation	Loss of power. Failure of the power supply, the supply cord connection or the fuse.	Check the connection of the main ACP power supply, the green LED of power supply must be ON. If not, unplug the system and check the integrity of the fuse and replace it if damaged. If able to solve the problem; check if the corneal tissue has not been damaged and start again the different steps described in paragraph VIII. If unable to solve the problem: do not use ACP system and contact your distributor or MORIA to have it fixed.

B. DIGITAL DISPLAY

Message	Reason	Intervention
Same displayed value whatever the position of the turning knob or 888 display or lines display or loss of digital display	Failure of display card	Do not use ACP system and contact your distributor or MORIA to have it fixed.

C. PRESSURE SETTING

Message	Reason	Intervention
The pressure doesn't increase nor decrease with the turning knob	No connection of ACP tubing set or leakage Failure of internal components	Check the integrity of ACP tubing set. If leakage is present, change it. Check the sealing of the artificial chamber too. Do not use ACP system and contact your distributor or MORIA to have it fixed.

D. ABSENCE OF AIR DELIVERY

Message	Reason	Intervention
No air delivery while green LED is ON.	ACP tubing set leaking or not well connected. Output connector malfunctioning. Internal components malfunctioning.	Disconnect ACP tubing set and place your finger on the output connector to check whether air is delivered outside or not: <ul style="list-style-type: none"> • if yes: then change ACP tubing set. • if unable to solve the problem: do not use ACP system and contact your distributor or MORIA to have it checked.

E. LEAKAGE AFTER CORNEAL TISSUE SETTING ON THE ARTIFICIAL CHAMBER

Message	Reason	Intervention
Leakage when the corneal tissue is on the artificial chamber	ACP tubing set not properly connected or artificial chamber leaking. Output valve malfunctioning. Internal components malfunctioning.	Check if the corneal tissue has not been damaged after collapsing. Change tubing and start again the different steps described in paragraph VIII.

XII. CARE & MAINTENANCE

A. CARE

- To clean ACP system, wipe the system with a lint-free cloth moistened with disinfecting solutions. Do not use solvents, abrasives, fluids or solutions on ACP system, ACP footswitch and ACP power supply.
- Avoid any liquid projection on ACP system, ACP footswitch and ACP power supply.
- Do not open neither ACP system nor its accessories.
- In cases of emergency, disconnect the power socket or the power supply connector.
- In cases of unusual vibrations or noises, instability, unusual performance, do not use ACP system and contact MORIA or your distributor for fixing.
- In order to assess any variation of the system accuracy, it is recommended to check on a routine basis the obtained pressure inside the artificial chamber with a 65 and/or 90 mmHg Barraquer tonometer and ACP system set up at 65 and/or 90 mmHg accordingly.

MORIA strongly recommends having the system thoroughly inspected by MORIA on a routine basis every year.

Please contact MORIA for any other servicing.

B. MAINTENANCE

1. REPLACING FUSE

Steps	What to do	 IMPORTANT WARNINGS	Related pictures
1	Switch OFF ACP system		E1
2	Unscrew the fuse on ACP rear panel		E2
3	Gently pull the fuse on ACP rear panel		E3
4	Remove the cartridge		E4
5	Change the fuse	T1.6AL250V (5x20mm)	E5
6	Gently push the fuse		E6
7	Screw the fuse on ACP rear panel		E7

XIII. WARRANTY

Warning: You may activate your warranty via MORIA website (<http://www.moria-surgical.com>) after being registered, or via the free Moria App app (to download from PlayStore or App Store).

A. SCOPE OF WARRANTY

DESCRIPTION	MORIA reference
ACP system	19188
ACP footswitch	19189
ACP power supply	19190

- The above items as well as labor necessary for their repair are covered by warranty. Any items returned must be sent in their original packaging, after having previously been decontaminated.
- The maintenance operations and the replacement of spare parts will be exclusively carried out by technicians authorized by MORIA.

B. NON-APPLICATION OF WARRANTY

The warranty will not be applicable under any of the following conditions:

- Defects or malfunction that occurs out of the warranty period (paragraph XIII.C).
- Normal wear and tear.
- Negligence or usage that does not comply with the specifications in the user manual.
- ACP system connected to an inappropriate power supply.
- The use of supply, tubing or accessories other than those supplied by MORIA.
- Disassembly, modification, or intervention carried out on the devices by a person not authorized by MORIA.

C. WARRANTY PERIOD

- The warranty takes effect on the date the material is dispatched.
- The duration of the warranty is 12 months from the date of effect.

D. LIABILITY

- The liability of MORIA is limited to the supply of the services mentioned in paragraph XIII.A. MORIA will not be held responsible for any direct or indirect damage notably financial, suffered by the client owing to the interventions within the scope of this warranty.
- For any dispute concerning the interpretation or the execution of the present contract or the present general terms and conditions, the Commercial Court of Nanterre, France will have sole jurisdiction.

XIV. DRAWINGS

A. ACP SYSTEM (FRONT)

Figure	Description
1	Cover
2	Pressure output
3	Digital display
4	Turning knob
5	Yellow LED – No pressure or overpressure = 250 mmHg
6	Green LED – Pressure ON

B. ACP SYSTEM (BACK)

Figure	Description
1	Power supply information
2	Switch ON/OFF
3	Alimentation connector
4	Fuse T1.6AL250V (5x20mm)
5	ACP footswitch connector
6	External air inlet.
7	MORIA reference and serial number

C. ACP FOOTSWITCH

Figure	Description
1	ACP Footswitch

D. POWER SUPPLY

Figure	Description
1	ACP power supply

E. CHANGING FUSE

Figure	Description
E1	Switch OFF ACP system
E2	Unscrew the fuse with a screwdriver
E3	Gently remove the fuse
E4	Change the fuse
E5	Replace the fuse
E6	Gently push the fuse
E7	Screw the fuse with a screwdriver

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MORIA S.A. – 27, rue du Pied de Fourche - 03160 Bourbon l'Archambault – France