EN COMPOSITES SQUEEZE FORCEPS Instructions for Cleaning, Sterilization and Maintenance of Moria Composites Squeeze Forceps 23, 25 G

GENERAL COMMENTS

The Composites instrument line has been designed for optimal operation, provided that the recommendations listed in these instructions are followed carefully. The most recent version of this document is available on the Moria website: <http://www.moria-surgical.com>.

It is the responsibility of the processor to ensure that the processing is performed using proper equipment, materials, and personnel to achieve the desired result. Validation and routine monitoring of the process are required. Any deviations from these instructions by the processor should be validated prior to use for effectiveness and potential adverse consequences.

DESCRIPTION AND INTENDED USE

Moria Composites instruments are reusable handheld devices for performing ophthalmic surgical procedures. The active elements of the forceps are surgical stainless steel and are attached to handles made from a composite polymer. The use of these instruments is apparent to qualified personnel.

PRECAUTIONS AND WARNINGS

⚠ Moria shall not be liable for any direct or indirect injury experienced by the user or his patients in the event of inappropriate use, care, or maintenance of the instruments.

⚠ Moria disclaims all liability in the event of an accident resulting from repair or modification of the instruments by persons not authorized by Moria.

⚠ The instruments should be used only by physicians and well-trained personnel.

⚠ **Caution USA only: Federal law restricts the use of these devices to physicians or licensed practitioners.**

⚠ Ensure that the selected instrument is appropriate for the intended medical procedure, is in working order, and that it will not interfere with other medical devices used during the procedure.

⚠ If any instrument is received in defective condition, notify Moria. Never use damaged or defective instruments.

⚠ Instruments are delivered non-sterile and must be cleaned and sterilized before use.

⚠ The instruments are supplied with rigid protectors. ⚠ Instruments should be immediately cleaned and sterilized after each use.

⚠ Incompatibilities: Do not use bleach.

⚠ For effective sterilization, check that the instruments are entirely clean and free of any residue.

⚠ After sterilization, and prior to each use, inspect instruments under a microscope (minimum 20 times magnification) for changes in physical characteristics that may affect their safe use, effectiveness, or performance. Misaligned forceps tips, bent instruments, dull instruments, and degradation, damage, or loosening of composite material handles indicate that an instrument is at the end of its product life.

In particular:

- Discard forceps with misaligned tips.
- Discard instruments that no longer have the requisite sharpness for satisfactory performance.
- Discard instruments in which the metal parts are separated or loose from the composite material handles.
- Discard instruments with degraded or damaged handles.
- Discard instruments with evidence of oxidation. The above list is not all-inclusive. Sound professional judgment must always be used in determining whether an instrument is in appropriate condition for use.
- ⚠ In case of doubt, consult the Moria website, the Moria catalogue, or your Moria representative to obtain further information.

Cleaning (see addendum 65100 on Moria website for pictures and material used)

A. LA TO DO IN POINT OF USE

- a. Right after use, place the protector on the instrument to protect the tip.
- b. Post procedure without delay (with the tip protector on the instrument); immerse the instrument in water*.

B. TO DO IN CLEANING SPACE

1. TAP WATER RINSING

- a. Rinse instrument (with its tip protector) with warm tap (or demineralized) water (30°C – 40°C / 85°F – 105°F) for 30 seconds minimum, rotating it to expose all surfaces and cavities to flowing water.

2. MANUAL CLEANING

- a. Immerse the instrument (with its tip protector) in warm neutral pH cleaning solution* (30°C –

- a. 40°C / 85°F – 105°F) to moisten and loosen soil.
- b. Remove the protector and scrub all surfaces of instrument (without tip protector) with the cleaning solution using a soft cleaning brush, without damaging the tip, until all visible soil has been removed.
- c. Visually inspect, and if necessary repeat step "b", until all visible soil has been removed.
- d. Activate mechanism of distally controlled instrument (with tip protector) 5 times.
- e. Clean the distal end of instrument (with tip protector) by flushing the cleaning solution 5 times per lumen, using a 5ml syringe and the special cannula provided.

3. TAP WATER RINSING

- a. Rinse the instrument (with tip protector) with warm tap water (30°C – 40°C / 85°F – 105°F) for 30 seconds minimum, rotating it to expose all surfaces and cavities to flowing water.

4. US BATH CLEANING

- a. Immerse the instrument (with tip protector) in an ultrasonic bath filled with warm fresh neutral pH cleaning solution (30°C – 40°C / 85°F – 105°F) and sonicate for 5 minutes*.
- b. Check that the instrument is cleaned and free from any visible soil (ideally under microscope). Repeat process if necessary.
- c. Recommended to repeat process as necessary.

5. TAP WATER RINSING

- a. Rinse the instrument (with tip protector) with warm tap water (30°C – 40°C / 85°F – 105°F) 30 seconds minimum, rotating it to expose all surfaces and cavities to flowing water.

6. RINSING

- a. Rinse the distal end of instrument by thoroughly flushing the lumens with warm demineralized water (30°C – 40°C / 85°F – 105°F) using the special cannula provided.
- b. Immerse instrument (with tip protector) in a bath of warm fresh demineralized water (30°C – 40°C / 85°F – 105°F).

7. DRYING

- a. Rinse the distal end of the instrument tip (with tip protector), each lumen, with 5 cc of 70° alcohol using the special cannula provided.
- b. Remove alcohol by flushing the distal end of the instrument thoroughly, each lumen, with filtered dry air using the special cannula provided.

- a. Wipe the instrument with a non-woven lint free disposable cloth.
- d. Air dry.

(*) in case of use of neutral pH cleaning solution , refer to manufacturer recommendation of use regarding concentration, duration, T°.

STERILIZATION

For effective sterilization, check that the instruments are entirely clean and free of any residue.

Sterilization with autoclave (steam sterilizer).

Steam sterilization cycle	Preparation	Exposure time ¹	Temperature	Minimum drying time, in room chamber ²
Pre-vacuum	wrapped	18 minutes	134°C / 273°F	20 minutes
Pre-vacuum (UK)	wrapped	3 minutes	134°C / 273°F	20 minutes
Pre-vacuum (USA) ³	wrapped	4 minutes	132°C / 270°F	20 minutes

- ¹ US only, Note: As per ANSI/AAMI ST79.
- ² Sterilization temperature must not exceed 140°C/284°F.
- ³ Use the autoclave drying cycle to avoid oxidation.

Moria recommends the above sterilization processes. The use of a cleaning and/or sterilization protocol, other than the one recommended by Moria, will be the responsibility of the end user. Users must validate their sterilization process.

The handle made of a composite material may change color after several cleaning and sterilization cycles but this change has no effect on its integrity.

Inactivation of UTA (unconventional transmissible agents) (only for patients at risk of Creutzfeldt-Jakob disease – CJD questionnaire / please refer to French Ministry instruction. DGS/R13/2011/449-December 2011).

STORAGE

- Store instruments in clean dry conditions at a temperature between 15°C and 35°C (59°F and 95°F) and relative humidity between 30% and 75%, in order to prevent risk of condensation inside the packaging.
- Do not store instruments near products or in an environment that may have a corrosive or magnetic effect.
- Avoid contact between instruments.

REPAIRS

These instruments are not intended for repair.

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