FIBEROPTIC PROBES
20, 23, 25 OR 27GAUGE, STERILE, SINGLE USE
(0.81mm, 0.64mm, 0.50mm, 0.41mm)

DESCRIPTION/INDICATION FOR USE (by physicians in the operating room)

Laser Probes are intended for use in vitreoretinal surgery to perform endo-ocular laser photocoagulation treatments at operating wavelengths of 500nm to 900nm.

Illuminated Laser Probes are intended for use in vitreoretinal surgery to perform endo-ocular laser photocoagulation treatments at operating wavelengths of 500nm to 900nm and to illuminate a surgical site.

Light Pipes and Chandeliers are intended for use in vitreoretinal surgery to illuminate a surgical site.

Replaceable Fibers are intended for use in vitreoretinal surgery to perform endo-ocular laser photocoagulation treatments at operating wavelengths of 500nm to 900nm.

CAUTION: FEDERAL (U.S.A) LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN

- INSTRUCTIONS FOR USE:
- Laser and Illuminated Laser Probes:
- Release the probe from pouches and remove protective cap.
 - Insert the laser coupler into the laser machine and – if applicable – the 19ga Needle into the illumination adapter that is plugged into the illumination machine.
 - Follow the recommendations of the laser and – if applicable – illumination source manufacturers (e.g.: calibration of the laser source as well as the need for protective eye wear for the corresponding wavelength of the laser model used).
- The Probe is now ready for use.
- Light Pipes and Chandeliers:
- Release the probe from pouches and remove protective cap.
 - Insert the 19ga Needle into the illumination adapter that is plugged into the illumination machine.
 - Follow the recommendation of the illumination source manufacturers.
 - Insert the supplied cannula into the eye.
 - Install the chandelier into the cannula.
 - The Light Pipe is now ready for use.
- Replaceable Fibers:
- The coupler, handpiece and fiber assembly are designed to be used together and all must be supplied by Katalyst Surgical, LLC to ensure proper function.
 - Autoclave sterilize the handpiece according to the IFU instructions. In the sterile field, identify the proximal side of the hand piece. Insert the fiber carefully from this side and thread the fiber until it is fully threaded, approximately 3/4 of a turn.
 - The coupler does not need to be sterilized. Insert the fiber carefully from this side and thread the fiber until it is fully threaded, approximately 3/4th of a turn. Connect the connector to the laser machine in accordance with the laser manufacturer's IFU..
 - The Probe is now ready for use.

- Additional instructions: Chandelier Infusion Cannula
- Catalog Numbers: V220X-YYZ
- Follow above bullets under "Light Pipes and Chandeliers"
 - Prime the cannula, removing all air bubbles.
 - Insert the supplied cannula into the eye.
 - Install the chandelier into the cannula. Push chandelier tubing over cannula body.
 - Fixate the infusion tubing near chin on side opposite the affected eye. Keep tubing as vertical as possible to direct illumination towards the posterior pole.
 - Adjust malleable portion to conform to the patient's facial contours.
 - The system is now prepared for specified use.

- FEATURES:
- Different couplers that suit most laser/illumination source models available on the market.
 - Sterilized by Ethylene Oxide, single-use, latex-free.
 - Dia 20ga, 23g, 25g or 27ga. Length 2.4m (8ft)
 - Illuminated Laser Probes, Chandeliers and Light Pipes come with a 19ga Needle connector. Adapters sold separately

RECOMMENDATION:

As the aiming beam passes down the same delivery system as the working beam, it provides a good means of checking the integrity of the delivery system. If the aiming beam is not present at the distal end of the delivery system, its intensity is reduced or it looks diffused, this is a possible indication of a damaged and malfunctioning delivery system.

It is recommended that the operator wear eye protection while using the Fiberoptic Probes during surgery.

CAUTION:

This device is to be used only by qualified physicians trained on the procedure and – if applicable – illumination machines.

Follow medical operating procedures during surgery.

This product is intended for single use, one patient only.

Reuse and re-sterilization may cause diminished product performance, contamination or infection.

Inserting, sharply bending, or improperly securing the fiberoptics may lead to damage of the fiber or delivery system and/or harm to the patient or machine operator. Please follow the recommendations from both Katalyst and the laser and/or illumination machine manufacturer to ensure proper fiberoptic handling.

PRECAUTION:

Do not use if package is open, damaged or wet.

Do not use if expiry date shown on the label is exceeded.

STORAGE:

Store in a dry area.

ADVERSE REACTIONS:

Phototoxicity. Complications are statistically. Assumed the user is adequately trained in the treatment and methods of avoidance to these complications..

WARNINGS:

Due to the potentially hazardous effects of phototoxicity to the retina from prolonged Endoillumination, this device should be used with an illumination source with filters that eliminate UV radiation (<400nm) and whenever possible, with filters that eliminate short wavelength blue light (<420nm). Care should be taken not to focus this device on a single point of the retina for prolonged, unnecessary periods of time.

ENGLISH

COMPATIBILITY WITH OTHER DEVICES:

The available laser machine connectors (A-J) defined by name of the compatible machine manufacturer and are able to connect to any machine made by the listed company.

Illumination machine connectors (R or V) define the specific machine compatible for connectivity.

MINIMUM EXPOSURE GUIDELINES:

Refer to the table below ISO 15752:2012 endoilluminator light guide maximum exposure guidelines.

Refer to Katalyst IFU 7507-0012-00 for other endoilluminator light sources.

Endoilluminator Light Source	Light Guide Optical Fiber Diameter (µm)	Field Angle ¹ (degrees)	Recommended Working Distance (mm)	Time to Reach Maximum Exposure Guideline			
				Without Retinal Protection Means at Maximum Light Output (minutes)	Without Retinal Protection Means at 50% of Maximum Light Output (minutes)	With Retinal Protection Means at Maximum Light Output (minutes)	With Retinal Protection Means at 50% of Maximum Light Output (minutes)
Synergetics Photon	250	60	15	27	28	27	28
		1. Maximum exposure times may be significantly reduced for field angles smaller than those shown. The light emitted from this instrument is potentially hazardous. See Notes 1 to 4. 2. Follow all exposure warnings.					
		Note 1: Maximum exposure times are for cumulative retinal exposure with a stationary distal tip of the light guide positioned at the specified distances from the retina. Changing the distance of the endoilluminator light guide from the retina will also significantly affect the risk factor. Note 2: Lower intensities increase the maximum exposure times in direct proportion to the decrease in intensities. Note 3: Movement of the light guide increases safe exposure time. Note 4: Maximum exposure times are given for clear media. Vitreous hemorrhage will increase these times.					

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Katalyst Surgical, LLC 754 Goddard Ave. Chesterfield, MO USA 63005