REVOLVER FIBEROPTIC PROBE REUSABLE PARTS
20, 23, 25 or 27 GAUGE, NON-STERILE

INDICATION FOR USE (by physicians in the operating room):

Revolver Laser Probes, Revolver Illuminated Laser Probes and Revolver Illuminated Probes are intended for use in vitreoretinal surgery to perform endo-ocular laser photocoagulation treatments. Revolver Laser Probes can only be used with a medical laser at operating wavelengths of 500nm to 900nm. The Revolver assemblies containing illumination should only be used with the light wavelength range of 425nm to 700nm.

DESCRIPTION:

The Revolver Fiberoptic Probe Reusable Parts are intended to be used in conjunction with Revolver Fiberoptic Replaceable Fibers. The Revolver Fiberoptic Replaceable Fiber will thread into the hand piece. The coupler is designed to connect the fiber assembly to the illumination and/or laser machine.

CAUTION: FEDERAL (U.S.A) LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN

The coupler, hand piece and fiber assembly are designed to be used together and all must be supplied by Katalyst Surgical, LLC to ensure proper function.

INSTRUCTIONS FOR USE:

- In the sterile field, identify the proximal side of the hand piece. Insert the fiber carefully from this side and thread the fiber until it is fully threaded, approximately 3/4 of a turn.
- Identify the proximal side of the coupler. Insert the fiber carefully from this side and thread the fiber until it is fully threaded, approximately 3/4th of a turn. Connect the connector to the laser machine in accordance with the laser manufacturer's IFU.
- The Probe is now ready for use.

CLEANING AND STERILIZATION OF THE HAND PIECE AND COUPLER:

- Initial Use of New Instruments** - Every instrument must be cleaned and sterilized before it is used for the first time
- Inspection and Functional Check** - It is very important to examine each surgical instrument for breaks, cracks or malfunctions before each use. Check insulation, cables and couplers for cuts, voids, cracks, tears, abrasions, etc. Do not use damaged instruments. Never attempt to make repairs yourself. Service and repairs should be referred to trained qualified persons only. Refer questions about repair to Katalyst Surgical, LLC customer service.
- Sterilization**
 - Sterilize instruments according to hospital guidelines.
 - Steam autoclaving has been validated in the wrapped configuration and can be sterilized wrapped or unwrapped.
- Standard autoclave cycle: Steam sterilize at 270°F (132°C) for 4 minutes or more. 20 minutes of dry time has been validated and was determined to be sufficient to completely dry the instruments. Other time and steam temperature cycles may also be used. However, user must validate any deviation from the recommended time and temperature.
- Contact the manufacturer of your steam autoclave to confirm appropriate temperatures and sterilization times.
- Autoclave temperatures should not exceed 280°F (137°C); handles, insulation or other non-metallic parts may be damaged.
 - Use of a plasma sterilization will compromise instrument life, and will void the warranty.
 - ETO sterilization is recommended due to low temperature and moisture levels experienced during processing.
- Cleaning**
 - Clean and rinse instruments immediately after each use. Failure to clean promptly may result in adherent particles that may resist cleaning and complicate or resist future sterilization. Instruments must be completely cleaned and rinsed of all foreign matter. Proper cleaning and maintenance procedures will ensure correct function and extend the useful life of the instrument
 - Remove blood, saline, tissue and other contaminants immediately by flushing with distilled/deionized water.
To flush: Attach a syringe to the back end of the handle and push fluid forward through the center of the instrument.
 - Clean instruments with a mild soap.
Do not use corrosive cleaning agents (i.e. bleach). Cleaning solutions and rinses at or near a neutral pH (7.0) are best.
Do not use abrasive cleaners.
Enzymatic cleaners should be used to remove protein deposits. If used, follow the enzymatic cleaner's instructions and rinse thoroughly.
 - Flush thoroughly with distilled/deionized water.
 - Flush with alcohol.
 - Dry instruments completely and carefully with compressed air (including inside channels and highly inaccessible areas).
 - Inspect and test each instrument carefully.
- Ultrasonic Cleaning and Drying**
 - Follow the manufacturer's instruction for proper use of your ultrasonic cleaner.
 - Use only detergents with distilled/deionized water (pH 7.0).
 - Place tip protectors on instruments before placing in ultrasonic cleaner. (DO NOT ALLOW the instruments to bump each other while in the tray.)
 - Immerse instruments together in a mesh type tray for the specified time period as described by the ultrasonic cleaner's instructions. This time period usually runs less than 5 minutes.
 - Dry instruments completely and carefully with compressed air (including inside channels and highly inaccessible areas).
- Storage and Sterilization**
 - Hand pieces should be stored in a clean, dry, moisture free area. Microsurgical instruments should be stored individually in their shipping packaging or in a protective tray with partitions. Always keep a protective cover on delicate tips when the instrument is not in use. Allow instrument to air cool to room temperature before use.

REORDE CODES:

Handle Configurations	Product Code	Coupler	Product Code
Straight Handle	RLP000-2X	Universal Laser Coupler	RLP-A
Curved Handle	RLP100-2X	Quantel Laser Coupler	RLP-C
Flexible Curved Handle	RLP111-2X	DORC Laser Coupler	RLP-J
Steerable Handle	RLP222-2X	Nidek Laser Coupler	RLP-E
Extendable Steerable Handle	RLP222E-2X	Constellation Illumination Adapter	RLP-Q
		Stellaris Illumination Adapter	RLP-S
		DORC Illumination Adapter	RLP-U
		Short 19ga Needle Illumination Adapter	RLP-V
		Zeiss Illumination Adapter	RLP-Z

FEATURES:

- Different couplers that suit most laser/illumination source models available on the market.
- Dia 20ga, 23g, 25g or 27ga.

CAUTION:

This device is to be used only by qualified physicians trained on the procedure. Follow medical operating procedures during surgery.

PRECAUTION:

Do not use if instrument is damaged.

ADVERSE REACTIONS:

There are no adverse reactions associated with the use of this product other than those related to the specific procedure being performed (to be assessed by the physician).

WARNINGS:

Due to the potentially hazardous effects of phototoxicity to the retina from prolonged Endoillumination, this device should be used with an illumination source with filters that eliminate UV radiation (<400nm) and whenever possible, with filters that eliminate short wavelength blue light (<420nm). Care should be taken not to focus this device on a single point of the retina for prolonged, unnecessary periods of time.

COMPATIBILITY WITH OTHER DEVICES:

Refer to the **product label to determine connector type and the machine(s) that this device may be connected to.**

MAXIMUM EXPOSURE GUIDELINES:

Refer to the **table below ISO 15752:2010 endoilluminator light guide maximum exposure guidelines.**

Endoilluminator Light Source	Light Guide Optical Fiber Diameter (µm)	Field Angle (degrees)	Recommended Working Distance (mm)	Time to Reach Maximum Exposure Guideline			
				Without Retinal Protection Means at Maximum Light Output (minutes)	Without Retinal Protection Means at 50% of Maximum Light Output (minutes)	With Retinal Protection Means at Maximum Light Output (minutes)	With Retinal Protection Means at 50% of Maximum Light Output (minutes)
Accurus High Brightness	250	60	15	N/A	N/A	62	192
Alcon Accurus	250	60	15	N/A	N/A	468	689
Alcon	250	60	15	N/A	N/A	26	140
B&L Millennium	250	60	15	N/A	N/A	379	758
B&L Stellaris PC	250	60	15	N/A	N/A	65	137
Synergetics Photon	250	60	15	N/A	N/A	27	28
DORC Eva	250	60	15	22	10	22	10
1. Maximum exposure times may be significantly reduced for field angles smaller than those shown. The light emitted from this instrument is potentially hazardous. See Notes 1 to 4. 2. Follow all exposure warnings.							
Note 1: Maximum exposure times are for cumulative retinal exposure with a stationary distal tip of the light guide positioned at the specified distances from the retina. Changing the distance of the endoilluminator light guide from the retina will also significantly affect the risk factor. Note 2: Lower intensities increase the maximum exposure times in direct proportion to the decrease in intensities. Note 3: Movement of the light guide increases safe exposure time. Note 4: Maximum exposure times are given for clear media. Vitreous hemorrhage will increase these times. Note 5: Since the focal light pipes produce the smallest beam angle of any Katalyst endoillumination product, they produce the most energy dense light beam, meaning the light provided by the illuminator source will be absorbed by a smaller area of the retina. Thus, the focal light pipe is considered to represent the worst case scenario product when compared to other endoilluminators.							

■ Katalyst Surgical, LLC 754 Goddard Ave. Chesterfield, MO USA 63005
Tel: 636-536-5950 or 888-452-8259 Fax: 636-787-0603 www.katalystsurgical.com

Rx only

Manufacturer

Lot Number

Catalog Number

Attention: See Instructions for Use

Quantity

Rx only

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

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🏠Katalyst Surgical, LLC 754 Goddard Ave. Chesterfield, MO 63005, USA