



Moria



SBK without compromise Clinical Data Sheet



One Use-Plus **SBK**






Predictible thin sub-Bowman flaps

| # | Author | Year | Nb of eyes | Clinical topic | Clinical proof |
|-------|-------------------|------|------------|----------------|---|
| 2 | Chen et al. | 2010 | 82 | | <p><i>"The central flap thickness was dramatically thinner in the OUP SBK group":</i></p> <ul style="list-style-type: none"> intended: 110 µm (Speed 1) accuracy: OD: 114.7 µm / OS: 109.4 µm predictability: OD: ± 10.1 µm / OS: ± 11.0 µm reproducibility: 5.3-µm difference between OD and OS |
| 3 | Lian et al. | 2010 | 20 | | <p><i>"The least variation in flap thickness is SBK":</i></p> <ul style="list-style-type: none"> intended: 100 µm (Speed 2) accuracy: OD: 97.50 µm / OS: 96.73 µm predictability: OD: ± 11.39 µm / OS: ± 10.45 µm reproducibility: <1.0-µm difference between OD and OS |
| 4 | Du et al. | 2011 | 60 | | <p><i>"The SBK head demonstrated the most accurate flap thickness":</i></p> <ul style="list-style-type: none"> intended: 100 µm (Speed 2) accuracy: OD: 97.50 µm / OS: 96.73 µm predictability: OD: ± 11.39 µm / OS: ± 10.45 µm reproducibility: <1.0-µm difference between OD and OS |
| 5 | Sun et al. | 2012 | 57 | | <p><i>"SBK is better than LASIK in creating much uniform corneal flap":</i></p> <ul style="list-style-type: none"> ±2 mm from apex: 91.09 ± 7.85 µm ±3 mm from apex: 89.58 ± 7.88 µm temporal: 89.30 ± 7.64 µm nasal: 89.20 ± 7.96 µm superior: 89.91 ± 8.38 µm inferior: 90.44 ± 7.69 µm |
| 6 | Zhang et al. | 2012 | 32 | | <p>Mean central flap thickness at Speed 1:</p> <ul style="list-style-type: none"> intended: 110 µm obtained: 108.15 µm → accuracy: <2.0 µm from intended |
| 7 | Zhai et al. | 2013 | 44 | | <p>Mean central flap thickness at Speed 1:</p> <ul style="list-style-type: none"> intended: 110 µm obtained: 113.85 ± 8.07 [97.50-130.00] µm → accuracy: <2. µm from intended |
| 8 | Zhang et al. | 2014 | 60 | | <p><i>"The flaps in the [...] SBK group were more regular, showing an almost planar configuration".</i></p> <p>Mean central flap thickness at Speed 1:</p> <ul style="list-style-type: none"> intended: 110 µm accuracy: OD: 110.6 µm / OS: 108.2 µm predictability: OD: ± 7.4 µm / OS: ± 6.1 µm reproducibility: 2.4-µm difference between OD and OS |
| 9 | Al-Thomali et al. | 2014 | 70 | | <p><i>"The One Use-Plus SBK is a reliable microkeratome with reasonable predictability for the creation of SBK flaps":</i></p> <p>→ mean central flap thickness: 88.74 µm</p> |
| 10 | Mimouni et al. | 2015 | 2560 | | <p><i>"The eyes in [SBK] group had statistically significantly thinner flaps" in this large-scale study:</i></p> <ul style="list-style-type: none"> intended: 110 µm (Speed 1) accuracy: OD: 110.1 µm / OS: 107.6 µm → 108.9 µm predictability: OD: ± 15.2 µm / OS: ± 15.4 µm → ± 15.3 µm reproducibility: 2.5-µm difference between OD and OS |
| 11 | Xu et al. | 2015 | 40 | | <p><i>"At each time period after SBK, flap thickness in the Moria group was significantly thinner:</i></p> <ul style="list-style-type: none"> at 1 day post-op: 95.8 ± 7.6 µm at 1 week post-op: 95.5 ± 7.8 µm at 2 weeks post-op: 96.5 ± 7.9 µm at 1 month post-op: 98.2 ± 8.2 µm |
| 12 | Katz et al. | 2015 | 344 | | <p>Mean central flap thickness at Speed 2:</p> <ul style="list-style-type: none"> intended: 100 µm accuracy: 96 µm <p>In a subgroup of 132 eyes of 66 patients:</p> <ul style="list-style-type: none"> accuracy: OD: 97.29 µm / OS: 93.77 µm reproducibility: <4.0-µm difference between OD and OS |
| TOTAL | | | 3369 | | |

CLINICAL DATA SHEET




Flap creation in less than 4 seconds

| # | Author | Year | Nb of eyes | Clinical topic | Clinical proof |
|-------|-----------|------|------------|---|---|
| 1 | Xu et al. | 2014 | 82 |  | "The application of the microkeratome suction ring has been shown in other studies to induce changes in the perilimbal conjunctiva, especially on goblet cell density, which contributes to the pathology of dry eye. <i>In Femto-LASIK and SMILE procedures, the control suction is longer and lighter. The effect on goblet cells needs to be investigated in future studies.</i> " |
| TOTAL | | | 82 | | |




Very fast visual recovery

| # | Author | Year | Nb of eyes | Clinical topic | Clinical proof |
|-------|-----------|------|------------|---|--|
| 1 | Xu et al. | 2014 | 82 |  | "No significant difference was found between groups [SBK vs SMILE vs LASIK with either a femtosecond laser or mechanical microkeratome]" |
| 11 | Xu et al. | 2015 | 40 | | "At 1 month postop, all the eyes had UCVA of 20/20 or better [SBK vs Femto-LASIK]." |
| TOTAL | | | 122 | | |



Excellent safety profile

| # | Author | Year | Nb of eyes | Clinical topic | Clinical proof |
|-------|-------------------|------|------------|---|--|
| 14 | Falcon et al. | 2016 | 2883 |  | "The Moria One Use-Plus SBK microkeratome is an excellent device that allows the easy creation of thin corneal flaps and regular and <i>smooth corneal beds</i> , safely and predictably even in extremely flat corneas without complications." "There were <i>no intraoperative nor immediate postoperative complications</i> : - intraoperatively: free cap 0%; incomplete flap 0%; button hole 0%; epithelial erosion 0%; bleeding 0%; irregular stromal bed 0% - postoperatively: flap displacement 0%; punctate keratitis: 8%; LASIK retreatments: 12%; safety: 100% (no loss of lines of visual acuity)." |
| 15 | Gauthier et al. | 2019 | 146 | | "The cornea was exposed and a flap was cut at 100 µm with a Moria microkeratome, <i>selecting the suction ring and the settings to achieve maximum flap diameters (> 9.5 mm).</i> " "High hyperopia can be treated with LASIK, if wide OZ, TZ, and TAZ are used, and <i>the flap is large enough to perform this extreme peripheral photoablation.</i> " |
| 16 | Kasetsuwan et al. | 2016 | 157 | | "In the femtosecond laser group, 2.0% of eyes had early postoperative diffuse lamellar keratitis (DLK) <i>compared with 0% in the [Moria SBK] microkeratome group.</i> " |
| 17 | Friehmann et al. | 2018 | 30574* | | "In this [large-scale retrospective cohort] study, the incidence rate (0.49%) of epithelial ingrowth was <i>somewhat lower than previously reported studies.</i> " |
| 18 | Pokroy et al. | 2016 | 9177 | | "Retreatment rates decrease as surgeon's experience is gained and technology improved. <i>The myopic LASIK retreatment rate in recent years improved to below 0.5%.</i> " |
| 19 | Mimouni et al. | 2018 | 1104 | | "In this study, the use of a Moria M2-90 microkeratome (as opposed to Moria SBK-90) was associated with a higher retreatment rate." |
| TOTAL | | | > 13467 | | |

* 30574 eyes of 15287 patients treated with two different automated microkeratomes (M2SU 90 then One Use-Plus SBK by Moria) between January 2006 and December 2014; split not provided by the authors

CLINICAL DATA SHEET



Excellent stromal surface smoothness

| # | Author | Year | Nb of eyes | Clinical topic | Clinical proof |
|-------|---------------|------|------------|----------------|--|
| 6 | Zhang et al. | 2012 | 32 | | "The flap margin after the [...] OUP SBK technique[s] appeared microscopically as a <i>clear-cut edge</i> ." |
| 13 | Duffey et al. | 2008 | - | | Laboratory study with human corneas not suitable for transplantation. |
| TOTAL | | | 32 | | |



Supported by in-vivo confocal microscopy

| # | Author | Year | Nb of eyes | Clinical topic | Clinical proof |
|-------|--------------|------|------------|----------------|---|
| 6 | Zhang et al. | 2012 | 32 | | "[...] in the OUP SBK group, it [the density of stromal keratocytes] was higher than that of the FS-LASIK [...] at 3 months postsurgery ($P<0.05$)." "The repairing velocity of subbasal nerve fibers in the OUP SBK group was a little faster than that of the FS and conventional LASIK groups." |
| TOTAL | | | 32 | | |



Excellent quality of vision

| # | Author | Year | Nb of eyes | Clinical topic | Clinical proof |
|-------|-------------------|------|------------|----------------|---|
| 9 | Al-Thomali et al. | 2014 | 70 | | "[...] the absence of correlation of flap thickness with postoperative CDVA and induced HOAs [Higher Order Aberrations] indicates that flap thickness does not influence visual outcomes. |
| 20 | Hassanin et al. | 2013 | 114 | | To evaluate changes in corneal HOAs (spherical, coma, trefoil, aberration coefficients, total aberration coefficient) as well as corneal asphericity (Q-value) following optimized LASIK ablation for moderate to highly myopic eyes. |
| 21 | McAlinden et al. | 2010 | 65 | | To compare the change in HOAs after LASIK and LASEK and to determine which method of flap creation induced a greater increase in HOAs. |
| 22 | McAlinden et al. | 2011 | 100 | | To investigate the internal HOAs following LASIK. |
| 23 | Malhotra et al. | 2015 | 50 | | "Flaps created using the Moria Evo 3 One Use-Plus SBK mechanical microkeratome were associated with significantly lower induction of total HOAs and spherical aberrations as compared with 150 KHz iFS flaps at the end of 3 months follow-up." |
| 24 | Wang et al. | 2013 | 67 | | "When comparing straylight values in the femtosecond laser group with those in the [OUP SBK] mechanical microkeratome group, differences were not statistically significant at any follow-up time point." |
| TOTAL | | | 466 | | |

REFERENCES

Refer to brochure #66065 available on Moria website:
<https://www.moria-surgical.com/resources-center>

Moria

MORIA SA
 27, rue du Pied de Fourche
 03160 Bourbon L'Archambault
 FRANCE
 Phone: +33 (0)1 46 74 46 74
 Fax: +33 (0)1 46 74 46 00
moria@moria-int.com
www.moria-surgical.com

MORIA Inc
 1050 Cross Keys Drive
 Doylestown, PA 18902
 USA
 Phone: (800) 441 1314
 Fax: +1 (215) 230 7670
orders@moriausa.com
www.moria-surgical.com

Moria Surgical Spain
 C. de Pablo Iglesias, 58,
 08908 L'Hospitalet de Llobregat,
 Barcelona - Spain
 Phone: +34932896270
oftaltech@oftaltech.com
www.oftaltech.com

Moria Japan K.K.
 Arcadia Building 6F
 1-12-3 Kanda
 SudachoChiyoda-Ku
 Tokyo 101-0041
 JAPAN
 Phone: 81-3-6260-8309
 Fax: 81-3-6260-8310
moria@moriajapan.com
www.moriajapan.com

AL.CHI.M.I.A. S.R.L.
 (a Moria company)
 Viale Austria, 14
 35020 - Ponte San Nicolò (PD)
 ITALY
 Phone: +39 049 89.62.074
 Fax: +39 049 89.62.071
info@alchimiasrl.com
www.alchimiasrl.com

Moria COMMERCIAL (CHINA) CO., LTD.
 上海目利亚贸易有限公司
 Room H 6F Kaili Building NO.432
 West Huai Hai Road, Changning
 district, Shanghai, 200052, P.R.C
 中国上海市长宁区淮海西路432号
 凯利大厦6楼08室
 Phone/Fax: +86 021 52586095
moria@moria-int.com

Moria Limited
 Suite 1904, 19th Floor,
 ST. George's Building,
 No. 2 Ice House Street,
 Central - HONG KONG
moria@moria-int.com
www.moria-surgical.com

