



# Moria



## Hyperopia & Flat Corneas Single-Use Large-Cut Kit

For LARGE to EXTRA-LARGE stromal beds




**Expert-Focus on large LASIK flaps  
for hyperopia & wide treatments**




With latest generation of excimer lasers offering **larger optical zones** and **smoother transition zones**<sup>1</sup>, hyperopic LASIK is now gaining more and more popularity even if **challenges still remain**:

- ▶ hyperopic eyes often present with an unusually large angle Kappa ( $\kappa$ ), resulting in a **nasally displaced light reflex**<sup>2</sup>
- ▶ **size of the transition zone** of latest excimer lasers varies according to the magnitude of hyperopia to be treated<sup>3</sup>



**Specifications for an ideal flap sizing<sup>4</sup>:**  
to be large enough to treat:


- a** - the optical zone
- b** - & the transition zones
- c** - & some additional peripheral tissue while not compromising flap edges



That's why Moria specifically proposes a **Hyperopic Single-Use Large-Cut Kit**:

- the -1 Large-Cut single-use suction ring for regular cases from 39 to 44  $\delta$
- the -1 Large-Cut reusable suction ring for very flat Ks or extra-large flaps

Suction ring		Single-use -1 Large-Cut ①	Reusable -1 Large-Cut ②
Moria nomogram		#65041	#65067
Calibrated Large-Cut head		110- $\mu$ m & 130- $\mu$ m	
Vertical flap diameter based on K readings (assumptions, mm)	37	-	9.3
	38	-	9.4
	39	9.3	>9.7
	40	9.5	>10.0
	41	9.8	>10.2
	42	10.0	>10.3
	43	10.3	-
	44	10.5	-




### Suction ring

- ① **UU -1LC** → single-use plastic ring, size -1 Large-Cut, coming from single-use pack (head & ring): **#19354/110**
- ② **RU -1LC** → reusable metallic ring, size -1 Large-Cut, with a purple-colored stopper: **#19519/-1**

### Cutting head

- 110 LC** → 110-micron calibrated Large-Cut head, coming from single-use pack (head & ring): **#19354/110**
- 130 LC** → 130-micron calibrated Large-Cut head, coming from single-use pack (head & ring): **#19354/130**



## My Technique

Dr Laurent Gauthier-Fournet  
(Saint-Jean de Luz, France)

## HYPEROPIA & LARGE TREATMENTS

### Rationale for One Use-Plus microkeratome to allow large treatment zones

The treatment of high hyperopia in corneal refractive surgery has always been a challenge<sup>1</sup>

With latest excimer laser platforms, such as the Amaris by Schwind (Germany), the size of the transition zone varies according to the magnitude of the hyperopia to be treated:


- +1.50  $\delta$  hyperopia: OZ = 6.70 mm → total treatment zone = 7.81 mm
- +6.00  $\delta$  hyperopia: OZ = 6.70 mm → total treatment zone = 9.20 mm\*

\* For such large treatment zones, **a LASIK flap size of at least 10 mm is recommended.**

The higher the amount of ablation in the periphery, the more important it is to optimize the optical quality and therefore to use a keratome to safely make so large LASIK flaps.

To achieve this, **we systematically use the One Use-Plus** (Moria SA, France) linear automated microkeratome<sup>2-5</sup> because:


- it comes with a **unique choice of Large-Cut heads & rings**
- a **nasal hinge is beneficial** in the way:
  - it allows to expose a larger corneal surface compared to the superior hinge
  - it is **less disturbed by the rolling movements of the eye**, more pronounced in hyperopic patients during surgery.



# MY NOMOGRAM FOR HYPEROPIA\*

choice of ring type and stop value based on Horizontal K

Horizontal Keratometry	110 LC head		Flap Diameter (mm)	Stroma available (mm)	Hinge size (mm)
	SU ring -ILC	RU ring -ILC			
39		8.5	assumptions	assumptions	assumptions
39.5		8.5	9.8	9.2	2.4
40		9	10.0	9.5	2.6
40.5		9	10.2	9.7	2.4
41		9	10.3	9.7	2.5
41.5		9	10.4	9.8	2.6
42		9	10.4	9.8	2.6
42.5	9		10.4	9.8	2.6
43	9		10.2	9.7	2.3
43.5	9		10.3	9.7	2.5
44	9		10.4	9.8	2.6
44.5	9		10.5	9.9	2.7
45	9		10.7	10.0	2.8
45.5	9		10.8	10.1	3.0
46	9		>10.8	>10.1	>3.0
46.5	9		>10.8	>10.1	>3.0
47	9		>10.8	>10.1	>3.0



\* Extract from marketing brochure #66105INT - My Nomogram - Gauthier - OUP (EN-FR)

## NOTES

- Both -1 Large-Cut single-use suction ring and -1 Large-Cut reusable suction ring :
  - can only be used with the 110- and 130-µm Large-Cut single-use calibrated heads coming from packs #19354/110 and #19354/130 respectively;
  - can not be used with the 90- and 130-µm single-use calibrated heads coming from boxes #19393/XXX (XXX = 90/130) and from packs #19336/XXX (XXX = 90/130) and #19337/XXX (XXX = 90/130).
- Both -1 Large-Cut single-use suction ring and -1 Large-Cut reusable suction ring can only be used by refractive surgeons having extended surgical experience with the use of the One Use-Plus microkeratome.

## BIBLIOGRAPHY

- Varley et al.** LASIK for hyperopia, hyperopic astigmatism, and mixed astigmatism. A report by the American Academy of Ophthalmology. *Ophthalmology* **2004**;111:1604-1617
- Malecaze et al.** A randomized comparison of pupil-centered versus vertex-centered ablation in LASIK correction of hyperopia. *Am J Ophthalmol.* **2011**;152:591-599
- Gauthier-Fournet L.** LASIK for high hyperopia with the Schwind AMARIS. Schwind website: <https://www.eye-tech-solutions.com/en/infocenter/userreports/dr-gauthier-lasik-high-hyperopia-with-schwind-amaris>. Accessed on April 17<sup>th</sup> 2024
- Gauthier-Fournet L.** Correction des très fortes hypermétropies par LASIK avec profils spécifiques et flaps extra-larges. Presented during: SAFIR 2019 annual congress, May **2019**, Paris, France
- Gauthier-Fournet et al.** Six-month outcomes after high hyperopia correction using LASIK with a large ablation zone. *Cornea* **2019**;38:1147-1153

**MORIA SA**  
27, rue du Pied de Fourche  
03160 Bourbon L'Archambault  
FRANCE  
Phone: +33 (0)1 46 74 46 74  
Fax: +33 (0)1 46 74 46 00  
[moria@moria-int.com](mailto:moria@moria-int.com)  
[www.moria-surgical.com](http://www.moria-surgical.com)

**MORIA Inc**  
1050 Cross Keys Drive  
Doylestown, PA 18902  
USA  
Phone: (800) 441 1314  
Fax: +1 (215) 230 7670  
[orders@moriausa.com](mailto:orders@moriausa.com)  
[www.moria-surgical.com](http://www.moria-surgical.com)

**Moria**

**Moria Japan K.K.**  
Arcadia Building 6F  
1-12-3 Kanda  
SudachoChiyoda-Ku  
Tokyo 101-0041  
JAPAN  
Phone: 81-3-6260-8309  
Fax: 81-3-6260-8310  
[moria@moriajapan.com](mailto:moria@moriajapan.com)  
[www.moriajapan.com](http://www.moriajapan.com)

**Moria COMMERCIAL (CHINA) CO., LTD.**  
上海目利亚贸易有限公司  
Room H 6F Kailli Building NO.432  
West Huai Hai Road, Changning district, Shanghai, 200052, P.R.C  
中国上海市长宁区淮海西路432号  
凯利大厦6楼08室  
Phone/Fax : +86 021 52586095  
[moriachina@moria-int.com](mailto:moriachina@moria-int.com)



MORIA #66107EN-B-04.2024