EN

Instructions for use RETINA FORCEPS – SINGLE USE

1- INTENDED USE, PERFORMANCE, AND INDICATION FOR USE

All of these instruments are used in ophthalmic surgeries and more specifically for what it is described in the table.

References
17456/23; 17456/25; 17456/27; 17457/23; 17457/25; 17457/27; 17458/23; 17458/25; 17458/27; 17459/23; 17459/25; 17459/27; 17460/23; 17460/25; 17460/27; 17461/25 ; 17461/27
Intended Use
To grasp and peel the vitreoretinal membranes
Indication for use
Vitreoretinal surgeries

To the best of our knowledge, these retina forceps do not have any restriction to the target population (age, sex and health status) and there are no contraindications.

2- DESCRIPTION/FEATURES

These instruments are for single use and ready for use (Ethylene Oxide sterilized).

The retina forceps have a polymer handle and a metallic active part. The active parts of the retina forceps have different designs to meet surgical needs.

3- PRE-REQUISITES BEFORE USE AND DIRECTIONS FOR USE

These instruments do not require specific direction for use because they are used safely by experienced surgeons who master these instruments and surgeries.

For reminder, these forceps are designed to be inserted in the posterior segment of the eye by the entrance hole. Insert the forceps in the eye and remove it from the eye in a closing position.

At the end of the surgical procedure, discard the single-use forceps into an appropriate container for destruction (according to appropriate regulations)

MORIA SA
27, rue Pied de Fourche
03160 Bourbon l'Archambault
FRANCE
Phone: +33 (0)4 70 67 09 04
Fax: +33 (0)4 70 67 06 61

Moria

MORIA #65166-B-04.2025

4- WARNING / PRECAUTION / ADVERSE EVENTS

This device must be used only by experienced ophthalmic surgeons. These instructions do not describe the steps necessary to perform the surgery. The surgical procedure remains the responsibility of the surgeon experienced in performing the surgery. Each instrument should be used for its purpose/design and no other.

- Store this device in a clean, dry location at room temperature and away from direct sunlight.
- Do not use if expiration date has passed.
- Do not use if the packaging is damaged or if the packaging has been opened unintentionally (check the integrity of the primary packaging and inspect the active part under a microscope before use).
- Do not have any physical contact with the active part to avoid injuries and/or eventual contaminations.
- Stop the use of the device in case of any malfunction during the procedure, including potential complication described in part 5.
- Do not reuse: its design is for 'single use only'; Reuse may affect its mechanical and biological features, and its clinical performances; and may cause device failures or may exposes the patient to adverse events such as contamination, bacterial infections, or inflammation, irritation, eye lesion.

5- POSSIBLE SIDE EFFECTS AND INFORMATION TO RELAY

To date, the identified risks by the manufacturer associated to retina forceps are:

- Functional issues

Any serious incident occurring in conjunction with these instruments must be reported to the manufacturer and the competent authority of the user's and/or patient's Member State.

6- LIST OF REFERENCES

17456/23	ILM 360 FORCEPS 23G
17456/25	ILM 360 FORCEPS 25G
17456/27	ILM 360 FORCEPS 27G
17457/23	END GRASPING 360 FORCEPS 23G
17457/25	END GRASPING 360 FORCEPS 25G
17457/27	END GRASPING 360 FORCEPS 27G
17458/23	ECKARDT 360 FORCEPS 23G
17458/25	ECKARDT 360 FORCEPS 25G
17458/27	ECKARDT 360 FORCEPS 27G
17459/23	ULTRA GRIP 360 FORCEPS 23G
17459/25	ULTRA GRIP 360 FORCEPS 25G
17459/27	ULTRA GRIP 360 FORCEPS 27G
17460/23	ASYM TANO 360 FORCEPS 23G
17460/25	ASYM TANO 360 FORCEPS 25G
17460/27	ASYM TANO 360 FORCEPS 27G
17461/25	YONEDA 360 FORCEPS 25G
17461/27	YONEDA 360 FORCEPS 27G

7- EXPLANATION OF SYMBOLS (General list of symbols used by Moria)

Symbol	Meaning
	CATALOG REFERENCE
	QUANTITY
	USE BY DATE
	BATCH CODE
	STERILIZED USING ETHYLENE OXIDE
	DO NOT REUSE
	MANUFACTURER
	DATE OF MANUFACTURE
	CONSULT OPERATING INSTRUCTIONS
	CAUTION: CONSULT ACCOMPANYING DOCUMENT(S)
	DO NOT USE IF PACKAGE IS DAMAGED
	KEEP DRY
	KEEP AWAY FROM SUNLIGHT
	CAUTION FOR USA ONLY: US FEDERAL LAWS RESTRICT THIS DEVICE TO SALE BY, OR ON THE ORDER OF, A PHYSICIAN.
	MEDICAL DEVICE
	UNIQUE STERILE BARRIER SYSTEM WITH INNER PROTECTIVE PACKAGING

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