ADDENDUM regarding the M2 LC (Large Cut) suction ring

Description

The M2 LC-1 suction ring #19379/-1 is a complement to the existing standard M2 suction rings #19325/XX. It allows the surgeon to cut larger diameter flaps than those obtained using standard rings.

Standard and large cut rings are used in the same manner. Thus the M2 LC-1 suction ring is compatible with both M2 (reusable) heads and M2SU (Single Use) heads.

For more information regarding set up, maintenance, sterilization etc., refer to the most recent version of the M2 instruction manual #65013 available on Moria's website: www.moria-surgical.com (Users' Club heading).

Performance

The corneal flap diameter obtained with the M2 LC-1 suction ring is 0.5 mm larger than that obtained on the same eye with the standard -1 M2 suction ring.

Indications

Moria recommends not to use the M2 LC-1 suction ring on corneas with steepest keratometry greater than 44.5D.

Nomogram

No specific stop adjustment is required with the M2 LC-1 suction ring compared to the -1 M2 suction ring; see the indicative nomogram.

Reminder regarding the nomogram

The nomogram provided by Moria is a guideline only, and the results obtained may differ from those mentioned in the nomogram. A variety of factors have an influence on the surgical result:

- Intra-ocular pressure
- Degree of myopia
- Age
- Gender
- Ethnic group
- Corneal hydration
- Suction time
- Pre-operative topography
- Surgical technique

Moria strongly recommends that surgeons develop their own nomogram based on analysis of their post-operative data.

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MORIA Inc