#65086

EN

Instructions for use for the ONE® Single-Use Adjustable Vacuum Trephine # 17202Dxxx

Indication

Moria Single-Use trephines are devices intended to incise a circular disc with partial or full thickness on the patient's cornea or on donor corneal tissue.

Device components

The device consists of a trephine with an adjustment ring (Fig. 1) and suction unit (Fig. 2). The adjustment ring (E) is used to set the depth of the cut.

The most recent version of this document is available on Moria website (<http://www.moria-surgical.com>) via the restricted access heading «Documentation/Library». Translations are also available in other languages under that same heading. Registration required.

Warning

- These instructions explain only how to use this medical device and do not describe the steps necessary to perform a corneal graft. The surgical procedure remains the responsibility of the surgeon experienced in performing corneal grafts.
- The trephine is sterile and ready for use. Do not use it if the pouch in which it is packaged is damaged.
- Do not dismantle the trephine, as this may damage the device.
- The blade is sharp. Do not touch it.
- The trephine is for single-use only: it is not designed to be resterilized or reused under any circumstances.
- The trephine is indicated for use with restriction (Child, pregnant and/or breast-feeding are not eligible for Ophthalmic surgery without doctor's consent)
- At-risk populations (pediatric population, pregnant and breastfeeding women etc...) eligibility is subject to surgeon's judgment and experience

Directions

1/ Select the appropriate diameter and assemble the device

The diameter is indicated on the packaging and on the trephine (Fig.3).
Connect the tubing (G, J) to the syringe (Fig.4).

2/ Check before using

The adjustment ring (E) is preset in «0» position. The blade is retracted at 250µm above this value (i.e. one full turn).

To control the blade: Put the trephine as in figure 6a, set the stop in maximal position (marked by the black zone). Turn over the trephine so it is positioned as in figure 7 (blade downward), then move the blade down by turning the collar, clockwise until the stop. Check under the microscope (minimum 20x) that the cutting edge of the blade has no visible defect. Defects can compromise the performance of the cut. Retract the blade fully in the trephine before placing it on the patient's eye.

To set the blade at the «0» position: check the position of the blade by moving it down until the stop (value «0») then retract it a maximum of one complete turn; this is to avoid touching the cornea while positioning the trephine on eye.

To adjust the depth of blade descent: first perform pachymetry on the patient's eye at the apex of the blade diameter (not only at the center). Pachymetry performed by the surgeon gives an approximate value and should be taken into account when the blade descends (Fig 5).

Turn the adjusting ring to display the depth in the viewing window (A) (Fig. 6a)

| Adjustment setting | 0 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | X | — |
|--------------------|---|-----|-----|-----|-----|-----|-----|-----|-----|------|-------|
| Depth* (µm) | 0 | 200 | 300 | 400 | 500 | 600 | 700 | 800 | 900 | 1000 | >1200 |

* Depth of blade descent in relation to 0. The depth of the cut depends on patient parameters (pressure exerted on eye, keratometry, pathologies, ...).

The setting « 10 » is indicated by « X » and the maximal position by the large black zone. The maximal value is reached as soon as the black zone is visible in the viewing window. Intermediate points between two settings represent increments of 50µm. (Fig. 6a and 6b)

3/ To use the device

Center the trephine on the cornea: Mark the center of the cornea, then position the trephine using the centration indicator (C).

Check that the device is clear of interference from other devices such as speculums. Check that the cornea is fully flattened under the centration indicator. In case of severe keratoconus with marked protrusion of the cone, the trephination can be improved by flattening the apex of the cone.

Apply suction: Depress the syringe piston completely (F), then press the trephine gently and uniformly onto the cornea. Slowly release the syringe piston.

Confirm suction: Pull upward on the trephine very gently. It should remain firmly attached to the cornea.

Trephining:

- hold the trephine between thumb and index finger on the grip area (H) (Fig.7)
- slowly turn the collar (I) clockwise (Fig.8) to cause the blade to descend. The blade descent will be stopped by the adjustment ring setting. If the adjustment ring setting exceeds the corneal thickness, the descending trephine will perforate the cornea.
- if desired, the reference on the collar (I) can be used to count the number of turns made (Fig.9). Each complete turn corresponds to a blade (D) descent of 250µm.

Removing the trephine: Depress the syringe piston to release suction and remove the trephine vertically.

The base of the trephine in contact with the cornea has 16 equidistant marks (B) that can be used for the remainder of the procedure. At the end of the procedure, retract the blade (D) and place the trephine and its accessories in an appropriate container for destruction (in compliance with applicable regulations).

Explanation of symbols (General list of symbols used by Moria)


| | | |
|--|---|--|
|  REF xxxxxx | CATALOGUE REFERENCE | |
|  | QUANTITY | |
|  YYYY/MM/DD | USE BY | |
|  LOT xxxxxx | BATCH CODE | |
|  STERILE EO | STERILIZED USING ETHYLENE OXIDE | |
|  | DO NOT REUSE |  DO NOT RESTERILIZE |
|  xxxxxxxx | MANUFACTURER | |
|  YYYY/MM/DD | DATE OF MANUFACTURE | |
|  MD | MEDICAL DEVICE | |
|  i xxxxxxx | CONSULT OPERATING INSTRUCTIONS | |
|  | DO NOT USE IF PACKAGE IS DAMAGED | |
|  | PROTÉGER DE LA LUMIÈRE DU SOLEIL | |
|  | KEEP DRY | |
|  | DISCARD IN A SEPARATE COLLECTION BIN | |
|  Rx only | CAUTION for USA only: US federal laws restrict this device to sale by, or on the order of, a physician. | |

List of references

| | |
|----------------------|---------------------------------|
| Reference: 17202Dxxx | |
| 17202Dxxx | ADJUSTABLE VACUUM TREPHINE Dxxx |
| Packaging | 1 per box/Sterile/Single-Use |

Accessories

| | |
|-------------|--|
| Reference | 17200Dxxx |
| Description | ONE® – Disposable corneal vacuum punch Diameters: 6.50 to 9.00 mm with 0.25 mm increments. Also 6.00 and 9.50 mm diameters |
| Packaging | 1 per box/Sterile/Single-Use |
| Reference | 17204 |
| Description | ONE® – Artificial Chamber for Single-Use Trephine This trephine is compatible only with the single-use artificial chamber #17204. |
| Packaging | 1 per box/Sterile/Single-Use |

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MORIA Inc