## Instruction for use of One Use Large-Cut (OU LC) heads (#19184/xxx) for Artificial Chamber

### Indication

OU LC heads (#19184/xxx) for artificial chamber are indicated for the creation of complete lamellar resection (without hinge) for donor graft preparation during keratoplasty surgery.

### Contents

OU LC heads (#19184/xxx) for artificial chamber have a pre-assembled blade, with different calibrations (fig. 1).  
OU LC heads (#19184/xxx) for artificial chamber can be used in combination with:

- a fully autoclavable and hand-driven system containing:
  - the ONE turbine (#19155) associated with its pneumatic hose (#19353)
  - either a reusable ALTK artificial chamber (#19161-19162) with its reusable OU LC guide ring (#19168)
  - or a single-use artificial chamber (#19182) with its single-use OU LC guide ring (#19186)
- or a fully automated system containing:
  - the One Use-Plus (OUP) handpiece for artificial chamber (#19175)
  - either a reusable ALTK artificial chamber (#19161-19162) with its reusable OUP ALTK guide ring (#19173)
  - or a single-use artificial chamber (#19182) with its single-use OUP LC guide ring (#19187)
- Evolution 3E control unit (#19380); to activate the ONE turbine (#19155) or OUP handpiece (#19175)
- infusion tubing set by MORIA (#19181); to establish pressure inside artificial chamber.

All elements to be used are MORIA components.

### Warning

Use the OU LC heads for artificial chamber exclusively with accessories manufactured by MORIA for this microkeratome; use of other accessories can damage the microkeratome and may cause significant injuries.

Read prior to use and depending the case:

- the user manual: - of the ALTK Artificial Chamber for ONE & CBM microkeratomes (#65009)  
- of the Evolution 3E control unit (#65060/INTL)
- addendums related to the guide rings:
  - OU LC guide ring for ALTK chamber (#65104)
  - OUP ALTK guide ring (#65105)
- addendum related to the OUP handpiece (#65106).

- This instruction for use only refers to the use of medical device but does not describe surgical steps to perform keratoplasty. Graft surgery must be performed only by experienced corneal surgeons.
- Check the integrity of the pouches.
- Check the expiration date on the labels.
- Inspect the blade under a microscope.
- The OU LC head for artificial chamber is sterile and for single use only. Do not use if packaging is damaged.
- **The reuse of single-use products, or the use of consumables other than those supplied by MORIA, may entail serious surgical consequences for the patient and damage the microkeratome.**
- **Design of this medical device does not allow its reuse. Reusing this medical device affects its clinical performances and exposes the patient to adverse events.**

For USA only: CAUTION : Federal law (USA) restricts these devices to sale by or on the order of a physician.

### Performance

Multiple parameters affect cap thickness and standard deviation: keratometry, pachymetry, graft quality, donor age, corneal edema, corneal storage medium, pressure inside the artificial chamber, speed of the pass, etc.

MORIA can only guarantee the dimensional characteristic of the head, not the surgical result.

Refer to document "One Use Large-Cut Head Selection for Creating Ultra-thin Corneal Lamellar Grafts" (#65101). That is a general indicative guide only. MORIA recommends that every user practice with non-transplantable corneal tissue to determine his own performance. Please read the user manual prior to performing actual surgery.

## Instructions for use

### 1/ Preparing the Evolution 3E control unit:

- Install the control unit according to instructions within user manual (#65060/INTL).
- Select manual mode (ONE turbine) or OUP (One Use-Plus handpiece).
- Proceed to vacuum test.
- Connect flexible hose or OUP handpiece to the control unit.
- In case the OUP handpiece is used: select Speed 2 (advance speed).

### 2/ Preparing the artificial chamber and connecting to infusion tubing :

- Connect the artificial chamber to short line from infusion tubing set (amortizing line)(fig 2 and 3).
- Raise the IV bottle to a height of 1.20 m between the drip chamber and the top of the artificial chamber (cornea) to ensure obtaining pressure of 90 mmHg (fig. 3 and 4).
- Attach the IV bottle to the long line (IV infusion line).
- Open the air valve. Purge the line, and use the roller to close the system. Check that there is no air bubble in the IV line up to the chamber.
- Attach the IV line to the short line (fig. 2 and 3) and release the roller of the IV line while keeping the chamber stopcock closed.

### 3/ Positioning the graft on the artificial chamber:

RECOMMENDATION: Use grafts with a rim having a minimum diameter of 16 mm and a well-centered cornea.

- If desired, apply cohesive viscoelastic to the endothelium.
- Open the IV line.
- Position the donor graft on the top of the piston, endothelium downwards.
- Slightly lift the graft with forceps to remove the possible air bubbles.
- Place the cover carefully in position (inserting the pin into the slot). Check that the graft is well centered. Seal the chamber:
  - single-use chamber: by placing the white locking ring (B) over the cover (A) (fig. 5)
  - ALTK chamber: by screwing the piston (C) to the cover (D) (fig. 8)
- Open the stopcock at least three seconds. Then, close the stopcock (between short line and IV line) to maintain pressure of 90 mmHg (fig 5).
- If desired, check the pressure with a tonometer.
- Place the guide ring on the cover and lock it:
  - single-use chamber:
    - with ONE turbine → screw the single-use guide ring (F) on with the blue screw (G) on the cover (A) (fig. 5 and 7)
    - with OUP handpiece for artificial chamber → screw the single-use guide ring (K) on with the blue screw (G) on the cover (A) (fig. 5 and 7a)
  - ALTK chamber:
    - with ONE turbine → screw the lateral screw (I) in order to lock the reusable guide ring (H) on the cover (D) (fig. 8a and 8b)
    - with OUP handpiece for artificial chamber → remove the diameter adjustment ring from the chamber, then screw the reusable guide ring (J) on the cover (D), and lock the guide ring (fig. 9a and 9b).

### 4/ OU LC head set up:

- Measure the pachymetry and select proper OU LC head according to cap performance. Refer to document "One Use Large-Cut Head Selection for Creating Ultra-thin Corneal Lamellar Grafts" (#65101).
- Screw the OU LC head into:
  - ONE turbine (fig. 6) then connect the turbine to pneumatic hose previously connected to the Evolution 3E control unit
  - OUP handpiece for artificial chamber (fig. 10) already connected to the Evolution 3E control unit.
- While assembled, inspect blade under microscope then activate the ONE turbine or OUP handpiece to check that blade is engaged and oscillates.

### 5/ Graft preparation:

- Open the stopcock at least three seconds. Then, close the stopcock (between short line and IV line) to maintain pressure of 90 mmHg (fig. 5).
- Irrigate both OU LC head and corneal tissue with BSS.
- Proceed to the cut:
  - ONE turbine: activate turbine and proceed to the cut of the graft by sliding the OU LC head into the guide ring for 4-5 seconds in a single motion (fig. 7 and 8)
  - OUP handpiece for artificial chamber: press vacuum footpedal then press forward footpedal to proceed to the cut of the graft.
- After use, unscrew the head.
- If willing to proceed to a second pass:
  - ONE turbine: rotate 180° the chamber and repeat actions §3-4
  - OUP handpiece for artificial chamber: unscrew the reusable guide ring (J), rotate 180° the guide ring, screw again and repeat actions §3-4.

### 6/ End of procedure:

After use, discard everything into an appropriate container for disposal (according to the regulation in force).

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MORIA Inc