
IMPORTANT: Check the integrity of the pouches, as well as the expiration date on the labels. The medical device is sterile and for single use only. It should not be resterilized nor reused. Design of this medical device does not allow its reuse. Reusing this medical device affects its clinical performances and exposes the patient to adverse events.

To create Ultra-Thin Corneal Lamellar Grafts, please refer to the indicative guide entitled «Creating Ultra-Thin Corneal Grafts with the Moria system Using Practice Donor Tissue» #65084EN.

To select the appropriate combination of disposable CB heads, please refer to the document entitled «CBSU Head Selection for Creating Ultra-Thin Corneal Lamellar Grafts» #65082.

Both documents are available on electronic versions on Moria website www.moria-surgical.com (restricted area - «Library» heading). An individual version is also available upon request (Moria's addresses indicated on back page).

Contents

19182 Artificial chamber base with 2 ports (one of which with a one-way stopcock), chamber cover, and white screw

19180 CBSU guide ring and blue screw

19181 Infusion tubing set consisting of one IV line with drip chamber, extension tubing, and one short line with one-way stopcock and syringe

Short line

IV line

Remove all the sterile parts from the pouches. Place the base of the chamber on the table and make sure it is stable.

Preparing the lines to obtain 90 mmHg pressure in the artificial chamber:

Attach the short line to the chamber (to the port without stopcock). Check that the 2 stopcocks are open.

Fill the syringe with the corneal storage medium or physiologic solution (BSS® or similar). Attach the syringe to the short line. Place a finger on the top of the piston and fill the line completely (from the short line to the end of the piston and fill the line completely (from the short line to the end of the chamber tubing). Close the chamber stopcock.

Height between the liquid in the drip chamber and the top of the artificial chamber: around 1.20 m


Raise the IV bottle to a height of 1.20 m between the drip chamber and the top of the artificial chamber (cornea) to ensure obtaining pressure of 90 mmHg.

Attach the IV bottle to the IV line. Open the air valve. Purge the line, and use the roller to close the system. Check that there is no air bubble in the IV line up to the chamber.

Attach the IV line to the chamber stopcock and release the roller of the IV line while keeping the chamber stopcock closed.

Positioning the graft on the artificial chamber:

If desired, apply cohesive viscoelastic to the endothelium.

Use the syringe to gently infuse the system with corneal storage medium or BSS®. Position the donor graft on the top of the piston, endothelium downwards. Slightly lift the graft with forceps to remove the possible air bubbles.

Precisely center the cornea on the stage and place the cover carefully in position (inserting the pin into the slot). Check that the graft is well centered. Seal the chamber by placing the white locking ring over the cover.

Open the chamber stopcock at least three seconds. Then, close the chamber stopcock to maintain pressure of 90 mmHg. **Check the 2 stopcocks are closed.**

If desired, check the pressure with a tonometer.

RECOMMENDATION: Use grafts with a rim having a minimum diameter of 16mm and a well-centered cornea.

Donor graft cut:

Main steps:

- remove the epithelium and dry the corneal surface
- measure corneal thickness
- if desired, check the pressure with a tonometer
- place the guide ring for CBSU head on the cover and screw it on with the blue screw
- make the cut(s) with the CB turbine and the appropriate CBSU heads (see document "CBSU Head Selection for Creating Ultra-thin Corneal Lamellar Grafts" #65082)


Removal of the graft and disposal of the devices:

Once the cuts are completed, unscrew the guide ring.

To prevent a possible collapse of the graft, open the chamber stopcock, unscrew the white screw, and gently lift the cover to grasp the graft.

The graft is now ready to be placed on a Moria punch to obtain the desired diameter.

At the end of the surgical procedure, discard the single-use artificial chamber and its accessories, the CBSU head, and the single-use tubing into an appropriate container for disposal (according to the regulations in force).

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MORIA Inc
(CHINA) CO., LTD.