

Single-Use
Instruments



NEW

Moria Single Use ROP Examination Kit

Now available in **EasySet** exam **Tray**



Moria ROP Exam Kit...
The Right Solution for the NICU



Supporting best practices with the ROP Kit

American Academy of Pediatrics (AAP) published a revised policy statement, *Screening Examination of Premature Infants for Retinopathy of Prematurity*¹. Included is a recommendation that "sterile instruments should be used to examine each infant to avoid possible cross-contamination of infectious agents."

The importance of using single-use instruments during the ROP examination has now been well proven²⁻³ in order to reduce the risk of infection and cross-contamination in premature infants.

Through a continuous innovation, Moria is introducing the EasySet for ROP Examination so as to support best practices.

Bringing simplicity and confidence



- ✓ Sterile
- ✓ Ergonomic and stable tray
- ✓ Good sterile support for retinal exam lens or eye drops
- ✓ High-quality surgical material. Not made with natural rubber Latex, DEHP or PVC


The tray is safe and helpful during examination.

17257X10



NEW Loop Depressor

- Gentle curve to facilitate scleral depression without slippage
- Ergonomic design for the infant
- High-quality finish



NEW Pediatric Speculum design

- Interdigitating blades for easy insertion between eyelids
- Grasping discs - well positioned for positive feel when wearing gloves
- Excellent spring tension to reduce speculum rotation
- Low profile design that can also be used with pediatric eye imaging (e.g., RetCam™)

■ Cost Effective for the NICU

Single-use ophthalmic ROP EasySet is in a sterile peel-pack. 10 kits per box.
For more information, please contact your local representative.

1. From the American Academy of Pediatrics: Policy Statement: Screening Examination of Premature Infants for Retinopathy of Prematurity AMERICAN ACADEMY OF PEDIATRICS Section on Ophthalmology, AMERICAN ACADEMY OF OPHTHALMOLOGY, AMERICAN ASSOCIATION FOR PEDIATRIC OPHTHALMOLOGY AND STRABISMUS, and AMERICAN ASSOCIATION OF CERTIFIED ORTHOPTISTS. Pediatrics 2013; 131:1 189-195; published ahead of print December 31, 2012, doi:10.1542/peds.2012-299615.
2. Hered RW. Reducing infection risk from retinopathy examinations: a new single-use sterile instrument set. Poster 2007, Washington, DC, USA
3. Hered RW. The retinopathy of prematurity screening examination: ensuring a safe and efficient examination while minimizing infant discomfort. Neonatal Network 2010;29(3):143-151

To obtain more information



MORIA SA
15, rue Georges Besse
92160 Antony
FRANCE
Phone: +33 (0) 1 46 74 46 74
Fax: +33 (0) 1 46 74 46 70
www.moria-surgical.com

MORIA Inc
1050 Cross Keys Drive
Doylestown, PA 18902
USA
Phone: (800) 441 1314
Fax: + 1 (215) 230 7670
www.moria-surgical.com

The Single Use ROP Examination Kit is patented by the Nemours Foundation (US Patent #6440065) and is manufactured and distributed solely by Moria under exclusive license from the Nemours Foundation.
RetCam™ is a registered trademark of Clarity Medical Systems, Inc.

Moria
Leading Innovations in Ophthalmology

