



LASIK Surgery

Open the door to
customized flaps using

One Use-*Plus* SBK



Unrivalled features...

✓ Unrivalled ease of use

Automated & linear motion

✓ Unrivalled ease of handling

Pre-assembled & one-handed usage possible

✓ Unrivalled safety

Fits on every patient's eye:

- size ring $\Rightarrow$ flap diameter
- stop value $\Rightarrow$ hinge length

✓ Unrivalled reliability

Two independent motors:

- one for head advancement
- one for blade oscillation

✓ Unrivalled simplicity

- single-use head
- single-use ring option

✓ Unrivalled convenience

Intraop visibility for confirmation of suction




✓ Ring size $\Leftrightarrow$ customization of flap diameter





✓ Stop value $\Leftrightarrow$ customization of hinge length

✓ Cutting head $\Leftrightarrow$ customization of flap thickness

✓ Speed 1 or Speed 2




... for unrivalled performance





One Use-Plus **SBK**

SBK with Moria



- ✓ Predictable thin flaps
- ✓ Planar architecture
- ✓ Beveled edges
- ✓ Customized flap diameters
- ✓ Unrivalled bed smoothness
- ✓ Accuracy and predictability equivalent to Femto-SBK...
- ✓ ...with faster visual recovery
- ✓ ...with greater patient comfort
- ✓ ...with a high level of safety
- ✓ ...at a fraction of the cost

To obtain more information



MORIA SA
15, rue Georges Besse
92160 Antony
FRANCE
Phone: +33 (0) 1 46 74 46 74
Fax: +33 (0) 1 46 74 46 70
www.moria-surgical.com

MORIA Inc
1050 Cross Keys Drive
Doylestown, PA 18902
USA
Phone: (800) 441 1314
Fax: +1 (215) 230 7670
www.moria-surgical.com

MORIA JAPAN KK
Inami Bldg. 6F
3-24-2 Hongo Bunkyo-ku
Tokyo 113-0033
JAPAN
Phone: +81 3 5842 6761
Fax: +81 3 5842 6762
www.moriajapan.com

MORIA in CHINA
Moria Commercial (China) Co. Ltd
Rm H. I, 6FL, Kai Li Building
N° 432 Huai Hai West Road
Shanghai 200052
CHINA
Phone: +86 21 5258 5066
Fax: +86 21 5258 5067
www.moria-surgical.com.cn

Moria
Leading Innovations in Ophthalmology