


**Single-Use
Instruments**

Moria Disposable KIT IVT Intravitreal injection


Speculum with tabs & Scleral Marker for
Intravitreal injection



MADE IN FRANCE




Moria
Ophthalmic Instruments



Speculum with Solid Blades/ Tabs

- Solid Blades to keep away eyelashes
- Curved blades for better comfort
- Discs/ Tabs for easy manipulation
- Good spring tension
- Full metal for a better rigidity
- Sterile, ready for use



Scleral marker

- Designed to mark the point for any pars plana incision or injection point
- Two calibrated markers 3.5mm/ 4.0mm
- Sterile, ready for use

References	Designations
17542X10	IVT KIT SPECULUM SOLID BLADE/MARKER
17538X10	SPECULUM SOLID BLADE/TABS
17539X20	SCLERAL MARKER

To obtain more information

MORIA SA
15, rue Georges Besse
92160 Antony
FRANCE
Phone: +33 (0) 1 46 74 46 74
Fax: +33 (0) 1 46 74 46 00
www.moria-surgical.com



MORIA Inc
1050 Cross Keys Drive
Doylestown, PA 18902
USA
Phone: (800) 441 1314
Fax: +1 (215) 230 7670
www.moria-surgical.com



Moria
Ophthalmic Instruments