


## STAIN CLEARLY WITH TWIN and viewILM

**Peel with Moria flat handle precision**  
**Confidence in every peel**



# Moria

**MORIA SA**  
27, rue du Pied de Fourche  
03160 Bourbon L'Archambault  
FRANCE  
Tél. : +33 (0)1 46 74 46 74  
Fax : +33 (0)1 46 74 46 00  
[moria@moria-int.com](mailto:moria@moria-int.com)  
[www.moria-surgical.com](http://www.moria-surgical.com)

**ALCHI.M.I.A. S.R.L.**  
(a Moria company)  
Viale Austria, 14  
35020 - Ponte San Nicolò (PD)  
ITALY  
Phone: +39 049 89.62.074  
Fax: +39 049 89.62.071  
[info@alchimiasrl.com](mailto:info@alchimiasrl.com)  
[www.alchimiasrl.com](http://www.alchimiasrl.com)

CE 0051  
CE 0477

**MORIA Inc**  
1050 Cross Keys Drive  
Doylestown, PA 18902  
USA  
Phone: (800) 441 1314  
Fax: +1 (215) 230 7670  
[orders@moriausa.com](mailto:orders@moriausa.com)  
[www.moria-surgical.com](http://www.moria-surgical.com)

**Moria Surgical Spain**  
C. de Pablo Iglesias, 58,  
08908 L'Hospitalet de Llobregat  
Barcelona - SPAIN  
Phone: +34932896270  
[info@moria-spain.com](mailto:info@moria-spain.com)  
[www.ofalteltech.com](http://www.ofalteltech.com)

**Moria Japan K.K.**  
Arcadia Building 6F  
1-12-3 Kanda  
SudachoChiyoda-Ku  
Tokyo 101-0041  
JAPAN  
Phone: 81-3-6260-8309  
Fax: 81-3-6260-8310  
[moria@moriajapan.com](mailto:moria@moriajapan.com)  
[www.moriajapan.com](http://www.moriajapan.com)

**ALCHI.M.I.A. S.R.L.**  
(a Moria company)  
Viale Austria, 14  
35020 - Ponte San Nicolò (PD)  
ITALY  
Phone: +39 049 89.62.074  
Fax: +39 049 89.62.071  
[info@alchimiasrl.com](mailto:info@alchimiasrl.com)  
[www.alchimiasrl.com](http://www.alchimiasrl.com)

**Moria COMMERCIAL  
(CHINA) CO., LTD.**  
上海目利亚贸易有限公司  
Room H 6F Kailli Building NO.432  
West Huai Hai Road, Changning  
district, Shanghai, 200052, P.R.C  
中国上海市长宁区淮海西路432号  
凯利大厦6楼08室  
Phone/Fax : +86 021 52586095  
[moria@moria-int.com](mailto:moria@moria-int.com)

**Moria Limited**  
Suite 1904, 19th Floor,  
ST. George's Building,  
No. 2 Ice House Street,  
Central - HONG KONG  
[moria@moria-int.com](mailto:moria@moria-int.com)  
[www.moria-surgical.com](http://www.moria-surgical.com)



MORIA #68921EN-A-06.2026

# Moria

**See it, Grip it, Peel it**

**The complete solution  
for a trusted ILM / ERM peel**



**1 Stain with precision**  
Stain the **ILM / ERM** with  
high **purity dyes** for clear,  
high-contrast visualization

**2 Peel with confidence**  
Use **retinal forceps** designed  
for optimal control and safety

RETINAL DYES

TWIN  
for ILM / ERM staining



viewILM  
for ILM staining



Maximum cohesiveness  
Enables precise targeting of stain



Excellent contrast  
Clearly highlights membranes for reliable identification



Cytotoxicity tested  
Batch level testing – ensuring safety for retinal surgeries  
Blulife® purity > 99%.



Ready to use Pre-filled syringes  
designed for convenience, precision, and safety

| COMPOSITION               | TWIN                           | viewILM    |
|---------------------------|--------------------------------|------------|
| Reference                 | RMB 004-00                     | RMB 003-00 |
| Blulife® concentration    | 0.03%                          |            |
| Trypan Blue concentration | 0.18%                          | -          |
| pH                        | 7.00 – 7.40                    |            |
| Osmolality                | 340 – 365 mOsm/kg              |            |
| Endotoxins                | < 0.5 EU/ml                    |            |
| Shelf life                | 36 months                      |            |
| Pack                      | 6 pre-filled syringes (0.7 ml) |            |

PRECISION TO REACH EVERY MEMBRANE

A complete range of disposable forceps designed to address every posterior retinal membrane peeling need



Precise & Controlled  
• Optimal rigidity for stability  
• Gentle, controlled peeling



Smooth & Safe  
• Atraumatic tip closure  
• Reliable grip and reactivity



Ergonomic & Practical  
• Dual-sided actuation  
• Colour coded gauge indicator  
• Low-profile, non-reflective shaft for better visibility

| Models                                     |     | ILM                               | End-Grasping            | Eckardt                 | Tano                    | Ultra-Grip                            | Yoneda         |
|--------------------------------------------|-----|-----------------------------------|-------------------------|-------------------------|-------------------------|---------------------------------------|----------------|
| Reference<br>(Flat-Handle Retinal Forceps) |     | 17450<br>(23 / 25 / 27)           | 17554<br>(23 / 25 / 27) | 17452<br>(23 / 25 / 27) | 17454<br>(23 / 25 / 27) | 17453<br>(23 / 25 / 27)               | 17455<br>(27)  |
| Grasping strength                          |     | 23G 25G ≥ 50 g ••<br>27G ≥ 30 g • |                         |                         |                         | 23G 25G ≥130G ••••<br>27G ≥ 80 g •••• | 27G ≥ 80 g ••• |
| Working length<br>(mm)                     | 23G | 34.6                              | 34.5                    | 35.3                    | 33.5                    | 33.8                                  | -              |
|                                            | 25G | 32.6                              | 32.5                    | 32.0                    | 31.4                    | 31.7                                  | -              |
|                                            | 27G | 29.5                              | 29.7                    | 30.0                    | 29.2                    | 29.4                                  | 29.4           |
| Pack                                       |     | 5 units / box                     |                         |                         |                         |                                       |                |