


OUR SOLUTIONS FOR BEST PRACTICES



The American Academy of Pediatrics (AAP) published a recommendation that "sterile instruments should be used to examine each infant to avoid possible cross contamination of infectious agents."¹

The following graph shows the significant reduction of the risk of infection that can be achieved by following best practice.²



1. From the American Academy of Pediatrics: Policy Statement: Screening Examination of Premature Infants for Retinopathy of Prematurity AMERICAN ACADEMY OF PEDIATRICS Section on Ophthalmology, AMERICAN ACADEMY OF OPHTHALMOLOGY, AMERICAN ASSOCIATION FOR PEDIATRIC OPHTHALMOLOGY AND STRABISMUS, and AMERICAN ASSOCIATION OF CERTIFIED ORTHOPTISTS. *Pediatrics* 2013; 131:1 189-195; published ahead of print December 31, 2012, doi:10.1542/peds.2012-299615.
2. Hered RH. Use of nonsterile instruments for examination for retinopathy of prematurity in the neonatal intensive care unit. *Journal of Pediatrics*. 2004 September; 145: 308-311




Disposable ROP Exam Kit with

Our support for your best practices




#17257X10
ROP exam kit W/WS-L
 Pediatric Wire speculum with tabs & Loop depressor

- Sterile
- Ergonomic and stable tray



#17265X10
ROP KIT S SPECULUM/LOOP DEPRESSOR
 Pediatric Wire speculum with tabs & Loop depressor

- Sterile & ergonomic design
- **New** compact packaging with less plastic




Pediatric speculum

- Easy to insert and remove
- Secure grasping
- Appropriate size and flex characteristics for premature infants
- Low profile design that can be used with pediatric eye imaging



Loop depressor

- Gentle curve to facilitate scleral depression
- Ergonomic design for infants
- High-quality finishing




Double ended depressor

- For 4 distinct orientations
- Ergonomic design for infants
- High-quality finishing


Reference	Designations
17257X10	ROP EXAM KIT W/WS-L Pediatric Wire speculum with tabs & Loop depressor
17265X10	ROP KIT S SPECULUM/LOOP DEPRESSOR Pediatric Wire speculum with tabs & Loop depressor
17263X10	ROP KIT S SPECULUM/2IN1 DEPRESSOR Pediatric Wire Speculum with tabs & Double ended Depressor

NEW



ROP kits with small neonatal speculum

- 4.1 mm wire blades
- Staggered closed wire blades to fit together, for easier insertion and removal
- Appropriate size and flex characteristics for premature infants
- Developed with surgeons



Reference	Designations
17264X10	ROP KIT XS SPECULUM/LOOP DEPRESSOR Small Wire Speculum & Loop Depressor
17262X10	ROP KIT XS SPECULUM/2IN1 DEPRESSOR Small Wire Speculum & Double ended Depressor

Single-Use Neonatal Eyelid Speculum & Scleral Depressor Kit used by leading NICUs in more than 37 states reduces the risk of cross-contamination and healthcare-associated infections in retinopathy of prematurity exams, other procedures and injections.



ROP kits with small solid blades speculum

- Solid small blades that limits contact with eyelashes during exam
- Optimal design for easy insertion between eyelids
- Appropriate size and flex characteristics for premature infants
- Low profile that can be used with pediatric eye imaging



Reference	Designations
17255X20	ROP EXAM KIT-L Speculum & loop depressor
17210	ROP EXAM KIT-D Speculum & double-ended depressor

✓ Safe for Infants:

- Reduces risks from improperly sterilized instruments, corneal toxicity and cross contamination
- Made of high-quality surgical materials. Not made with natural rubber Latex, DEHP, or PVC.

✓ Designed by a pediatric ophthalmologist for premature infants with an appropriate speculum size easy to insert and remove

✓ Promoted to ensure best practices

- Facilitates and ensures compliance with AAP, CDC and APIC guidelines
- Enables hospitals and clinics to ensure that sterile instruments are used for every ROP exam

Disposable ROP Exam kit

Our support for your best practices

Reference	Designations
ROP Exam Kit with Pediatric wire Speculum with tabs – Open Blades	
17257X10	ROP EXAM KIT W/WS-L Pediatric Wire speculum with tabs & Loop depressor
17265X10	ROP KIT S SPECULUM/LOOP DEPRESSOR Pediatric Wire speculum with tabs & Loop depressor
17263X10	ROP KIT S SPECULUM/2IN1 DEPRESSOR Pediatric Wire Speculum with tabs & Double ended Depressor
ROP Exam Kit with Small Pediatric Wire Speculum – open blades	
17264X10	ROP KIT XS SPECULUM/LOOP DEPRESSOR Small Wire Speculum & Loop Depressor
17262X10	ROP KIT XS SPECULUM/2IN1 DEPRESSOR Small Wire Speculum & Double ended Depressor
ROP Exam Kit with Small Solid Blades Speculum	
17255X20	ROP EXAM KIT-L Speculum & loop depressor
17210	ROP EXAM KIT-D Speculum & double-ended depressor

 **MORIA SA**
15, rue Georges Besse
92160 Antony
FRANCE
Tel.: +33 (0) 1 46 74 46 74
Fax: +33 (0) 1 46 74 46 00
www.moria-surgical.com

MORIA Inc
1050 Cross Keys Drive
Doylestown, PA 18902
USA
Phone: (800) 441 1314
Fax: +1 (215) 230 7670
www.moria-surgical.com

**Moria**
Ophthalmic Instruments



MORIA #67090US-A-02.2022