


Hand-Held Instruments

New



For easy grasping, insertion, and positioning of donor lamellar button.




Busin Spatula #20006

- For easy insertion of **conventional thickness** donor lamellar button

Mini Busin Spatula #20007

- For easy insertion of **Ultra-thin** donor lamellar button



- Especially designed for "pull-through" technique
- Wing design aids in scooping up graft directly from the punch
- Facilitates unfolding of graft
- Simplifies centration of donor button in the anterior chamber
- Minimizes intraoperative manipulation of graft and possibility of endothelial cell loss



Busin Forceps 20G #20004

- Micro-incision forceps, 20G, distal action
- Tips only contact the periphery of the graft



Busin Forceps 23G #20011

- Micro-incision forceps, 23G, distal action
- Atraumatic ridges to safely grasp the donor button



- Designed to position graft in the Busin spatula and pull it from the spatula into the anterior chamber

To obtain more information



MORIA SA
15, rue Georges Besse
92160 Antony
FRANCE
Phone: +33 (0) 1 46 74 46 74
Fax: +33 (0) 1 46 74 46 70
www.moria-surgical.com

MORIA Inc
1050 Cross Keys Drive
Doylestown, PA 18902
USA
Phone: (800) 441 1314
Fax: +1 (215) 230 7670
www.moria-surgical.com

MORIA JAPAN KK
Inami Bldg. 6F
3-24-2 Hongo Bunkyo-ku
Tokyo 113-0033
JAPAN
Phone: +81 3 5842 6761
Fax: +81 3 5842 6762
www.moriajapan.com

MORIA in CHINA
Moria Commercial (China) Co. Ltd
Rm H, I, 6FL, Kai Li Building
N° 432 Huai Hai West Road
Shanghai 200052
CHINA
Phone: +86 21 5258 5066
Fax: +86 21 5258 5067
www.moria-surgical.com.cn

Moria
Leading Innovations in Ophthalmology