



Moria



Master your rhexis

**Reusable Cross-Action
Capsulorhexis Forceps**



Actual size of the forceps

A COMPLETE RANGE FOR YOUR PREFERENCE



4 tip designs

to match your
preferred technique



7 lengths

from 8.5mm
to 11.5mm



**Curved &
straight**

to shaft options

AVAILABLE TIPS DESIGNS

Sharp pointed tips



Corydon

≈ 90° pointed tip

Ref: **20050/xxx**



Crozafon

≈ 45° pointed tip

Ref: **20051/xxx**

Blunt tips



Utrata

≈ 45° pointed tip

Ref: **20052/xxx**



Serrated

≈ 45° pointed tip

Ref: **20053/xxx**

BENEFITS THAT MATTER



Protect Every Incision: Cross-action design

Minimizes stress on micro-incisions as small as 1.8 mm.



Comfort in Every Movement:

Lightweight bird-cage handle offers excellent balance and handling.



Extended Reach, Precise Control:

Wide jaw opening delivers superior access without enlarging the incision.



Built to Last:

Durable, integrated design supported by Moria repair service.



Gentle by Design:

Rounded reinforced arms help ensure smooth, atraumatic performance.

KEY ENGINEERING FEATURES



Reinforced rounded arms



Precision tip coaptation



Dedicated protective cap included



Premium stainless steel construction

Moria

MORIA SA
27, rue du Pied de Fourche
03160 Bourbon L'Archambault
FRANCE
Tél. : +33 (0)1 46 74 46 74
Fax : +33 (0)1 46 74 46 00
moria@moria-int.com
www.moria-surgical.com

CE 0459

Moria Inc
1050 Cross Keys Drive
Doylestown, PA 18902
USA
Phone: (800) 441 1314
Fax: +1 (215) 230 7670
orders@moriausa.com
www.moria-surgical.com

Moria Surgical Spain
C. de Pablo Iglesias, 58,
08908 L'Hospitalet de Llobregat
Barcelona - SPAIN
Phone: +34932896270
info@moria-spain.com
www.oftaltech.com

Moria Japan K.K.
Arcadia Building 6F
1-12-3 Kanda
SudachoChiyoda-Ku
Tokyo 101-0041
JAPAN
Phone: 81-3-6260-8309
Fax: 81-3-6260-8310
moria@moriajapan.com
www.moriajapan.com

ALCHIMIA S.R.L.
(a Moria company)
Viale Austria, 14
35020 - Ponte San Nicoló (PD)
ITALY
Phone: +39 049 89.62.074
Fax: +39 049 89.62.071
info@alchimiasrl.com
www.alchimiasrl.com

**Moria COMMERCIAL
(CHINA) CO., LTD.**
上海目利亚贸易有限公司
Room H 6F Kaiili Building NO.432
West Huai Hai Road, Changning
district, Shanghai, 200052, P.R.C
中国上海市长宁区淮海西路432号
凯利大厦6楼08室
Phone/Fax : +86 021 52586095
moria@moria-int.com

Moria Limited
Suite 1904, 19th Floor,
ST. George's Building,
No. 2 Ice House Street,
Central - HONG KONG
moria@moria-int.com
www.moria-surgical.com



MORIA #67092EN-B-062026