

The New State-of-the-Art

Introducing
360° Vitreo-Retinal
Peeling Forceps



360° handle

- ▲ **Low pressure** required to actuate the handle
- ▲ Progressive closure & short stroke for an **optimal control**
- ▲ Color coding for easy **23G-25G-27G identification**

Refined active part

- ▲ **Optimal grasping strength** to prevent from slipping
- ▲ **Fixed proximal part** to prevent backward movement when in contact with trocar

**PATENT
PENDING**

Optimal shaft stiffness

- ▲ **System+** to further enhance it with 27G

Made in
 **France**



► 16 new references available

REFERENCE	DESIGN AND SIZE	USEFUL LENGTH	GRASPING FORCE
17456/23	ILM 360 forceps 23G	34.6 ± 0.5 mm	≥ 50 g
17456/25	ILM 360 forceps 25G	32.5 ± 0.5 mm	≥ 50 g
17456/27	ILM 360 forceps 27G	29.5 ± 0.5 mm	≥ 30 g
17457/23	END-GRASPING 360 forceps 23G	34.5 ± 0.5 mm	≥ 50 g
17457/25	END-GRASPING 360 forceps 25G	32.5 ± 0.5 mm	≥ 50 g
17457/27	END-GRASPING 360 forceps 27G	29.7 ± 0.5 mm	≥ 30 g
17458/23	ECKARDT 360 forceps 23G	35.3 ± 0.5 mm	≥ 50 g
17458/25	ECKARDT 360 forceps 25G	32.0 ± 0.5 mm	≥ 50 g
17458/27	ECKARDT 360 forceps 27G	30.0 ± 0.5 mm	≥ 30 g
17459/23	ULTRA GRIP 360 forceps 23G	33.8 ± 0.5 mm	≥ 130 g
17459/25	ULTRA GRIP 360 forceps 25G	31.7 ± 0.5 mm	≥ 130 g
17459/27	ULTRA GRIP 360 forceps 27G	29.4 ± 0.5 mm	≥ 80 g
17460/23	TANO 360 forceps 23G	33.5 ± 0.5 mm	≥ 50 g
17460/25	TANO 360 forceps 25G	31.4 ± 0.5 mm	≥ 50 g
17460/27	TANO 360 forceps 27G	29.2 ± 0.5 mm	≥ 30 g
17461/27	 YONEDA 360 forceps 27G	29.4 ± 0.5 mm	≥ 80 g

► 6 tip designs



ILM



ECKARDT



TANO



END-GRASPING



ULTRA GRIP



YONEDA



Micro-serrated, blunt, rounded wedge with large grasping surface



MORIA SA
27, rue du Pied de Fourche
03160 Bourbon-L'Archambault
FRANCE
Phone: +33 (0)1 46 74 46 74
Fax: +33 (0)1 46 74 46 00
moria@moria-int.com
www.moria-surgical.com

MORIA Inc
1050 Cross Keys Drive
Doylestown, PA 18902
USA
Phone: (800) 441 1314
Fax: +1 (215) 230 7670
orders@moriausa.com
www.moria-surgical.com

Moria

Moria Japan K.K.
Arcadia Building 6F
1-12-3 Kanda
SudachoChiyoda-Ku
Tokyo 101-0041
JAPAN
Phone: 81-3-6260-8309
Fax: 81-3-6260-8310
moria@moriajapan.com
www.moriajapan.com

Moria COMMERCIAL (CHINA) CO., LTD.
Room H 6F Kaili Building NO.432
West Huai Hai Road,
Changning district, Shanghai,
200052, P.R.C
Phone/Fax: +86 021 52586095
moriachina@moria-int.com

