

Moria



See it, Grip it, Peel it

The complete solution
for a trusted ILM / ERM peel



1

Stain with precision

Stain the **ILM** / **ERM** with high **purity dyes** for clear, high-contrast visualization

2

Peel with confidence

Use **retinal forceps** designed for optimal control and safety

RETINAL DYES

TWIN

for ILM / ERM staining



view ILM
for ILM staining



Maximum cohesiveness

Enables precise targeting of stain



Excellent contrast

Clearly highlights membranes for reliable identification



Cytotoxicity tested

Batch level testing – ensuring safety for retinal surgeries
Bluelife® purity > 99%.



Ready to use Pre-filled syringes

designed for convenience, precision, and safety

COMPOSITION	TWIN	view ILM
Reference	RMB 004-00	RMB 003-00
Bluelife® concentration	0.03%	
Trypan Blue concentration	0.18%	–
pH	7.00 – 7.40	
Osmolality	340 – 365 mOsm/kg	
Endotoxins	< 0.5 EU/ml	
Shelf life	36 months	
Pack	6 pre-filled syringes (0.7 ml)	

PRECISION TO REACH EVERY MEMBRANE

A complete range of disposable forceps designed to address every posterior retinal membrane peeling need



Precise & Controlled

- Optimal rigidity for stability
- Gentle, controlled peeling



Smooth & Safe

- Atraumatic tip closure
- Reliable grip and reactivity




Ergonomic & Practical

- Dual-sided actuation
- Colour coded gauge indicator
- Low-profile, non-reflective shaft for better visibility

Models		ILM	End-Grasping	Eckardt	Tano	Ultra-Grip	Yoneda
Reference (Flat-Handle Retinal Forceps)		17450 (23 / 25 / 27)	17451 (23 / 25 / 27)	17452 (23 / 25 / 27)	17454 (23 / 25 / 27)	17453 (23 / 25 / 27)	17455 (27)
Grasping strength		23G 25G ≥ 50 g •• 27G ≥ 30 g •				23G 25G ≥ 130G •••• 27G ≥ 80 g ••••	27G ≥ 80 g •••
Working length (mm)	23G	34.6	34.5	35.3	33.5	33.8	-
	25G	32.6	32.5	32.0	31.4	31.7	-
	27G	29.5	29.7	30.0	29.2	29.4	29.4
Pack		5 units / box					

STAIN CLEARLY WITH TWIN and view ILM

Peel with Moria flat handle precision
Confidence in every peel



Moria

MORIA SA
27, rue du Pied de Fourche
03160 Bourbon L'Archambault
FRANCE
Tél. : +33 (0)1 46 74 46 74
Fax : +33 (0)1 46 74 46 00
moria@moria-int.com
www.moria-surgical.com

ALCHIMIA S.R.L.
(a Moria company)
Viale Austria, 14
35020 - Ponte San Nicolò (PD)
ITALY
Phone: +39 049 89.62.074
Fax: +39 049 89.62.071
info@alchimiasrl.com
www.alchimiasrl.com

CE 0051
CE 0477

MORIA Inc
1050 Cross Keys Drive
Doylestown, PA 18902
USA
Phone: (800) 441 1314
Fax: +1 (215) 230 7670
orders@moriausa.com
www.moria-surgical.com

Moria Surgical Spain
C. de Pablo Iglesias, 58,
08908 L'Hospitalet de Llobregat
Barcelona - SPAIN
Phone: +34932896270
info@moria-spain.com
www.oftaltech.com

Moria Japan K.K.
Arcadia Building 6F
1-12-3 Kanda
SudachoChiyoda-Ku
Tokyo 101-0041
JAPAN
Phone: 81-3-6260-8309
Fax: 81-3-6260-8310
moria@moriajapan.com
www.moriajapan.com

ALCHIMIA S.R.L.
(a Moria company)
Viale Austria, 14
35020 - Ponte San Nicolò (PD)
ITALY
Phone: +39 049 89.62.074
Fax: +39 049 89.62.071
info@alchimiasrl.com
www.alchimiasrl.com

**Moria COMMERCIAL
(CHINA) CO., LTD.**
上海目利亚贸易有限公司
Room H 6F Kaili Building NO.432
West Huai Hai Road, Changning
district, Shanghai, 200052, P.R.C
中国上海市长宁区淮海西路432号
凯利大厦6楼08室
Phone/Fax : +86 021 52586095
moria@moria-int.com

Moria Limited
Suite 1904, 19th Floor,
ST. George's Building,
No. 2 Ice House Street,
Central - HONG KONG
moria@moria-int.com
www.moria-surgical.com



MORIA #6892IEN-A-06.2026